



**Example Protocol for (1,000ml)  
Nanosized Oil in Water Emulsion: CBD**  
  
(Using customer developed surfactant mix)



**Model #Q700CA Sonicator**

## Materials and Equipment:

- CBD isolate
- Purified water
- Mixed surfactant: average HLB of ~11-14 (see Workflow document for details)
- Beakers of various volumes: 100-2000ml
- Calibrated balance
- Weigh boats
- Spatulas
- Stir bar
- Heated stir plate
- Stir bar retriever
- Overhead homogenizer
- Calibrated thermometer / thermometer gun
- #Q700CA Sonicator system
- Compressed air with inlet/outlet
- Ice bath: composition can vary based on user preference
- .45-micron filter
- .22-micron filter
- Pipettes
- Flask weights (optional)
- Parafilm or equivalent (optional)
- Laser pointer (optional)
- Small glass bottles (optional)
- Paper with black print text (optional)

## Formula:

<b>Purified Water</b>	<b>92.5%</b>
<b>Surfactant</b>	<b>5.0%</b>
<b>CBD Isolate</b>	<b>2.5%</b>


**This protocol creates a 25mg/ml concentration.**

\*Common compounding practice allows  $\leq 5\%$  error

### General Preparatory Work:

- Gather appropriate materials and calibrate any balances to be used
- Place 925g of purified water into a 2,000ml beaker
- Place 50g surfactant into a weigh boat or beaker (pre-melt if solid at room temperature)
- Place 25g CBD isolate into a weigh boat or beaker
- Begin heating hot plate to  $\geq 60^{\circ}\text{C}$
- Place stir bar into the beaker containing 925g purified water
- Prepare ice bath and store in cool place until needed

Note: Ice bath should be composed of  $\geq 70\%$  water. It is optional to use alcohols and sodium chloride in conjunction with ice water for a greater cooling effect. Operator should account for volume displacement of ice as well as potential vapor pressures from alcohol evaporation.



925g Purified Water



25g CBD Isolate




50g Surfactant Blend


### Coarse Emulsion:

- Place purified water onto heated stir plate and activate the magnet
- Heat water to  $\geq 60^{\circ}\text{C}$
- Place pre-weighed surfactant into the heated water solution
- Place CBD isolate into the water: surfactant solution  
(Optional step - pre-mix CBD isolate with surfactant then add to heated water)
- Stir on heat for  $\geq 30$  minutes
- Cool to room temperature
- Remove stir bar


Note: do not exceed  $100^{\circ}\text{C}$  as water will begin to vaporize



**CBD added to surfactant before mixing**



**Completed coarse emulsion**



**CBD added to surfactant after mixing**




**CBD added to surfactant after mixing**

## Microemulsion:

- Place homogenizer/mixer into the coarse emulsion
- Homogenize at maximum speed permitted for  $\geq 5$  minutes
- Pipette aliquot into small glass bottle for comparison (optional)



**Overhead homogenization to form microemulsion**



**Completed microemulsion**



**Sample aliquots:  
unfiltered (left), .45 micron (middle), and .22 micron**


## Nanoemulsion Preparatory Work:

- When using 1" probe (#4205), operator must use compressed air to cool the Converter
- See Converter cooling instructions to prevent voiding warranty
- Place ice bath vessel into Sound Enclosure
- Secure beaker in the ice bath with flask weights or clamp (red weights used in photo)
- Close Sound Enclosure door

Note: For best results, ensure ice water volume corresponds to the microemulsion

- Only use a 1,500mL or 2,000mL beaker for a 1,000mL sample volume
- If the probe is too shallow or deep it will not deliver consistent results
- Place cover onto top of beaker to prevent splashing or aerosol (optional)
- Determine appropriate processing parameters

Note: If solution is overheating or foaming excessively, decrease ice bath temperature or increase pulse OFF time. Check probe depth and confirm beaker size.



Note air cooling hoses - one enters the Converter and the other exits the enclosure to exhaust air to the room



Ensure ice bath volume corresponds to sample volume




Program screen for current protocol




## Nanoemulsion:


- Power on the generator using the switch on the back of the unit
- Indicate you are not using a microtip (“Are you using a microtip?” - select NO)
- Select: “To Select or Modify a Program or Sequence, Press HERE”
- Set amplitude to 90
- Set process time to 1 hour (Note: lower concentrations will require less processing time)
- Set pulse ON time to 4 minutes
- Set pulse OFF time to 1 minute
- Save program
- Select RUN
- Closely monitor solution temperature and refresh ice bath as needed
- Once run is completed, open door and allow solution to cool to room temperature




Sample secured in ice bath



Sample appearance at T=30min





**Sample aliquot appearance after .22 micron filter (volume = ~10mL)**

#### **Custom Hardware Processing Parameters:**

While each sample matrix is individual, and most relationships in colloidal chemistry are not linear, below is a table with general guidelines for customizing processing time in relation to oil load assuming 90% amplitude/intensity setting.

<b>Oil Amount</b>	<b>Run Time</b>
1%	~30 min
1.5%	~36 min
2%	~48 min
2.5%	~60 min
3%	~72 min



## Filtering:

Most nanoemulsions must be filtered at least once in order to achieve translucency. Additionally, filtering prevents the presence of probe titanium in the solution.

- Ensure filter matrix is compatible with solution
- Filter sample through .45-micron filter until appearance does not change-it is not uncommon to filter multiple times.
- Filter sample through .22-micron filter until appearance does not change -it is not uncommon to filter multiple times.

When selecting proper filtering matrix, most vendors include a table which lists solvent compatibility.

## Syringe Filters Solvent Compatibility Chart

Group of Substance & Chemical Reagents	Cellulose Acetate	Nylon	PES	PTFE	PVDF
<b>ACIDS</b>					
Acetic, 5%	L	R	R	R	R
Acetic, 10%	L	R	R	R	R
Acetic, 25%	N	L	R	R	R
Acetic, Glacial	N	N	R	R	R
Boric	-	L	-	R	-
Formic 25%	L	N	-	R	-
Hydrochloric 15%	L	L	R	R	L
Hydrochloric 25%	N	N	R	R	-
Hydrochloric concentrated	N	N	L	R	N
Hydrofluoric 10%	N	N	-	-	-
Hydrofluoric 35%	N	N	-	R	-
Nitric 25%	N	N	R	R	-
Nitric 6N, 38%	N	N	L	R	R
Nitric concentrated	N	N	N	R	N
Phosphoric 25%	L	N	R	R	-
Sulfuric 25%	N	N	N	R	-
Sulfuric 6N, 29%	N	N	N	R	-
Sulfuric concentrated	N	N	N	R	N
Trichloroacetic 10%	N	N	-	R	R
<b>ALKALINES</b>					
Ammonium Hydroxide 25%	N	R	R	R	L
Formalin 30%	L	L	R	-	-
Sodium Hydroxide 3N, 12%	N	R	R	R	R



### Top 10 Reasons to Use Restek Syringe Filters

- 1 Protect any analytical system.
- 2 Extend LC column lifetime.
- 3 Achieve more reproducible analyses.
- 4 Variety of membranes, porosities, and diameters available.
- 5 Luer lock inlet provides strong, leak-tight syringe connection to withstand filtration pressure.
- 6 Rugged polypropylene construction— autoclavable to 121 °C for 15 minutes (75 psi).

### Sample Filter Chart from Restek

#### Full Chart:

<https://www.restek.com/pdfs/GNTS2122-UNV.pdf>

## Particle Size Estimation (Optional):

While the most credible way to determine particle size is dynamic light scattering (DLS), there are a few economic methods to check your work.

### Option 1: Printed text test


- Place up to 5mLs of filtered aliquot into a small glass bottle/beaker-beakers often yield best results
- Place beaker/bottle over printed text
- If the particles are nanosized, the text should be extremely easy to read

### Option 2: laser-pointer method


- Place up to 5mLs of filtered aliquot into a small glass bottle/beaker
- Place beaker/bottle in front of a plain colored background such as a lab notebook cover. Use plain white or plain black for best results
- Turn off lights (optional)
- Shine laser beam through solution and check for penetration onto the background.
- This may be done from the anterior or posterior aspect of the solution
- Nanoemulsions will allow light to pass through them




**Unfiltered printed text test**



**.45 micron printed text test**




**.22 micron printed text test**



**Unfiltered laser pointer test**




**.45 micron laser pointer**



**.22 micron laser pointer test**

### **Determine API recovery (Optional):**

The conditions of sonicating in order to form a nanoemulsion can cause side reactions or degradation for a multitude of reasons. The most common being poor condition monitoring (such as inappropriately high temperatures). If you are nanosizing a product for commercial retail, all label claims must be verified with a potency assay via instrumentation. Common methods are liquid and gas chromatography although liquid chromatography is the most used.



**.22 micron laser pointer test in the dark**

**Qsonica, LLC.**  
53 Church Hill Rd.  
Newtown, CT 06470

Phone: 203.426.0101  
[info@sonicator.com](mailto:info@sonicator.com)