

## Cell Lysis

*Example protocols and results are based on customer feedback.*

### Sonicator Model: Q700MPX system with Chiller



#### Protocol:

**Sample Type:** *Bacterial cell suspension in 8M urea 4% SDS*

**Sample Volume:** *100 $\mu$ l*

**Sonicator Amplitude Setting:** *50%*

**Sonication Pulse Rate:** *10 seconds ON, 10 seconds OFF*

**Total Sonication ON Time:** *5-10 minutes*

**Sample Process Temperature:** *8C - controlled by Chiller*

#### Customer Notes:

Bacterial cell samples can be fully lysed in 300  $\mu$ l 96-well PCR plates as well as processed in 200 $\mu$ l volumes in 1.7 ml EP tubes. We had no difference between samples processed using the Microplate Horn when compared to using a Q125 sonicator with single probe.

We recommend rotating the microplate a few times during the procedure to ensure all wells are processed evenly.

Note: The Microplate Horn was also applicable to resuspend proteins during metaproteomic sample processing. This replaces pipettes that work manually. 10% amplitude for approximately 5 seconds was sufficient. 1 ml acetone was used to resuspend around 500  $\mu$ g purified proteins. We also use 100-150  $\mu$ l 6M urea buffer to resuspend the 500  $\mu$ g proteins.