



## Avian Blood and Tissue Genomic DNA Shearing

Example protocols and results are based on customer feedback.


350bp




450bp



750bp



1000bp



### Protocol Information

**Cell Type:** Avian Blood and Tissue samples

**Total Sample Volume:** 100ul

**Sonicator Amplitude Setting:** 25%

**Sonication Pulse Rate:** 15 seconds On, 15 seconds Off

**Total Sonication On Time:** Desired fragment size dependent. See Reference guide below

**Sample Process Temperature:** 3°C

### Fragment Size Reference Guide:

Fragment size (bp)	1000	750	450	350
Amplitude (25)	25	25	25	25
Pulse (on:off)	(15:15)	(15:15)	(15:15)	(15:15)
Duration (Seconds)	60	105	150	180

**Customer Notes:**

- 1) Turn on chiller and set to 3 °C
- 2) Adjust water level in cup horn to match the sample level in the tube.
- 3) Aliquot 100 µl high molecular weight DNA (TE buffer or sterile water) into individual 0.2 ml PCR tubes
- 4) Chill the DNA aliquots to 4 °C
- 5) When chilled, spin down PCR tubes for 30 seconds at 1000 RPM
- 6) Carefully place tubes into sample cradle, making sure not to splash sample onto the tube walls
- 7) Carefully attach loaded sample cradle to the rotator, place into instrument and close the lid
- 8) Sonicate for desired duration, amplitude and pulse
- 9) Sheared DNA may now be used as starting material for library preparation

**Concentration Independent shearing to 750bp**

