

GAIAM WELLBEING™

HOT & COLD NECK CRADLE

Setup, Workout Tips and Care & Safety Guide

Please read before setting up or using your Hot & Cold Neck Cradle.

SETUP

1. Remove Hot & Cold Neck Cradle from box and check for possible shipping damage (scratches, punctures or gouges).
2. Your Hot & Cold Neck Cradle is now ready for use.

INSTRUCTIONS FOR USE

To soothe neck pain, lie down on your back and place the Hot & Cold Neck Cradle under your neck. Adjust placement until your head and neck are resting comfortably on the cradle.

COOLING THERAPY

Place the Hot & Cold Neck Cradle in the freezer for 4-6 hours to freeze gel; product will stay cold 15-20 minutes. Refreeze as necessary.

HEAT THERAPY

Place the Hot & Cold Neck Cradle in 120-degree (very hot) tap water for 7-10 minutes; test before applying to skin. Roller stays heated approximately 10-15 minutes. Reheat as necessary.

DO NOT USE OTHER HEATING METHODS—conventional or microwave ovens or stove tops, for example—as overheating may cause burns or damage the Hot & Cold Neck Cradle. We do not recommend placing the Hot & Cold Neck Cradle in boiling water, as this may damage the seals on the product and/or cause burns.

GAIAM WELLBEING™

HOT & COLD NECK CRADLE

Setup, Workout Tips and Care & Safety Guide

Please read before setting up or using your Hot & Cold Neck Cradle.

CARE

1. Dry the Hot & Cold Neck Cradle with a towel after use. Check for nicks or tears.
2. Spot clean Hot & Cold Neck Cradle with a damp cloth. Air dry.
3. Avoid exposing the Hot & Cold Neck Cradle to sharp or abrasive surfaces, or to direct heat or excessive sunlight.

IMPORTANT SAFETY CAUTIONS!

1. Not all exercise is suitable for everyone. To reduce the risk of injury, consult your doctor before beginning this or any exercise program or before using this or any other exercise equipment. The instruction presented herein is in no way intended as a substitute for medical counseling.
2. If you suffer from heart disease, high blood pressure or any other disease or condition, consult your doctor before beginning this workout.
3. Perform exercises in a slow and controlled manner.
4. Stop and rest if you feel dizzy or short of breath.
5. This product is sold without any warranties or guarantees of any kind and manufacturer and distributor disclaim any liability, loss or damage caused by its use.

CAUTION: Use only as directed. Consult with your doctor or physical therapist before beginning any self-treatment program. Suitable for ages 13 and over. Keep out of reach of younger children. Do not use if you have diabetes, or if you have open sores or cracked or blistered skin in the treatment area. If symptoms worsen, decrease pressure; if adverse symptoms persist, discontinue use and contact your doctor.

Learn more at <http://life.gaiam.com>

©2016 Gaiam and the "Flower of Life" logo design are registered trademarks and Wellbeing is a trademark of Gaiam Americas, Inc. or its subsidiaries. Manufactured and distributed by Gaiam Americas, Inc., Louisville, CO 80027-2452.

MADE IN CHINA.

SKU 05-5210/CARE