## EquiPlus

**Flüssiges Konservierungsmedium für Hengstfrischsamen**  
**Liquid Culture Medium for Chilled Stallion Semen**  
**Medio de Conservación Líquido para Semen Equino**  
**Жидкая культуральная среда для охлажденной спермы**  
**жеребцов**

**DE**

**EN**

**ES**

**RU**

**13570/0301**

31357/0311

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**EquiPlus** is a culture medium in which stallion semen can be preserved at +4°C to +6°C for approx. 3 days. **EquiPlus** is also recommended for the centrifugation of stallion semen. **EquiPlus** liquid is delivered in a 100 ml glass bottle and is whitish; no further components must be added.

## 1. Composition

**EquiPlus** consists of the following components:

- Caseinates
- Sugars
- Buffer
- Antibiotic

## 2. Package size

**EquiPlus** liquid, 100 ml, with Gentamicin sulphate  
13570/0301

## 3. Application

### 3.1. Preparation

The required amount of **EquiPlus** must be warmed to +30°C to +35°C.

### 3.2. Water quality

**EquiPlus** liquid is manufactured with ultra-pure water. Addition of water is not necessary.

### 3.3. Recommended application

#### 3.3.1. Preparation of chilled semen

The freshly collected stallion semen is maintained at a temperature of +30°C to +33°C and analysed. After removing the gel, its volume, concentration and motility are determined. Depending on the desired number of progressive motile sperm cells in the insemination dose, the required amount of extender is calculated. The extender and the semen must have the same temperature at the moment of the first dilution (+/-1°C). The extender is slowly added in several steps to the filtered, gel-free ejaculate and the 2 fractions are gently but thoroughly mixed.

#### Calculating the final volume of extended semen:

##### Dilution without centrifugation

The optimum dilution rate should be at least 1:3 (semen:extender), or 25 million motile semen cells/ml after final dilution. The recommended motile sperm per insemination dose amounts 600 million. The inseminated volume should be 10-30 ml.

Example:

Data of the freshly collected ejaculate:

Volume of ejaculate = 45 ml

Concentration of Semen = 145 million/ml

Total of sperm cells in ejaculate = 6,525 billion ( $6,525 \times 10^9$ )

Motility = 60%

→ motile sperm in the ejaculate = 3,915 billion

From this ejaculate 6 insemination doses can be obtained, each containing 600 million sperm (3,915 billion/600 million = 6.5 doses).

To reach an optimum dilution with 25 million motile sperm/ml, the final volume must be 156 ml (3,915 billion/25 million = 156.6 ml). Each semen dose will amount 26 ml (156 ml/6 doses = 26 ml).

#### Dilution with centrifugation:

The optimum dilution rate and the number of motile sperm per dose remain equal. The volume of the doses may be reduced (10-15 ml).

Example:

Data of the freshly collected ejaculate:

Volume of ejaculate = 45 ml

Concentration of semen = 145 million/ml

Total of sperm in ejaculate = 6,525 billion sperm

Motility = 60%

→ motile sperm in the ejaculate = 3,915 billion of sperm

After centrifugation and taking off the supernatant (leaving approx. 1 ml seminal plasma with the centrifuged pellet), 3-5 ml of extender are added to re-suspend the pellet by careful movement.

Then volume is determined, motility evaluated and based on counting with the counting chamber, the new number of sperm cell of the ejaculate is calculated.

After the dilution, the semen is no longer maintained in warm conditions. The recommended cooling rate is 0.05° to 0.1°C per minute, respectively within 2.5 to 5 hours from ambient temperature to +5°C.

EquiPlus is recommended for semen preservation for up to 3 days at approx. +5°C. A daily quality analysis until the application is highly recommended.

#### **3.3.2. Centrifugation**

The prepared extender and the gel-free ejaculate are mixed 1+1. Centrifugation is performed at 500-800 g for 5-10 minutes. After centrifugation the supernatant is removed and the pellet resuspended with extender for chilled or freezing semen. To achieve a higher yield of sperm cells it is possible to work with CushionFluid (Ref. 13580/0001).

## 4. Storage, shelf life and important hints

### 4.1. Storage, shelf life

The bottle with **EquiPlus** liquid must be stored in cool (below +10°C), dry and dark conditions. The minimum shelf life is printed on each bottle.

### 4.2. Important hints

**Keep the product out of the reach of children.  
The product is not suitable for treatment or consumption by human beings or animals.**

Before each use of the extender, make sure to analyse precisely, whether the extender shows the described properties, especially regarding its consistency and aspect. If the extender shows a different consistency or aspect as normal, it must not be used or Minitüb must be contacted first.

### 4.3. Trouble shooting

A couple of factors are crucial for good semen conservation and insemination success, among others:

- Use of exclusively sperm friendly materials and surfaces for semen and extender (collection gloves, semen and extender bottles, glassware free of mineral residues)
- Good temperature management, particularly avoiding temperatures below +30°C for undiluted semen and below +4°C for diluted semen, as well as temperature fluctuations
- Avoiding bacterial contamination (appropriate stallion management and preparations of collection; semen collection in clean ambient, appropriate hygiene protocol in the lab)
- Avoiding sperm toxic factors (water, UV light, dust, detergent)
- Quality of the purified water (see point 3.2)

### 4.4. Consulting and training

You can rely on the Minitube scientific experts for advice and technical support concerning all aspects of your semen production and/or artificial insemination process. Workshops and seminars are offered regularly and also on an individual basis. Please refer to our web site for information regarding our current program.

## 5. Production standard and quality control

### 5.1. Raw material quality

All raw materials used in Minitube media are produced according to GMP and DIN ISO 9001:2015 norms, and are certified for meeting Ph Eur, BP or USP standards. They are tested and comply with the current international guidelines regulating the quality of pharmaceutical substances.

### 5.2. Production standard

The mixing and bottling is performed under GMP conditions in a controlled atmosphere room and documented with weighing records. Under the supervision of a team of veterinary

specialists, the final product is subjected to a series of chemical and physical tests designed to ensure that Minitube media can consistently provide semen conservation with high spermatological efficacy over a wide range of climatic and process conditions.

### **5.3. Certificates**

The general quality certificate and the batch certificate are available upon request.

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## Literatur/ Literature/ Literatura/ Литература

Aurich, C. (2005): **Factors affecting the plasma membrane function of cooled-stored stallion spermatozoa.** Animal Reproduction Science, Vol. 89, Issues 1-4, 65-75.

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Pagl, R., Aurich, J.E., Müller-Schlösser, F., Kankofer, M., Aurich, C. (2006): **Comparison of an extender containing defined milk protein fractions with a skim milk-based extender for storage of equine semen at 5°C.** Theriogenology, Vol. 66, Issues 1115-1122.

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All rights reserved, including the translation.  
Reprint, even in extracts only with a prior written permission. No part of this manual may be reproduced or transmitted in any form or by any means electronic, mechanical, photocopying, without the prior written permission of Minitüb GmbH.  
All brand product names are trade marks or registered trade marks of the respective title holder.  
Errors and technical alterations excepted.  
Use product only in perfect condition and in compliance with the manual.  
Keep safe for future consulting (according to EN ISO 12100).  
Translated from German.

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