INTENDED USE

The Vcheck Feline NT-proBNP is an *in vitro* diagnostic test kit for the quantitative measurement of NT-proBNP concentration in feline serum. N-terminal pro-B type natriuretic peptide (NT-proBNP) is cleaved from BNP which is produced by the muscle cells of the heart and increases with excessive stretching of the cells. NT-proBNP concentration reflects the degree of cardiac activation secondary to stimulus, such as stretching, allowing this marker to be used to assess the magnitude of cardiac muscle stretching. NT-proBNP is a valuable biomarker for differentiating cardiac and respiratory causes of dyspnea and can be used for screening occult heart disease in asymptomatic cats. The BIONOTE Vcheck Feline NT-proBNP is designed to be used only by veterinarians.

PRINCIPLE

The Vcheck Feline NT-proBNP Test kit is a fluorescent immunoassay for the quantitative measurement of feline NT-proBNP concentration. The Vcheck Feline NT-proBNP Test kit uses specific anti-feline NT-proBNP antibodies that bind to feline NT-proBNP. The test procedure involves dissolving fluorescence microparticles conjugated to these specific antibodies into assay diluent that react with feline NT-proBNP in the sample.

MATERIALS PROVIDED

Reagent	5 Tests/ Kit	10 Tests/ Kit	20 Tests/ Kit
① Vcheck Feline NT-proBNP Test device	5	10	20
② Assay diluent tube	5	10	20
③ Disposable pipette tip	5	10	20
④ Instructions for use	1	1	1

MATERIALS REQUIRED, BUT NOT PROVIDED

- BIONOTE Vcheck Analyzer
- 100 µl pipette

STORAGE AND STABILITY

- Store the test kit at 1~30°C. **DO NOT FREEZE.**
- Do not store the test kit in the direct sunlight.
- The test kit is stable until the expiry date that is marked on the package label.

PRECAUTIONS

- The test kit is for feline use only. Do not use for other animals.
- The test device is sensitive to humidity and heat. Perform the test immediately after removing the test device from the foil pouch.
- Do not reuse the test components.
- Do not touch the membrane in the result window of the test device.
- Do not use the test kit beyond the stated expiry date marked on the label.
- Do not use the test kit if the pouch is damaged or the seal is broken.
- Do not mix components from different lot numbers; the components in this kit have been quality control tested as a standard batch unit.
- All samples should be handled as being potentially infectious. Wear protective gloves while handling samples. Wash hands thoroughly afterwards.
- Decontaminate and dispose of all samples, used kits, and potentially contaminated materials safely in accordance with national and local regulations.

COLLECTION AND PREPARATION OF SAMPLE

- Feline serum should be used with this test.
- [Serum]**
 - Collect the whole blood using a plain tube **ONLY**. Leave to settle for 30 minutes for blood coagulation and then centrifuge to obtain a serum supernatant.
 - The samples should be tested immediately after collection.** If serum samples are not tested immediately, serum can be frozen (-20°C or colder). Frozen samples should be brought to room temperature (15~30°C) prior to use.

TEST PROCEDURE

All reagents and frozen sample must be at room temperature(15~30°C) before use.

[Coding]

- Turn on V200 Analyzer and select "Standard Test".

- Remove the test device from the foil pouch. Once the "Insert Device" is displayed in the screen, insert the test device.


- After checking Feline NT-proBNP item name and test procedure on the display window, proceed as follows.

[Dilution of sample & Measurement]

- Remove the sealer of assay diluent tube, draw 100 µl sample using 100 µl pipette and add the sample into an assay diluent tube. And then, Mix the sample with diluent by pipetting for 5~6 times.

- Add the mixed sample(100 µl) into the sample hole of the test device. Then, press the "Start" button.

- The V200 analyzer will display the test result on the screen after 10 minutes.
- Remove the test device.

INTERPRETATION OF THE RESULT

- Read the concentration value of feline NT-proBNP appearing on the display of the BIONOTE Vcheck Analyzer. (50~1500 pmol/L)
- If "< 50 pmol/L" appears on the display, it means the concentration of feline NT-proBNP in the specimen is less than 50 pmol/L.
- If "> 1500 pmol/L" appears on the display, it means the concentration of feline NT-proBNP in the specimen is greater than 1500 pmol/L.

REFERENCE RANGE

< 100 pmol/L	≥ 100 pmol/L
Normal	Abnormal Additional diagnostics are recommended

* An abnormal NT-proBNP test result should always be interpreted in combination and other diagnostic findings.

** In cats with respiratory signs, if the NT-proBNP is ≥ 270 pmol/L, congestive heart failure is the most likely cause of the clinical signs.

SCREEN MESSAGES AND TROUBLE SHOOTING

[V200]

Error message	Error description
Contaminated Device	The test device is damaged or inserted improperly. Solution: Discard the test device and retest with a new test device and a new specimen.
Insufficient Sample	An insufficient amount of blood has been applied. Solution: Retest with a new test device with enough specimen, ensuring that blood is placed in to the narrow channel in the top edge of the test device.
Expired Device	The test devices are expired. Solution: Retest with a new test device that is not expired.
Temperature Error	The environmental temperature is above or below the operating range of the analyzer. Solution: Move to an area in the acceptable temperature range for the analyzer and perform the test. Do not heat or cool the analyzer artificially.
Printer Connection Fail	The communication between analyzer and barcode or printer has failed. Solution: Reconnect the analyzer and external device. If the error continues after turning ON/OFF the analyzer, please contact BioNote, Inc.
Barcode Error	The measured total hemoglobin is out of the range of 7 to 23 g/dL. Solution: This error occurs when a specimen has a total hemoglobin in the abnormal range. If the error continues after turning ON/OFF the analyzer, please contact BioNote, Inc.
Extremely High Total Hemoglobin	The test is invalid. Solution: Retest with a new test device and a new patient specimen. If the error continues after turning ON/OFF the analyzer, please contact BioNote, Inc.
Result: Invalid	The calibration is overdue. Solution: If the error continues after turning ON/OFF the analyzer, please contact BioNote, Inc.
Calibration Overdue	A test device that is not supported by the analyzer has been loaded. Solution: Check whether the test device is manufactured by BioNote, Inc.
Not Supported Device	Internal error has occurred. Solution: If the error continues after turning ON/OFF the analyzer, please contact BioNote, Inc.
EEE	

PERFORMANCE CHARACTERISTIC

- Measuring Range
50~1500 pmol/L of feline NT-proBNP is measurable in serum. For > 1500 pmol/L results, samples may be diluted with negative serum and measured again for a numerical result. To calculate the final feline NT-proBNP concentration, the measured result must be multiplied by the dilution factor.
- Interfering substances study
No interference was observed for each substance up to the concentration presented in the following table.

Interfering substances	Concentration
Hemoglobin	≤ 37.5 mg/dL
Intralipid	≤ 2,500 mg/dL
Cholesterol	≤ 800 mg/dL
Bilirubin (total)	≤ 15 mg/dL
Vitamin C	≤ 1,000 mg/L

Severely hemolyzed samples can have an effect on the results due to interferences of hemoglobin.

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