## INTENDED USE

Vcheck cTSH is an *in vitro* diagnostic test kit for the quantitative measurement of canine TSH (Thyroid Stimulating Hormone) in serum. TSH is secreted by the anterior pituitary and promotes release of thyroid hormones. Measurement of TSH provides information on the interaction between the pituitary and thyroid gland. Serum cTSH concentration is usually increased in dogs with hypothyroidism. Vcheck cTSH can assist in diagnosing canine hypothyroidism. It is designed to be used only by veterinarians.

## PRINCIPLE

The Vcheck canine TSH test kit is based on a sandwich immunoassay method for the quantitative measurement of cTSH concentration. When the sample is delivered to the sample hole of test device, cTSH in the specimen reacts with the colloidal gold-labeled anti-TSH monoclonal antibody in the conjugate pad. The mixture then migrates along the nitrocellulose membrane and forms complexes with the anti-TSH monoclonal antibody coated on the membrane. The density of the test line reflects the canine TSH concentration in the canine serum. The BIONOTE Vcheck Analyzer reads the density of the test line and calculates the TSH concentration from the calibration curve data. The control line is a reference line that indicates the test has been performed correctly.

## MATERIALS PROVIDED

Reagent	5 Tests/Kit
① Vcheck cTSH Test device	5
② Assay diluent tube	5
③ Pipette tip	5
④ Instructions for use	1

## MATERIALS REQUIRED, BUT NOT PROVIDED

- BIONOTE Vcheck Analyzer
- 100 µl pipette

## STORAGE AND STABILITY

- Store the test kit at 2~8°C. **DO NOT FREEZE.**
- Do not store the test kit in the direct sunlight.
- The test kit is stable until the expiry date that is marked on the package label.

## PRECAUTIONS

- The test kit is for canine use only. Do not use for other animals.
- The test device is sensitive to humidity and heat. Perform the test immediately after removing the test device from the foil pouch.
- Do not reuse test components.
- Do not touch the membrane in the result window of the test device.
- Do not use the test kit beyond the stated expiry date marked on the label.
- Do not use the test kit if the pouch is damaged or the seal is broken.
- Do not mix components from different lot numbers, the components in this kit have been quality control tested as a standard batch unit.
- All samples should be handled as being potentially infectious. Wear protective gloves while handling samples. Wash hands thoroughly afterwards.
- Decontaminate and dispose of all samples, used kits and potentially contaminated materials safely in accordance with national and local regulations.
- Hemolyzed samples or serum with precipitate such as fibrin may give erroneous results

## COLLECTION AND PREPARATION OF SAMPLE

- Canine serum should be used with this test. A method for preparing the serum is as follows.
- [Serum]** Collect the whole blood into a blood collection tube containing NO anticoagulant. Leave to settle for 30 minutes for blood coagulation and then centrifuge to obtain the serum supernatant.
- If serum samples are not tested immediately, they should be refrigerated at 2~8°C and used within 7 days. For longer storage, serum can be frozen (-20°C or colder). Frozen samples should be brought to room temperature (15~30°C) prior to use.

## TEST PROCEDURE

Allow all kit components and sample to reach room temperature (15~30°C) prior to testing.

### [Coding]

- Turn on V200 Analyzer and select "Standard Test".


- Remove the test device from the foil pouch. Once the "Insert Device" is displayed in the screen, insert the test device.

### [Dilution of sample & Measurement]

- Using a 100 µl pipette, draw 100µl of the sample and add to the assay diluent tube.

- Mix the sample with diluent by pipetting for 5~6 times.

- Add 100 µl of the sample mixture into the sample hole of the test device using a 100 µl pipette and press **[START]** to initiate testing.

- The V200 analyzer will display the test result on the screen after 15 minutes.

- Remove the test device.

## INTERPRETATION OF THE RESULT

- Read the concentration value of canine TSH appearing on the display of the BIONOTE Vcheck Analyzer. (0.25~5.00 ng/ml)
- If "<0.25 ng/ml" appears on the display, it means the concentration of canine TSH in the specimen is less than 0.25 ng/ml.
- If ">5.00 ng/ml" appears on the display, it means the concentration of canine TSH in the specimen is greater than 5.00 ng/ml.

## REFERENCE RANGE

Normal	High
< 0.5 ng/ml	≥ 0.5 ng/ml

\* It is recommended that each laboratory establishes its own reference range.

## SCREEN MESSAGES AND TROUBLE SHOOTING

[V200]

Error message	Error description
Contaminated Device	The test device is damaged or inserted improperly. <b>Solution:</b> Discard the test device and retest with a new test device and a new specimen.
Insufficient Sample	An insufficient amount of blood has been applied. <b>Solution:</b> Retest with a new test device with enough specimen, ensuring that blood is placed in to the narrow channel in the top edge of the test device.
Expired Device	The test devices are expired. <b>Solution:</b> Retest with a new test device that is not expired.
Temperature Error	The environmental temperature is above or below the operating range of the analyzer. <b>Solution:</b> Move to an area in the acceptable temperature range for the analyzer and perform the test. Do not heat or cool the analyzer artificially.
Printer Connection Fail	The communication between analyzer and barcode or printer has failed. <b>Solution:</b> Reconnect the analyzer and external device. If the error continues after turning ON/OFF the analyzer, please contact BioNote, Inc.
Barcode Error	The measured total hemoglobin is out of the range of 7 to 23 g/dL. <b>Solution:</b> This error occurs when a specimen has a total hemoglobin in the abnormal range. If the error continues after turning ON/OFF the analyzer, please contact BioNote, Inc.
Extremely Total hemoglobin	The test is invalid. <b>Solution:</b> Retest with a new test device and a new patient specimen. If the error continues after turning ON/OFF the analyzer, please contact BioNote, Inc.
Result: Invalid	The calibration is overdue. <b>Solution:</b> If the error continues after turning ON/OFF the analyzer, please contact BioNote, Inc.
Calibration Overdue	Loading a test device that is not supported by the analyzer. <b>Solution:</b> Check whether the test device is manufactured by BioNote, Inc.
Not Supported Device	Internal error has occurred. <b>Solution:</b> If the error continues after turning ON/OFF the analyzer, please contact BioNote, Inc.
EEE	

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