

User Manual

MELAtherm[®] 10

Washer-Disinfector

From software version 1.311



EN

Dear doctor,

We thank you for your confidence demonstrated by the purchase of this MELAG product. As an owner-run and operated family concern founded in 1951, we have a long history of successful specialization in hygiene products for practice-based use. Our focus on innovation, quality and the highest standards of operational reliability has established MELAG as the world's leading manufacturer in the instrument treatment and hygiene field.

You, our customer are justified in your demand for the best products, quality and reliability. Providing "**competence in hygiene**" and "**Quality – made in Germany**", we guarantee that these demands will be met. Our certified quality management system is subject to close monitoring: one instrument to this end is our annual multi-day audit conducted in accordance with ISO 13485 and ISO 9001. This guarantees that all MELAG products are manufactured and tested in accordance with strict quality criteria.

The MELAG management and team.

CE 0197

Contents

1 General Guidelines	5
User group and validity	5
Symbols used.....	5
Formatting rules	5
Symbols on the device	6
Disposal	6
2 Safety	7
3 Performance specifications	9
Intended use	9
User benefits.....	9
Program sequence.....	10
Approved process agents	11
4 Description of the device	13
Scope of delivery.....	13
Views of the device	13
Operating panel and acoustic signals	16
Menu structure	17
Water softening unit	18
5 First steps	19
Setup and Installation.....	19
Switching the device on and off	19
Opening and closing the door	19
Emergency release of the door	20
Inserting the basis basket	21
Filling the regeneration salt	21
Regenerating the water softening unit	22
Metering process agents.....	23
Holding process agents ready.....	23
Bleeding the metering hoses.....	26
6 Cleaning and disinfection	27
The nature of the load	27
Arranging the load.....	27
Decontaminating hollow-bodied instruments	28
Decontaminating dental transfer instruments.....	29
Decontaminating ophthalmological instruments.....	29
Overview of programs	31
Selecting, starting and following the program	32
Manual program abort.....	33
Removing the load after program end.....	34
7 Logging	35
Batch documentation	35
Output media.....	35
Outputting logs immediately and automatically	37

Subsequent log output	38
Deleting the saved logs	38
Determining the format for the program logs	38
Finding the logs	41
8 Settings	42
SETUP menu	42
Setting the water supply	42
Setting automatic logging	42
Setting date and time	43
Setting the display contrast	44
Selecting the language	44
Setting the water hardness	45
9 Function tests	46
Automatic and manual functional check	46
Measuring conductivity	46
10 Maintenance	47
Regular checks and cleaning	47
Cleaning on demand	49
Avoiding staining	49
Changing the filter in the drying fan	50
Maintenance	51
(Process) validation	51
11 Pause times	52
Decommissioning	52
Transport within the practice	53
12 Malfunctions	54
General events	54
Warnings	56
Malfunction messages	61
13 Technical data	71
14 Accessories and Spare Parts	72
Accessories	72
15 Documentation and approval	73
Glossary	74

1 General Guidelines

Please read this user manual carefully before commissioning the device. The user manual includes important safety information. The functionality and value-retention of this device depends on the care accorded to it. Please store these user manual carefully and in close proximity to your device. It represents a component of the product.




Should the user manual no longer be legible, damaged or lost, please obtain a new copy from MELAG. State the device type and your address in an e-mail.

The device type is specified on the type plate on the rear of the device.

User group and validity

This manual applies for the devices MELAtherm 10 DTA and MELAtherm 10 DTB. This manual is addressed to doctors, their assistants and service departments.

Symbols used

Symbol	Explanation
	Indicates a dangerous situation, which if not avoided, could entail slight to life-threatening injuries.
	Draws your attention to a situation, which if not avoided, could result in damage to the instruments, the practice fittings or the device.
	Draws your attention to important information.

Formatting rules

Example	Explanation
see Chapter 2	Reference to another text section within this document
Settings	Words or phrases appearing on the display of the device are marked as display text.

Symbols on the device



Manufacturer of the medical device



Date of manufacture of the medical device



Medical device serial number from the manufacturer



Article number of the medical device



Refers to the lowest and highest water temperature to which the device can be safely subjected.



Flow pressure on the water inflow connected from min. to max.



Internal device fuse, rated in amperes [A]



The user manual includes important safety information. Failure to comply with these instructions can result in injury and material damage.



Please read this user manual carefully before commissioning the device.



In affixing this CE mark, the manufacturer declares that this medical product fulfils the basic requirements of the Medical products directive. The four-digit number confirms that this is monitored by an approved certification agency.



The device may not be disposed as domestic waste. The vendor is responsible for appropriate disposal of the device - it must be delivered to the vendor to be disposed of. By the designation of an apparatus with this symbol, the manufacturer furthermore declares that he satisfies all requirements of the law concerning the release, redemption and environmentally sound disposal of electric and electronic appliances.

Disposal

MELAG devices are synonymous for long-term quality. When you eventually need to decommission your MELAG device, we offer a special disposal service. Simply contact your stockist.

The packaging protects the device against transport damage. The packaging materials have been selected for their environmentally-friendly and recycling properties and can be recycled. Returning the packaging to the material flow reduces the amount of waste and saves raw materials. Dispose of all non-required packaging materials at the collection points of the dual system.

Dispose of waste from process agents in accordance with the specifications from the manufacturer of the process agents. Information regarding this topic is provided by the safety data sheets or can be obtained directly from the manufacturer of the process agents.

Dispose of accessories and consumption media which you no longer require (e.g. used filters) in the appropriate manner. Comply with all relevant disposal specification in terms of possibly contaminated waste.

2 Safety



When operating the device, comply with the following safety instructions as well as those contained in subsequent chapters. Use the device only for the purpose specified in these instructions. Failure to comply with the specifications of the safety instructions can result in injury and/or damage to the device.

Set-up, installation and commissioning

- Check the device for any damage suffered during transport after unpacking.
- The device should only be set-up, installed and commissioned by MELAG authorized persons.
- The connections for electrical provision and water supply and effluent must be set-up by trained personnel.
- Use of the optional leak detector (water stop) minimizes the risk of water damage.
- In accordance with current VDE specifications, the device is unsuitable for operation in explosive atmospheres.
- Install and operate the device in a frost-free environment.
- The device is conceived for use outside the patient area. The device should be located a minimum of 1.5 m radius away from the treatment area.
- The documentation media (computer, CF card reader etc.) must be placed in such a way that they cannot come into contact with liquids.
- Observe all the information contained in the technical manual during commissioning.

Power cable and power plug

- Comply with all legal requirements and locally-specified connection conditions.
- Never operate the device if the plug or power cable are damaged.
- The power cable and plug should only be replaced by authorized personnel.
- Never damage or alter the power plug or cable.
- Never unplug by pulling on the power cable. Always take a grip on the plug.
- Ensure that the power cable does not become jammed in.
- Never lead the cable along a source of heat.
- Never fix the power cable with sharp objects.

Daily operation

- Use only those instruments designed by their manufacturer for automatic treatment in a washer-disinfector. Comply with the instructions issued by the instrument manufacturer in accordance with DIN EN ISO 17664. It is especially important to observe the manufacturer's information regarding cleaning instruments for the first time whilst purchasing new instruments.
- Use only original MELAG accessories or those from other suppliers authorized for use by MELAG.
- When using non-MELAG accessories for the mounting of instruments (especially hollow-bodied instruments) comply with the manufacturer's information.
- Comply with the specifications of the national standards pertaining to the decontamination of instruments and directives, the manufacturer's decontamination instructions and those from the AKI.
- Instrument decontamination may only be performed by trained personnel.
- The fore ventilation slits may not be covered.
- Only ever operate the device using with the basis basket provided for this purpose.
- Never operate the device unattended. Unsupervised operation of the device can result in damage to the device or your facility. In such a case, MELAG does not accept any liability.

Process agents

- Handle all process agents with care. The cleaning, neutralization and rinsing aids contain irritants and even caustic substances.
- Use only those process agents approved by MELAG for this device. Observe the operating and safety information from the process agent manufacturer. If, despite observation of the manufacturer's information, the process agents have a negative effect on the material of the instruments or the device, liability lies with the manufacturer of the process agents.
- The use of process agents not approved by MELAG absolves MELAG of all liability whatsoever for any damage to the device or the instruments.
- Should you have any questions concerning the compatibility of the process agents with the instruments, please consult the manufacturer. MELAG provides information for the application of the process agents in the device but does not take any responsibility for their effects on the instruments.
- Any fluids in the drawer and the floor tank underneath can also contain process agents in case of damage. Ensure that you observe the information of the respective process agent manufacturer.

Transport and storage

- Store and transport the device in a frost-free environment.
- Avoid strong shocks/vibrations.

Maintenance

- Have the maintenance done only by authorized persons.
- Maintain the specified servicing intervals.

Malfunctions

- Should the device issue the same malfunction message repeatedly, turn off the device and if necessary, inform your stockist.
- Only have the device repaired by authorized persons.

3 Performance specifications

Intended use

In accordance with DIN EN ISO 15883-1 and 2, this device is a washer-disinfector intended for use in a medical context such as clinics and medical and dental practices. You can subject thermostable medical instruments (i.e. instruments which are heat resistant to a temperature of 95 °C) to automatic decontamination as long as they are suitable for this purpose and have been approved for such treatment by their manufacturer. The cleaning is undertaken via the use of water and a chemical cleaning agent. Subsequent disinfection is thermal disinfection. This device is not intended to be used in a patient environment.

This device is NOT suitable for the decontamination of:

- Thermo-unstable instruments e.g. flexible endoscopes
- Laboratory waste requiring disposal
- Crockery
- Bedpans

User benefits

Universal use

The device both cleans and disinfects. The disinfection phase is conceived so as to reach an A0 value of at least 3000 is achieved. This kills vegetative bacteria, fungi and their spores and viruses (inc. HBV, HCV) so that the effective range AB is reached in accordance with the specifications of the Robert-Koch-Institut.

Active drying

The device is equipped with active drying. An integrated drying fan dries the instruments from outside and in after cleaning and disinfection. The additionally fitted HEPA filter guarantees drying with contamination-free air. The automatic decontamination of hollow-bodied instruments is also possible. This protects the instruments from stain accretion and rusting. The geometry of some hollow-bodied instruments mean that they require additional drying.

Automatic filter recognition

The device recognizes before a program start whether the fine sieve has been inserted in the base of the washing chamber. The fine sieve avoids a situation in which instrument components enter the opening of the drain pump or the circulation pump during cleaning, thereby compromising the function of the pumps, rinse arms and the injector rail.

Internal water softening

The device is equipped with an internal water softening unit. The water hardness of the local drinking water is set on the device. The internal water softening unit then automatically adjusts itself to the most suitable performance. This ensures best decontamination results.

Monitoring the speed of the rinse arm

The rotation speed of the rinse arms is subject to permanent monitoring during a program run. This ensures that the cleaning process proceeds without hindrance and the rinse arms do not become blocked e.g. by protruding instruments in the washing chamber.

Monitoring cleaning pressure

The rinse pressure is monitored by a pressure sensor during the program run. This ensures an effective cleaning performance. The device aborts a current program if too much foam is generated.

Metering monitoring

The required amounts of cleaning agent and neutralizer are measured out using a hose metering pump. A measurement turbine performs flow monitoring. The rinse aid is metered using a hose pump subject to monitoring for rotation speed.

Drawer for process agents

The drawer for the process agents is located in the lower area of the device in which the cleaning agent, neutralizer and rinse aid containers are stored.

Automatic conductivity measurement

If the MELAtherm 10 is provided with DI water in the final rinse, the DI feed water is subject to automatic, internal conductivity measurement.

Door emergency release

The door can be opened manually via the emergency release following a power outage or malfunction.

Program sequence

The following program steps are indicated on the display during the program run.

Pre-cleaning

Pre-cleaning is performed with cold (min. 22 °C), softened or de-ionized water without process agents. Instrument soiling and contamination (proteins and/or coarse organic deposits) will be mechanically removed to avoid denaturation from excessively-high water temperatures.

Cleaning

The pre-cleaning is followed by the actual cleaning phase with warm process water and the addition of a mildly-alkaline cleaning agent. As part of the cleaning process, a defined volume of process water is let into the washing chamber and heated to 40 °C. The mildly-alkaline cleaning agent is automatically metered to the process water as soon as the set temperature has been reached. The temperature is maintained (holding time). This is followed by the actual cleaning.

The process water is heated further for the actual cleaning process. Any soiling or organic deposits still on the instruments are removed at a temperature of 55 °C.

Neutralization

The main wash is followed immediately by neutralization. Neutralization reduces the alkalinity introduced during cleaning and removes acid-soluble deposits such as limescale and extraneous rust etc.

A metered amount of neutralizer is added during this phase. The neutralized process water is subject to continual circulation for a defined period. Finally, the process water is pumped out of the washing chamber completely.

Intermediate rinsing

Intermediate rinsing is the preparatory step for thermal disinfection in which the instruments are rinsed without process agents. The rinsing serve to reduce the residual process agents to a defensible level. In this phase, a metered amount of cold process water is let into the washing chamber and subject to continual circulation for a defined period. Finally, the process water is pumped out of the washing chamber.

Disinfection

The actual disinfection is performed after intermediate rinsing. The disinfectant effect of MELAtherm 10 is achieved using thermal disinfection and subsequent drying.

A rinse aid added to the process water at a temperature of 70 °C. The actual thermal disinfection is performed via continual circulation and heating of the process water to 90 °C and a holding time at this temperature of min. 5 minutes. Following the disinfection, the process water is pumped out of the washing chamber completely.

Drying

Drying is effected by drawing ambient air through a class H 13 HEPA filter. The instruments are dried inside and out with hot, filtered air. This prevents rust accretion on the instruments.

Displaying the batch counter

The display shows the last batch number run and the total batch counter after a program end or the end of a program abort.

Approved process agents



NOTICE

Observe the operating and safety information from the process agent manufacturer. Process agents from different manufacturers may not be mixed. Any change to another authorized combination may only be performed by trained service partners. The metering concentration must be adapted to local conditions. This is to be performed by the service technician during setup in accordance with the manufacturer's information.



NOTICE

Observe the operating and safety information from the process agent manufacturer. If, despite observation of the manufacturer's information, the process agents have a negative effect on the material of the instruments or the device, liability lies with the manufacturer of the process agents. The use of process agents not approved by MELAG absolves MELAG of all liability whatsoever for any damage to the device or the instruments.

The following combinations of process agents were tested for suitability for use with this device, brought into circulation on 1.7.2012:

Manufacturer	Cleaner	Neutralizer	Rinse aid
MELAG	MEtherm 50 - Mildly-alkaline cleaning agent	MEtherm 55 - Neutralizer C MEtherm 56 - Neutralizer P	MEtherm 60 - Rinse aid
Alpro*)	TR-3	TR-Neutralizer	TR-Clear
B. Braun	Helimatic cleaner alkaline	Helimatic neutralizer C	Helimatic rinse neutral
Bode	Dismoclean 21 clean	Dismoclean 25 acid	Dismoclean 64 neutra-dry
Borer*)	deconex 28 ALKA ONE-x deconex 22 LIQ-x	deconex 25 ORGANACID deconex 26 MINERALACID	deconex 64 NEUTRADRY
Dr. Schumacher	Thermoton Cleaner	Thermoton N	Thermoton clear
Dr. Weigert	neodisher MediClean forte neodisher MediClean Dental	neodisher N neodisher Z neodisher N dental neodisher Z dental	neodisher MediKlar neodisher MediKlar dental neodisher MediKlar special
Schülke*	Thermodont alka clean Thermosept RKF forte Thermosept alka clean forte	thermodent neutralizer thermosept NKP thermosept NKZ	Thermodont clear Thermosept BSK
Henry Schein	Eurosept Thermo Cleanser	Eurosept Thermo Neutralizer	Eurosept Thermo Rinse

*) When using process agents of manufacturers Alpro or Borer, it is absolutely essential and regardless of the maintenance message after twelve month to replace all hoses of the dosing unit during maintenance. Ensure that the process media selected are suitable for your instruments. Please address any questions to the manufacturer of the process agents. Please observe and comply with the notes specific to the process agents contained in the maintenance record.

Preset metering concentrations¹⁾**PLEASE NOTE**

The following setting values are effective from software version 1.311.

Program	Cleaner	Neutralizer	Rinse aid
Universal-Program	6 ml/l	1.5 ml/l	0.3 ml/l
Quick-Program	6 ml/l	1.5 ml/l	0.3 ml/l
Intensive-Program	10 ml/l	1.5 ml/l	0.3 ml/l
Ophthalmic-Program	6 ml/l	1.5 ml/l	-----

¹⁾ For using drinking water (hardness range medium - hard) we recommend the dosage concentrations specified above. The metering concentration must be adjusted individually depending on the local conditions. The pre-set metering quantities can be used for various combinations of process agents. Observe the operating and safety information of the process agent manufacturer. This information can be read-off from the canister in the units "ml/l" or on separate data sheets provided by the manufacturer. The process-relevant parameters must be adapted depending on the nature of the soiling, the quality of the tap water and other framework conditions. Please contact the manufacturer of the process agents regarding this question.

4 Description of the device

Scope of delivery

Please check the scope of delivery before setting up and connecting the device.


Standard scope of delivery

- MELAtherm 10 washer-disinfector
- User manual
- Technical manual
- Record of installation and setup
- Manufacturer's inspection report and declaration of conformity
- Guarantee
- Instructions for the use and care of the accessories
- MELAflash CF card for documentation purposes
- 1 litre storage container for rinse aid
- Filling funnel for the regeneration salt
- Starter package of regeneration salt
- Open-end wrench for the injector rail
- Container tap for a 5 l and 10 l container
- Effluent hose with hose clamp Ø 16-25/9

Optionally

- Accessories in accordance with the delivery slip

Views of the device



- 1 Operating and display panel
- 2 Cover plate (optional)
- 3 Power switch
- 4 Cover cap for CF card slot and Ethernet data connection (for service technician)
- 5 Door handle
- 6 Hinged door, opens forwards
- 7 Ventilation slots for air egress
- 8 Grip for the drawer
- 9 Drawer for process agents
- 10 Device foot

Fig. 1: View from front




Fig. 2: Cover cap CF card slot open

- 11 Slot
- 12 LED
- 13 Ejection button
- 14 Ethernet data connection




Fig. 3: View from rear

- 15 Demineralized water (DI water) connection
- 16 Cold water connection
- 17 Effluent connection
- 18 Transport rollers
- 19 Ethernet data connection for permanent network connection
- 20 Mains cable




Fig. 4: View inside

- 21 Injector rail connection fitting
- 22 Cold water inflow (CW) and demineralized water (DI)
- 23 Salt container
- 24 Rinse arm down
- 25 Coarse and fine sieve




Fig. 5: Process agent drawer open

- 26 Drying filter
- 27 Assignment of the process agents
- 28 Suction lance bracket
- 29 Container for rinse aid with suction lance
- 30 Container for neutralizer with suction lance
- 31 Container for cleaning agent with suction lance

Operating panel and acoustic signals

The operating panel consists of a two-row LED display and four membrane keys.














Fig. 6: 2-row LED display

Table 1: Key functions

	Key	Function / Explanation
1...2	 or 	Navigation: BACKWARDS, FORWARDS to adjust a value: SMALLER, LARGER
3		Unlock door BACK, ABORT Leave menu
4		Confirming messages (ENTER, OK, YES, SELECTION). QUIT with warning or malfunction messages
2+3	 and 	The system status is displayed with information relating to the device e.g. serial number, software version, daily and total batches etc.
1+3	 and 	QUIT + DOOR, i.e. to acknowledge the program abort and unlock the door
3...4	 or 	To delete all logs located in the internal log memory

Acoustic Signals

The device issues acoustic signals for information purposes.

Table 2: Acoustic Signals

Signal / tone	Meaning
1 x	Confirmation, warning or message
3 x	Refill with salt soon; program abort; end of the abort after drying abort reached
5 x	Program completed successfully
10 x	Malfunction

Menu structure

MAIN MENU


- | P01 Universal-Program
- | P02 Quick-Program
- | P03 Intensive-Program
- | P04 Ophthalmic-Program
- | Z01 Rinsing
- | Z02 Emptying
- | Z03 Conductivity DI
- | Z04 Air removal
- | Z05 Regeneration
- | Z06 Time metering 60s


M01 → DOCU MENU (output of a saved log via the following output media)

Select output medium : Automatic, CF card, MELAprint, PC

- | 01 Log list
- | 02 Last log
- | 03 Logs of day
- | 04 Logs of week
- | 05 Logs of month
- | 06 All logs
- | 07 Last malf. log
- | 08 Malf. logs of day
- | 09 Malf. logs of week
- | 10 Malf. logs of month
- | 11 All malfunction logs
- | 12 Caption log
- | 13 Status laog
- | 14 System log
- | 15 Format CF card

M02 → SETUP MENU

- | 01 DI water
- | 02 Autom. logging
- | L 
- | 03 Date
- | 04 Time
- | 05 Display contrast
- | 06 Language
- | 07 Water °dH
- | 08 → DIAGNOSIS+SERVICE

- | ACOUT AC outputs
- | DCOU DC outputs
- | AIN Analog. inputs
- | DINZ count. inputs
- | DIN Digital inputs
- | SERVICE MENU
- | L 
- | Maint. Counter Date
- | DEMO Mode

Water softening unit

The tap water is processed in a water softening unit to produce optimal cleaning results.

- ▶ Use coarse-grain regeneration salt (NaCl) to regenerate the water softening unit.



NOTICE

The fitted water softening unit has been optimized for a degree of hardness of 0–40°dH.

For higher degrees of hardness than 40°dH, you will require a dedicated water softening unit.



NOTICE

An incorrectly set degree of hardness can result in a higher salt consumption or limescale deposits on the instruments.

When using a dedicated water softening unit, the residual hardness of the dedicated water softening unit must also be set in the setup menu.

Table 3: Water hardness conversion table

°dH	mmol/l	°f	°e	°dH	mmol/l	°f	°e	°dH	mmol/l	°f	°e
1	0.2	2	2	15	2.7	27	19	28	5.0	50	36
2	0.4	4	3	16	2.9	29	20	29	5.2	52	37
3	0.5	5	4	17	3.1	31	22	30	5.4	54	38
4	0.7	7	5	18	3.2	32	23	31	5.6	56	39
5	0.9	9	7	19	3.4	34	24	32	5.8	58	41
6	1.1	11	8	20	3.6	36	25	33	5.9	59	42
7	1.3	13	9	21	3.8	38	27	34	6.1	61	43
8	1.4	14	10	22	4.0	40	28	35	6.3	63	44
9	1.6	16	12	23	4.1	41	29	36	6.5	65	46
10	1.8	18	13	24	4.3	43	31	37	6.7	67	47
11	2.0	20	14	25	4.5	45	32	38	6.8	68	48
12	2.2	22	15	26	4.7	47	33	39	7.0	70	49
13	2.3	23	17	27	4.9	49	34	40	7.2	72	51
14	2.5	25	18								

5 First steps

Setup and Installation



PLEASE NOTE


Comply with the specifications of the technical manual during set-up and installation. This contains all building-side requirements.

Record of installation and set-up

The responsible stockist is to complete the record of installation and setup as proof of the correct setup, installation and commissioning. A copy is to be sent to MELAG Medizintechnik oHG and the stockist. This is a constituent part of any guarantee claim.

Switching the device on and off

- ▶ Switch the device on or off at the power switch.



Opening and closing the door

The door is automatically closed via a motor. For this reason, it is important that the device is connected to the power supply and is switched on. The door unlocks automatically after a successful program run. The door cannot be opened following a power outage. In such a case, activate the [door emergency release](#) [▶ page 20].




NOTICE

The door can only be opened during a program run using a program abort.

The door is unlocked after the program abort has been acknowledged and sufficient cooling has been performed.

Opening the door

1. Switch on the device at the power switch.
2. Press the  key.
 - ↳ The door is unlocked.
3. Open the door forwards to open.

Closing the door

- ▶ Shut the door upwards and press it until the motorized lock sets in.

Emergency release of the door



CAUTION

Danger of scalding from hot steam.

Hot steam could be released upon opening the door.

- Never operate the emergency release while a program is active.
- Wear suitable protective clothing.



DANGER


Danger of contamination from insufficient disinfection

A program aborted by emergency release is classed as not completed successfully. Process water can remain in the washing chamber; the instruments have not been disinfected completely.

- The instruments must be subject to renewed treatment.

Proceed as follows to activate the emergency release:

1. Pull out the process agent drawer.




➔ An emergency release grip is located at the front left-hand side of the device.

2. Press down on the grip until you hear a clicking sound.
3. To do so, pull the door forwards strongly using its handle.

Inserting the basis basket

A port for the connection of the injector rail is located on the right-hand rear side of the washing chamber of the washer-disinfector.

- ▶ Slide the basis basket with the injector rail opening into the washing chamber until it connects to this port.



Filling the regeneration salt



CAUTION

Danger of injury from insufficient protective measures!

Performing work without first taking the corresponding protective measures can result in injuries.

- Comply with the working safety measures required by the respective tasks.



NOTICE

Malfunctions of the water softening unit from unsuitable regenerating salt.

Fine grain regeneration salt can cause device malfunctions. We do not recommend the use of pellets, as the salt dissolves too slowly.

- Use only special, coarse grain regeneration salt (additive-free NaCl).
- Never use cooking salt, table salt, de-icing salt, cattle salt or road salt. These salts usually contain insoluble components.
- Never pour cleaning agent or other process agents in the salt container.

Filling the regenerating salt for the first time (by the service technician upon first installation)

1. Fill enough water in the salt container **once** until it overflows. This enables dissolution of the salt.
 - ↳ No further water is required for all subsequent salt fillings.
2. Fill the salt container with 1 kg of regenerating salt.
 - ↳ The device can only be operated if the salt container contains sufficient regenerating salt.


Re-filling regenerating salt



Insufficient regenerating salt will result in the display of the corresponding display notification.

- Display of *Refill with salt soon* means that you can perform up to ten washing runs depending on the program selected and the hardness of the water.
- Display of *Salt exhausted. Please refill* means that you should refill with regenerating salt immediately. Otherwise you will be unable to start a further program.

You can refill the regenerating salt at any time without the corresponding display notification. After filling, start the supplementary “Rinse” program to remove overflowing brine and salt residue from the washing chamber.

In order to refill the regeneration salt, proceed as follows:

1. Acknowledge the notification `Refill with salt soon` with the  key.
2. Open the door.
3. Remove the basis basket.
4. Open the lid of the salt container by turning it anti-clockwise.
5. Insert the filling funnel for the regenerating salt in the opening.
6. Fill the salt container with regenerating salt using the filling funnel.
7. Clean the edge of the filling opening of salt residue.
8. Remove the filling funnel and screw on the lid of the salt container.
9. Insert the basis basket.
10. Wait three minutes and then start the rinsing program without (instrument) load.



Regenerating the water softening unit

The internal water softening unit regenerates automatically in certain intervals. The program run time is extended by a number of minutes. You can regenerate the water softening unit manually after e.g. having filled it with salt without a warning having previously been issued.

- ▶ To do so, start the “Regeneration” program.

Metering process agents

The concentration of the process agents is set once during the initial device setup performed by the customer services technician (see technical manual). During a program run, the preset concentration of the relevant process agents is metered automatically.

Holding process agents ready



DANGER

Danger of acid burns by irritant substances!

Improper handling of the process agents may cause caustic burns and health damage.

- Protect your eyes, hands, clothing and all surfaces from contact with the process agents.
- Comply with the information from the manufacturer of the process agents.
- Every type of fluid (e.g. in the drawer, in the device floor tank or fluid emerging from the device) issued from the device as the result of damage could potentially contain aggressive process agents.



NOTICE

Damage to instruments and the device from unsuitable process agents!

- Use only those process agents approved by MELAG for use in this device.
- Comply with the information from the process agent manufacturer.

An insufficient filling level of a certain process agent means will trigger the display of the corresponding display message. In this case, replace or refill the process agent container.

Comply with the following provisions pertaining to the use of the process agents:


- When re-filling, please use the same process agents, set on the device during installation (see the tie-on label on the can in the drawer).
- Changing the combination of process agents may only be performed by a qualified and trained service technician.
- Every change of a process agent in a validated device necessitates revalidation.
- Use a citric acid based neutralizer for the treatment of dental transfer instruments.
- Use a mild alkali cleaner wherever possible on technical, hygienic and ecological reasons.
- Check that your instruments are compatible for use with rinse aid.
- All air should be removed from the metering hoses before commissioning or after changing a container. See [Bleeding the metering hoses](#) [▶ page 26].

Containers for process agents

Every process agent has its own container and a suction lance with screw-on lid:


- Cleaning agent: 5 L container with a blue suction lance screw-on lid
- Neutralizer: 5 L container with a red suction lance screw-on lid
- Rinse aid: 1 L container with a black suction lance screw-on lid

- ▶ Place the container in the drawer in accordance with with the process agent assignment. A container can only be locked correctly if the colour of the process agent matches that of the suction lance screw-on lid.



Change the containers for cleaning agent and neutralizer


1. Unscrew the suction lance from the container and place it the suction lance bracket.



2. Place the new container in the process agents drawer and screw on the suction lance.

➤ The screw-on lid of the suction lance points forwards.

3. Vent the metering hose (see [Bleeding the metering hoses](#) [▶ page 26]).



Refilling rinse aid



PLEASE NOTE

Rinse aid may not be used to decontaminate ophthalmological instruments, see [Decontaminating ophthalmological instruments](#) [▶ page 29].




PLEASE NOTE

A streaky instrument surface could be caused by too much rinse aid.

Proceed as follows to fill an empty container with rinse aid:

1. Unscrew the suction lance from the container and place it the bracket behind it.




2. Transfer the rinse aid from the original packaging into the MELAG container.
 - ↳ Use the container tap (for a 5 l and 10 l container) included in the scope of delivery for ease of filling.
 - ↳ Fill the container with rinse aid $\frac{3}{4}$ full, otherwise the rinse aid will overflow during insertion of the suction lance.
3. Screw the suction lance onto the container.
4. Vent the metering hose (see [Bleeding the metering hoses](#) [▶ page 26]).


Bleeding the metering hoses

The process agent hoses must be vented after filling the container, performing a product change or removing the suction lances. Air removal removes all air pockets from the metering hoses and permits good metering.

-
- ✓ *The rinse aid suction lance not used for ophthalmic instruments must be inserted head-first in the suction lance bracket during the “Air removal” program run*
 - ✓ *The “Air removal” program must be started twice after initial filling of the container.*
-

1. When processing ophthalmic instruments, insert the suction lance of the non-used process agent head-first in the suction lance bracket **before the program start**.



2. Press the  button repeatedly to navigate in the main menu to Z04 Air removal.
3. Start the “Air removal” program by pressing the  key.

6 Cleaning and disinfection

The nature of the load

Comply with the specifications of the document *Instructions for the use and care of accessories* when loading the device. Only use the loading pattern specified and approved within the scope of the validation. This device can clean and disinfect max. 10 kg of the following type of load:

- Massive instruments
- Hollow-bodied instrument e.g. aspirator tips, which are fixed to injector nozzles or
- transfer instruments e.g. handpieces by using the adapter

Further accessories may be required when **reprocessing ophthalmological instruments** (not available from MELAG). The operator is responsible for validating the procedure in combination with special load accessories. It is especially important that feed lines to hollow-bodied instruments are maintained without kinking and as short as possible.

Arranging the load



DANGER

Danger of injury from sharp and pointed instruments.

Using inappropriate procedures to load the device with sharp and pointed instruments can result in injuries. Baskets and sieve cassettes with wire meshes or other openings do not provide any protection against penetration by sharp instruments.

- Wear protective gloves.
 - Proceed with caution when sorting sharp and pointed instruments.
 - Arrange sharp and pointed instruments in such a way as to prevent the danger of injury.
-



NOTICE

Use only those instruments designed by their manufacturer for automatic treatment in a washer-disinfector.

- Some brands are only authorized for thermal disinfection after a specific year of manufacture.
 - Comply with the information from the relevant instrument manufacturer.
-

To arrange the load, the basis basket must be used as a minimum (see chapter [First steps](#) [▶ page 19], [Inserting the basis basket](#) [▶ page 21]). This is used to accommodate insert racks, instrument baskets, wash trays and sieve cassettes. The basis basket with an injector rail is available for the decontamination of hollow-bodied instruments.

Further accessories and their user instructions such as insert racks for wash trays, sieve cassettes and instrument baskets etc. are listed in the document *Instructions for the use and care of the accessories*.

Please comply with the following specifications when arranging the load:

- ▶ Empty all residual liquids from containers before arranging them in the device. Rinse away any liquids (e.g. disinfectant solutions) thoroughly.
- ▶ Never place any individual instruments directly in the basis basket. Use baskets or trays to this end.
- ▶ Ensure that instruments do not protrude from the sides of the instrument basket or the basis basket. Protruding instruments can damage the seal and the surface of the door or the side walls of the washing chamber. The instruments can break.
- ▶ Place hollow-bodied instruments in the device in such a way as to ensure safe rinsing. If necessary, use the accessories developed especially for the decontamination of hollow-bodied instruments such as injector nozzles, Luer connections, adapters etc. See the information for the use and care of accessories.

- ▶ Avoid blockages of the rinse arm from instruments protruding upwards or downwards. The rinse arms must be able to rotate freely.
- ▶ Avoid spray shadow. A good cleaning outcome depends on the correct arrangement of the instruments.
- ▶ Arrange all containers such as glasses, basins etc. with their opening pointing downwards.
- ▶ Place components with openings or compressions at an angle, so that the water can run off them.
- ▶ Only use thermostable instruments approved by their manufacturer for decontamination.

Decontaminating hollow-bodied instruments



DANGER

Danger of contamination from insufficient disinfection

- Check the hollow-bodied instrument for free passage before decontamination.
- Residue on the hollow-bodied instruments can hinder water pass through and thus impair their disinfection.



PLEASE NOTE

All openings must be connected to an accessory piece when using a multiway terminal blocks or the injector rail. Only then can correct functioning be guaranteed.

- Seal non-used openings with sealing screws.



PLEASE NOTE

Use a filter insert for hollow-bodied instruments with an inside diameter ≤ 0.8 mm.

- Do not use any reusable filter screen and central filters in the treatment of ophthalmological instruments.

Comply with the following specifications before automatic decontamination:

- ▶ Flush all hollow-bodied instrument after use with patients / before automatic decontamination.
- ▶ Treat only those hollow-bodied instruments which guarantee sufficient and reproducible rinsing. Remove instruments with a recognisably reduced throughflow.
- ▶ Use only MELAG adapters (accessories list) for the injector rail to treat hollow-bodied instruments. The suitability of the hollow-bodied instruments for the respective adapter and the sufficient flushing of the instruments can only be proven by validation.
- ▶ Check the connections between the adapter and the hollow-bodied instrument for stability both before and after decontamination. Should a connection work loose after decontamination, the instruments must be decontaminated again.
- ▶ Comply with the cleaning and replacement intervals when using filter inserts. The cleaning and replacement intervals are listed in the separate document *Instructions for the use and care of the accessories*.
- ▶ When treating dental and ophthalmologic transfer instruments, comply with the specifications of the special decontamination instructions in the chapter [Decontaminating dental transfer instruments](#) [▶ page 29] und [Decontaminating ophthalmological instruments](#) [▶ page 29].

Rules for the use of filters / filter discs:

Diameter of the inner lumen	Use of a filter
≤ 0.8 mm	Filter required e.g. a triple distributor with a filter disc
> 0.8 mm	No filter required, direct connection of the adaptor to the injector rail possible.

Decontaminating dental transfer instruments

Comply with the following specifications before automatic treatment:

- ▶ The exterior surfaces of the handpieces should be free of all residue e.g. dental cement.
- ▶ The air and spray channels must be entirely clear.
- ▶ Avoid the drying-on of soiling, especially on the handpieces.
- ▶ Use a citric acid based neutralizer for the treatment of dental transfer instruments.
- ▶ In case of any residual moisture after decontamination, e.g. after a drying abort, dry hollow-bodied instruments with medical compressed air.

Care of the instruments and adapters

Dry the spray / air / water channels immediately after the successful cleaning and disinfection using clean (medical) compressed air and then maintenance with suitable maintenance and care products and oils.

The adapter for handpieces should be checked for dirt at regular intervals. Clean the individual adapter components under running water where necessary. The silicon inserts can be rubbed with a damp, non-fuzzing cloth.

Decontaminating ophthalmological instruments

Please comply with national recommendations for the cleaning of medical products under the aspect of decontamination of infectious prion proteins (vCJK).



DANGER

Danger of contamination from biological interactions.

Devices used to treat ophthalmologic instruments may only be used exclusively for this purpose.

- Do not treat any instruments which are used for operations of the posterior segment surgery (contact with retinal tissue, sub-retinal fluid and the optical nerve).
- These devices should be fitted with a suitable filter system: Ophthalmological instruments may only be treated using the filter discs (MELAG art. no. 64375).



NOTICE

Use only those instruments designed by their manufacturer for automatic treatment in a washer-disinfector.


- Some brands are only authorized for thermal disinfection after a specific year of manufacture.
- Comply with the information from the relevant instrument manufacturer.



NOTICE

Do not use a rinse aid to treat ophthalmological instruments.

1. If present, remove the rinse aid container from the process agents drawer and hang the black suction lance in the suction lance bracket so that the cover is positioned at the top.



2. The suction lance need only be placed head-first in the suction lance bracket if the metering hoses are to be bled. See [Bleeding the metering hoses](#) [▶ page 26].



PLEASE NOTE

Use demineralized water to treat ophthalmological instruments.

- To this end, connect a mixed-bed resin cartridge.

Suitable program

Treat ophthalmological instruments in the Ophthalmo-Program. Only this program enables monitoring of the conductivity during the disinfection phase; this ensures an uncritical residual conductivity.

Comply with the following specifications before automatic decontamination:

- ▶ Cleaning should be undertaken with a mildly alkali cleaning agent. Neutralization should be effected with a citric acid based neutralizer.
- ▶ Flush all hollow-bodied instrument after use with patients / before automatic treatment with DI water.
- ▶ Treat only those hollow-bodied instruments which guarantee sufficient and reproducible rinsing. Remove instruments with a recognisably reduced throughflow.
- ▶ All hollow-bodied instruments should be connected with the rinse bar deigned for the purpose.
- ▶ Ensure that plugs and / or cables from phaco handpieces are not able to slip through the basis basket, otherwise the rinse arm can become blocked.
- ▶ Try to prevent dirt from drying or encrusting on and in the instruments.
- ▶ Dry any residual moisture, e.g. after a drying abort or hollow-bodied instrument after decontamination, with medical compressed air.
- ▶ When using the rinsing system, individual outlets which are not connected can be sealed with suitable accessories.

Instrument care

Comply with the manufacturer's instructions regarding the care and maintenance of the instruments / the load accessories.

Routine controls

Perform a routine check of the pH value after decontamination of the hollow-bodied instruments.

1. Blow through the hollow-bodied instruments with medical compressed air onto indicator paper (e.g. from Macherey-Nagel: PEHANON pH 4.0-9.0). The measurement accuracy must amount to or exceed 0.5.
2. Compare the value displayed on the indicator paper with the pH value of the final rinse water from the previous performance qualification.
3. Should you discover any deviations, contact the customer services.

Overview of programs

- ▶ Select the program according to the level of soiling of the load. Comply with the specifications from the validation:
- ▶ The Universal-Program is sufficient for every-day general cleaning and disinfection. The Quick-Program is designed for lightly-soiled instruments.

The following table lists the correct program for each load.

Table 4: Programs and operating times

Program	Nature of the instruments / degree of soiling	Operating time* inc. drying time	
		DTA	DTB
Universal-Program 90 °C, 5 min. ¹	<ul style="list-style-type: none"> ▪ For normal-strongly soiled instruments. This complies with the general hygiene-related requirements of DIN EN ISO 15883-1. 	40 min.	59 min.
Quick-Program 90 °C, 5 min. ¹	<ul style="list-style-type: none"> ▪ For unsoiled or only lightly-soiled instruments 	36 min.	53 min.
Intensive-Program 90 °C, 5 min. ¹	<ul style="list-style-type: none"> ▪ For especially heavily soiled instruments ▪ Like the Universal-Program but with a longer cleaning time. 	51 min.	64 min.
Intensive-Program 90 °C, 5 min. ¹	<ul style="list-style-type: none"> ▪ For ophthalmological instruments ▪ Like the Universal-Program but with a longer cleaning time, double intermediate rinsing without subsequent rinse aid. 	42 min.	59 min.

*The operating times represent average values and apply for the recommended running water pressure at a cold water temperature of 15 °C.

Table 5: Additional programs

Additional programs	Application	Operating time*
Rinsing, 3 min. no disinfection, without process agents	<ul style="list-style-type: none"> ▪ For rinsing strongly-soiled instruments (e.g. blood). A disinfection program must then be started ▪ To rinse out the washing chamber after adding salt; without process agents, no disinfection 	3 min.
Emptying	<ul style="list-style-type: none"> ▪ Pumping out residual water from the washing chamber 	1 min.
Conductivity measurement DI	<ul style="list-style-type: none"> ▪ For measuring the conductivity of the DI water 	2 min.
Air removal	<ul style="list-style-type: none"> ▪ After filling / changing the process agents, i.e. product change etc. ▪ With decommissioning and commissioning 	5 min.
Regeneration	<ul style="list-style-type: none"> ▪ Regenerating the internal water softening unit 	8 min.
Time of metering	<ul style="list-style-type: none"> ▪ Only for technicians 	--

*The operating times represent average values and apply for the recommended running water pressure at a cold water temperature of 15 °C.

¹⁾ In accordance with the A0 concept from EN ISO 15883-1 thermal disinfection is performed with 90 °C (+ 5 °C, - 0 °C) and an application times of 5 minutes (min. A0-3000).

Selecting, starting and following the program



NOTICE

Unsupervised operation of the device can result in damage to the device or your facility. In such a case, MELAG does not accept any liability.


- Never operate the device unattended. Unattended operation is performed at the operator's risk.


Ensure compliance with the following prerequisites in order to secure the optimal cleaning performance before every program start:


- ▶ The process agents containers are sufficiently full.
- ▶ The injector rail nozzles / adaptor are clean.
- ▶ The rinse arms can be turned freely.
- ▶ The load must be arranged correctly.
- ▶ Baskets and inserts are inserted correctly.

Selecting and starting a program:

- Select a program in accordance with the [Program overview](#) [▶ page 31].

- Navigate to the desired program using . The display shows the program names, the temperature and the holding time.




- Start the selected program with .

Following the program run on the computer

You can follow the current progress of a program run on every computer in the practice network. To do so, an IP address must be issued for the device and it must be incorporated in the practice network.

- Open a web browser window in the practice PC (we recommend Mozilla Firefox or Internet Explorer).
- Enter the device IP address in the address bar of the web browser, e.g. 192.168.70.206 and confirm with [Enter].



- ➔ The program run and the device information such as e.g. serial number, device software version will be shown.

Manual program abort



NOTICE

Cancellation of a current program by deactivation at the power switch may cause damage at the device.

- Never abort a program by switching off at the power switch.
-

Abort the program during drying



CAUTION

Danger of burns from hot instruments and surfaces.

The instruments and chamber can remain hot even after the device has been switched off.

- Allow the device to cool before removing the instruments.
 - Wear suitable protective gloves.
-






WARNING

Nucleation because of residual dampness.

If a program is aborted during drying, residual dampness will remain on the instruments.

- Only abort a current program in exceptional reasons.
 - Dry the instruments manually.
-

If a program is aborted during drying, the program is classed as having been ended successfully. Proceed as follows to abort the program during drying:

1. Wait until the display shows the message **CANCEL DRYING ●4**.
 2. Press the  button to abort the program and confirm the abort with **Yes**.
 3. To open the door, press  and  simultaneously.
-

Abort the program before the start of drying



CAUTION

Danger of burns from hot instruments and surfaces.

The instruments and chamber can remain hot even after the device has been switched off.

- Allow the device to cool before removing the instruments.
 - Wear suitable protective gloves.
-




WARNING

Danger of infection from program abort.

If a program is interrupted before the start of a drying phase, the load is taken as not having been cleaned and disinfected. This endangers the health of the patient and the practice team.

- Only abort a current program in exceptional reasons.
 - Never open the door after a program start.
 - Treat the instruments again after a program abort.
-

Press the  buttons to abort a current program before drying begins and follow the information on the display.

Removing the load after program end



CAUTION

Danger of burns from hot instruments and surfaces.

The instruments and chamber can remain hot even after the device has been switched off.

- Allow the device to cool before removing the instruments.
- Wear suitable protective gloves.



PLEASE NOTE

Open the door immediately after the end of the program to prevent the accretion of condensation.



PLEASE NOTE

Dry the hollow-bodied instruments with clean (medical) compressed air after treatment in order to remove any residual moisture.




NOTICE

Soiling can become encrusted even after a rinsing program; this could result in instrument corrosion.

- Remove the instruments from the washing chamber after every program run.
- Do not leave the instruments in the washing chamber overnight.

The display message indicates whether and when a program has ended successfully. The display shows the last batch number run and the total batch counter after a program end or the end of a program abort.

1. Press  and open the door.
2. Remove the load whilst complying with all the hygiene and working safety regulations.
 - Do not leave the instrument in the washing chamber over night, otherwise this can result in corrosion of the instruments.

The load is only classed as having been cleaned and disinfected if all the following points have been satisfied. Otherwise, treat the load afresh.

- ▶ The program has been performed without interruption or malfunction.
- ▶ Instruments must be completely clean and dry.
- ▶ Hollow-bodied instruments are fixed.
- ▶ The interior of the hollow-bodied instruments are accessible.
- ▶ The injector rail still sits snugly on the connection tube in the washing chamber.
- ▶ The nozzles and connections to the basis basket are still connected.

7 Logging

Batch documentation

The batch documentation acts as proof of the successful conclusion of the sterilization process and represents an obligatory part of quality control. The data, such as type of program as well as batch and process parameters of the completed programs, are stored in an internal log memory of the device.

To obtain the batch documentation, you can read out the internal log memory and transfer its data to various output media. This can be performed immediately at the end of every program or at a later point, such as at the end of the day.


As delivered, the MELAflash CF card is set as the output medium. Setting a different output medium or activating the internal log memory is outlined under [Settings](#) [▶ page 42].

Capacity of the internal log memory

The device is equipped with an internal log memory. This saves all the data regarding the sterilization program automatically. The capacity of the internal log memory is sufficient for 15-20 logs. If the internal log memory is full, the display will show the warning `Internal program log memory full, do not output all logs`. If this warning is issued, provide the specified output medium (see [Settings](#) [▶ page 42] menu → Logging) and output the affected log (→ [Logging](#) [▶ page 42] menu). If the program is continued, the logs are deleted automatically; the last ten saved logs remain in the log memory.

We recommend outputting logs immediately.

Outputting logs immediately and automatically

- ▶ The text log is issued on the selected output medium after the end of the program run. At the same time, this text log is saved in the internal log memory and marked as “outputted”.
- ▶ If multiple output media have been activated, all activated output media must be connected to the device. Otherwise, the text logs are saved in the internal memory and are classed as “not outputted”.
- ▶ If the internal log memory is full, the device will register all the text logs which are classed as “not outputted”. The warning 386 appears after the program start. You can acknowledge this warning with the  key in order to continue the program run.
- ▶ With warning 372, logs not yet outputted must be outputted manually. Only then is a program start possible. The log memory is deleted automatically after manual issue; the last ten logs remain in the log memory. The manual outputting of logs is outlines under [Subsequent log output](#) [▶ page 38].

Output media

You can output the logs of the finished programs via the following output media:

- MELAflash CF card
- A computer via the practice network (LAN)
- MELAprint 42/44 log printer with network adaptor

The output media can be combined in any fashion. Thus it is possible both to save logs on the CF card (included in the scope of delivery) and also to print them on the log printer.



PLEASE NOTE

Further information about the log printer (e.g. the duration of the readability of the log print-out) is specified in the appendant operating manual.

Using the CF card as an output medium



NOTICE


Premature removal of the CF card from the card slot or its inappropriate handling can result in data loss, damage to the CF card, the device and/or its software.

- Never push the CF card in the slot with force.
- Never remove the CF card from the slot whilst it is being written or read. The red LED next to the card slot on the right will illuminate red in short irregular intervals during reading and writing access.
- Hold the CF card cover cap closed during operation so that the CF card ejection button is not actuated by mistake.

Inserting the CF card

The card slot for the CF card is located behind the cover cap on the right, adjacent to the door below the power switch. When inserting the CF card in the slot, ensure that it is aligned correctly.

1. Open the CF card cover cap.
2. Insert the CF card in the slot with the contacts at the front. The MELAG logo on the CF card points towards the LED.



3. Slide the CF card in the card slot until it clicks. Do not use force. When the CF card has been placed correctly, the red LED will illuminate shortly.
4. Close the cover cap.

Removing the CF card

1. Open the CF card cover cap.
2. Press the ejection button and remove the CF card.
3. Close the cover cap.

Using the computer as an output medium

You can either connect a computer directly to a device or via a network if the following conditions are fulfilled:

- ▶ The computer is fitted with a network card with a RJ45 bushing (LAN).
- ▶ An FTP server or an FTP service is installed on the computer (when the log is issued via FTP).
- ▶ A suitable program, e.g. MELAtrace/MELAview, is installed (when the log is issued via TCP).

Outputting logs immediately and automatically

As delivered, the MELAflash CF card is set as an output medium in the setup menu and thus the automatic log output at the end of a program (immediate output = YES) is thus activated. Log output on multiply activated media is performed successively. You can select or add an alternative output medium for automatic log output.

Text logs

The following requirements must be fulfilled in order to output text logs immediately after the end of a program.

- ▶ In the setup menu → 02 Autom. logging → immediate output is set to YES.
- ▶ In the setup menu → 02 Autom. logging → at least one output medium is selected and 02 Autom. logging is thus set to ACTIVE.
- ▶ The activated output medium is available (e.g. The MELAprint 42/44 log printer or a CF card).

Graphic Logs (optional)

The following additional requirements must be fulfilled in order to record graphic logs:

- ▶ In the SETUP MENU → 02 Autom. logging → Graphic logs, at least one output medium is set to YES.
- ▶ At least one of the output media selected corresponds to an output medium for the text logs. This means that at least the computer or the CF card must be activated as an output medium for both log types.
- ▶ The selected output medium has been connected.



NOTICE


Graphic logs cannot be saved in the internal log memory and cannot be outputted via the log printer MELAprint 42/44.

- Save the graphic logs on the CF card or the computer.

The following settings can be made to record graphic logs:


Graphics & CFC

One of the selected output media must conform with the selected output medium for text logs.




CFC interval

CFC interval or PC interval indicates the time intervals in which the program curve is recorded on the CF card or computer. The smaller the time interval, the more exact the curve. In the example, the time interval is set at 1 second.



PC backup

PC backup indicates the time interval in which the graphic logs are to be saved on the computer by the device. In the example, the backup interval is set to 30 seconds.



Subsequent log output

The Docu menu provides the option of issuing logs subsequently and independently of the point of the program end. Proceed as follows:

1. Press or to navigate to the `DOCU MENU`.
2. Press to open the `DOCU MENU`.
3. Press repeatedly to select an output medium. If you want to assume the settings from the menu `Automatic logging`, select the option `automatic`.
4. Press to get to the option `Log type`.
5. Press to choose between the log types, e.g. `last log`, `log of the day`, etc.
6. Press the key in order to start the log output.

Deleting the saved logs

Save the logs on an output medium before deleting them.

1. Press or to navigate to `DOCU MENU`.
2. Press to open `DOCU MENU`.
3. Press again.
4. Press to navigate to option `06 All logs`.
5. Press and for a short period. A security query is displayed: .
6. Hold and depressed to delete all logs.

Determining the format for the program logs

The log format enables you to determine which of the data saved is to be outputted. You can choose between the format (0001) and the more comprehensive format (0002). The log format (0002) is the standard format. Working in the setup menu, you can select the log format for the program log (see [Logging](#) [▶ page 35]).

Log types

In addition to logs for successfully completed programs, there are many other types of log. These can be issued via the list in the Docu menu. You can identify the log type by the ending of its file name.

Table 6: Overview of log types

Ending	Stands for	Explanation
PRO	Program log	Log of a successfully completed program
GPD	Graphic log	The log in which the process is recorded graphically
STR	Malfunction log	Log of an aborted program
STB	Malfunction in standby	Log with malfunctions without a program having run
LOG	System log	List of all the malfunctions and changes to the system in order of time (log book)
STA	Status log	Summary of all the important settings and system states (counter, measured values etc.) + a list of all procedure-relevant parameters (VRP)
LEG	Legend log	Contains all step abbreviations used in the program log
DEM	Demo log	Log of a successfully completed simulated program in DEMO mode (only for presentation purposes)
DES	Demo malfunction	Log of a program simulated as interrupted (presentation)

Example of a program log for a successful finished Program

10	MELAG MELAtherm 10-DTA	Header: Device name
15	Program: Ophthalmo-Program	Program name
20	TARGET °C min	Target values of the temperature and holding time of the partial cycles
21	Pre-clean: 22.0 01:30	Date
23	Cleaning : 55.0 10:00	Daily and total batches
28	Disinfect: 90.0 05:00	Control display
30	Date : 17.02.2016	Program abort (appears if program unsuccessful)
35	Batch : Day:10 Total:00010	ACTUAL-value Temperature (range) in C°, holding time of the partial cycles
40	Program successfully ended	ACTUAL-value Temperature conditions of the disinfection, A0 value
50	ACTUAL °C +/- K min	Conductivity of the DI water for the final rinse
51	Pre-clean: 29.1 +0.2/-0.2 01:30	Time upon program start
53	Cleaning: 57.2 +1.2/-0.3 10:00	Time upon program end
58	Disinfect: 91.9 +0.2/-0.2 05:00 3962	Device serial number
60	Conductivity : 9.8 (12.6) µS/cm	Installed firmware version
65	Start time: 11:19:21	Installed parameter version
70	End time : 12:49:32 (90:11 min)	Installed user interface
80	SN:201610-DTA1288	Time (mm:ss), since the program start Duration (mm:ss), of a min.
81	Firmware : V1.309 01.03.2013	°C Temperature of the rinse liquor in the washing chamber in Celsius
82	Parameter: V1.314 25.02.2015	ml Volume of CW/DI water, the process agent consumed during a process step
83	UI : V1.305 25.01.2013	mbar Rinse pressure
92		Up to 5 warnings as required
95		If necessary, malfunction numbers of the program aborts
		Proof of authenticity: Authenticity proof; should never be changed; permits inference that the data was created and not changed on an autoclave from MELAG. Sensor measurement values are displayed in the case of a malfunction. The values are helpful for a technician.
	Step Start End Time °C ml mbar	
	--> Process start	
	S1X01 0:00 0:00 0:00 52.3 ---- --	
	--> Pre-cleaning	
	V1E01 0:00 0:03 0:03 52.3 ---- --	
	V1P01 0:03 0:09 0:06 52.2 ---- --	
	V1S12 0:09 0:13 0:04 52.2 350 --	
	...	
	--> Cleaning	
	R1E02 4:19 4:19 0:00 35.1 ---- --	
	R1P01 4:19 4:25 0:06 35.3 ---- --	
	R1U03 4:25 4:54 0:29 35.4 ---- --	
	R1S02 4:25 5:16 0:51 25.2 5555 --	
	...	
	--> Neutralizing	
	N1E02 16:14 16:14 0:00 56.8 ---- --	
	N1P01 16:14 16:20 0:06 56.8 ---- --	
	N1U01 16:20 16:48 0:28 56.3 ---- --	
	N1S03 16:20 17:05 0:45 35.9 5055 --	
	...	
	-.> Intermediate rinsing	
	Z1E02 18:03 18:04 0:01 35.1 ---- --	
	...	
	--> Disinfecting	
	D1E02 19:47 19:48 0:01 24.3 ---- --	
	D1P01 19:48 19:53 0:05 24.5 ---- --	
	D1U03 19:53 20:29 0:36 25.5 ---- --	
	D1S05 19:53 20:52 0:59 25.8 5035 --	
	...	
	--> Drying	
	T1E01 35:15 35:28 0:13 89.0 ---- --	
	T1T04 35:15 38:25 3:10 79.8 ---- --	
	T1K02 35:15 38:25 3:10 79.8 704 --	
	...	
	--> Process end	
	E1Z01 60:52 61:00 0:08 61.0 ---- --	
	>> Never change code in following row << 490000ED0077000A6DFB0209030004C200000000 >> Proof of authenticity batch log <<	
	Voltage max/min: 241/217	
	KW:21.5 VE:5035	
	0.0 0.0 -0.00 0.0 -0.00	
	--et1---et2---eps---etu---epx---END--	


Finding the logs



PLEASE NOTE

Do not rename the directory, otherwise logs will be stored in both the renamed directory as well as the new device directory generated automatically by the device.

All storage media (CF card or computer) contain a directory with the encoded serial number of the device concerned following log output. The folder name consists of five characters identical with the first five characters of every log, e.g. CR0ZH. This directory contains a sub-directory with the month of log generation e.g. 01_2013 for January 2013. This contains all logs generated by the device this month. The device directory is entered in the main directory on the CF card.



The device checks the storage medium after every type of log output (immediate output after a completed cycle or the transfer of multiple logs simultaneously). Should a directory not exist, it creates one of device and a month. If logs are issued on the same storage medium more than once, a "duplicated" directory will be created under the device directory in which these logs will be saved only once.






Given direct log transfer to a computer, set the storage location in the program (FCP, FTP) used on your computer.

8 Settings

SETUP menu

The setup menu contains settings for the date, time and display contrast.







Navigate in the setup menu as follows:

1. Press  to navigate in the main menu to M02->SETUP-MENU.
2. Press  to open the setup menu.
3. Press  to leave the setup menu.
4. Press  to save changes or hold  depressed to discard the changes.

Setting the water supply

If the device is connected to a DI water supply e.g. MELAdem 53/MELAdem 53 C or another water treatment unit, this must be set on the device. In its delivery state, the water supply has been set to DI water YES.

To alter this setting proceed as follows:

1. Press  to open the setup menu.
 - ↳ This display registers the option 01 DI-Water YES.
2. Press  in order to change the option.
 - ↳ The value YES flashes.
3. Press  or  to switch between YES and NO.
4. Press  to accept Yes o No.
 - ↳ The value no longer flashes.
5. Press  to leave the setup menu.
 - ↳ The selected value is automatically saved when leaving the setup menu.

Setting automatic logging


Settings relating to log output are made in Menu 02 Automatic logging. The settings made here are saved for the respective output medium. The display image shows whether the option for log issue is ACTIVE. Detailed information regarding logging is provided in chapter [Logging](#) [▶ page 35].








Determining the output medium

You are able to output the logs of the completed programs on various media: Comply with the specifications of the manufacturer's operating manual of the respective device.

The example shows how to use the CF card as an output medium. Proceed in a similar manner to set a different output medium.

Working in the [SETUP menu](#) [▶ page 42] set the output medium as follows:

1. Press  to navigate to 02 Automatic logging.

2. Press  to open the menu 02 Automatic logging.
 - ↳ The selectable output media are displayed consecutively.
3. Press  repeatedly to navigate in the setup menu to 01 CF card YES .
 - ↳ The display of YES indicates that the log is to be saved on the CF card.
4. Press  if this value is to be changed.
 - ↳ The value YES flashes.
5. Press  or  to switch between YES and NO.
6. Press  to save the new value.
 - ↳ The value no longer flashes.
7. Press  to leave setup menu 02 Automatic logging.
 - ↳ The selected value is automatically saved when leaving the setup menu.

Determining the log format

You can find more detailed information regarding the log formats 0001 and 0002 here: [Determining the format for the program logs](#) [▶ page 38].

Setting date and time

Date and time of the device must be correctly set for proper batch documentation.












PLEASE NOTE

The time is not set automatically.

- The time setting to summer or winter time must be performed manually.








Setting the date

Working in the [SETUP menu](#) [▶ page 42] set the output medium as follows:

1. Press  to navigate to 03 Date.
2. Press  to change the date.
 - ↳ The display changes to 03 Change date.
3. Press  repeatedly to choose between day, month and year.
4. Press  to activate the selected parameters (day, year).
 - ↳ The current value flashes.
5. Press  or  to reduce or increase the value.
6. Press  to save the new value.
 - ↳ The value no longer flashes.
7. Press  to change the month next. Proceed in a similar fashion here.
8. Press  to leave the setup menu.
 - ↳ The selected value is automatically saved when leaving the setup menu.







Setting the time

Working in the [SETUP menu](#) [▶ page 42] set the time as follows:

1. Press  repeatedly to navigate to 04 Time.
2. Press  in order to change the language.
↳ The display changes to 04 Change date.
3. Press  to activate the selected parameters.
↳ The current value flashes.
4. Press  or  to reduce or increase the value.
5. Press  to save the new value.
↳ The value no longer flashes.
6. Press  to leave the setup menu.
↳ The selected value is automatically saved when leaving the setup menu.





Setting the display contrast


Working in the [SETUP menu](#) [▶ page 42] set the display contrast as follows:

1. Press  repeatedly to navigate to 05 Display contrast.
2. Press  to activate the selected parameters.
↳ The current value flashes.
3. Press  or  to reduce or increase the value.
4. Press  to save the new value.
↳ The value no longer flashes.
5. Press  to leave the setup menu.
↳ The selected value is automatically saved when leaving the setup menu.

Selecting the language

You can choose between two languages. Language 0001 is usually the local language or English, Language 0002 is English or other required language Working in the [SETUP menu](#) [▶ page 42] set the language as follows:

1. Press  repeatedly to navigate to 06 Language.
2. Press  to activate the selected parameters.
↳ The current value flashes.
3. Press  to change to language 0002.
4. Press  to save the new value.
↳ The value no longer flashes.


5. Press  to leave the setup menu.


↳ The selected value is automatically saved when leaving the setup menu.

Other languages can also be installed. To this end, the corresponding language update file must be downloaded on the device from the CF card. Please consult your MELAG customer services or stockist for this.



Setting the water hardness


Working in the [SETUP menu](#) [▶ page 42] set the water hardness as follows:

1. Press  repeatedly to navigate to 07 Water °dH.


2. Press  to activate the selected parameters.

↳ The current value flashes.

3. Press  or  to reduce or increase the value.

4. Press  to save the new value.

↳ The value no longer flashes.

5. Press  to leave the setup menu.

↳ The selected value is automatically saved when leaving the setup menu.

You can find a water hardness conversion table in the chapter [Description of the device](#) [▶ page 13] under [Water softening unit](#) [▶ page 18].

9 Function tests

Automatic and manual functional check

Automatic

The device components are monitored and checked automatically for their functionality and interplay. Should the parameter thresholds be exceeded, the device will issue warnings or malfunction messages. If necessary, it will abort a program with the relevant notification. The device will also display messages when a program has been completed.


Manual

You can follow the program run on the display and use the log recorded to check the success of a program. Further information is provided in chapter [Logging](#) [▶ page 35].

Measuring conductivity

You can access the water quality of the DI water on the device on the display at any time providing, that it is switched on.

- ▶ Press  to start the "Conductivity measurement DI" program.



10 Maintenance



DANGER

All servicing work, especially that performed in the washing chamber may only be performed after a successful disinfection program.

- Comply with the working safety regulations.

Regular checks and cleaning



NOTICE

Incorrect cleaning can damage the surfaces and sealing faces. Scratched or damaged surfaces and leaking sealing faces favour soiling deposits and corrosion in the washing chamber.

- Comply with all information regarding cleaning of the parts affected.



NOTICE


When the coarse and fine sieves are missing, residue may enter the flushing circuit and impair the device function.

- Ensure that the coarse and fine sieves are always in place before program start.


Check the sieve in the washing chamber

A coarse and a fine sieve are fitted in the washing chamber. The sieves are designed to hold back dirt particles or soiling from the instruments. They can become blocked over time.

1. Inspect the coarse and the fine sieves for soiling and small components which have fallen from the load.
2. Turn the grip of the coarse sieve leftwards to its fullest extent and remove it upwards.



3. Turn the knurled nut on the fine sieve leftwards and remove the fine sieve upwards.




4. Inspect the coarse and the fine sieves for soiling.
5. Rinse the soiled sieve under running water. Do not use any dish-washing detergent. Remove any deposits with a soft brush.


Check the rinse arm

Soiling can block the nozzles of the rinse arms. Check both rinse arms regularly and rinse the nozzles if necessary.

1. Check that the coarse and the fine sieves are installed.
2. Turn the knurled nut on the rinse arm leftwards and remove the rinse arm.



3. Clean blocked nozzles with a thin pointed object.



4. Return the rinse arms and ensure their free movement.

Check the door seal

Check the door seal for impurities, deposits or damage on a daily basis. If necessary, clean the door seal with a moist, non-fuzzing cloth and conventionally-available neutral liquid cleaning fluid.

Check the injector rail nozzles for free passage

We recommend checking the injector rail nozzles for their free passage on a monthly basis.

To test whether the injector rail nozzles are blocked, they must be held upright under running water. If the water flows freely through the nozzles, they are not blocked.

Check the accessories

Check the accessories in particular with regard to the plastic parts (e.g. inserts) for damage, deposits and soiling on a monthly basis, unless the document *Instructions for the use and care of the accessories* indicates otherwise.

Cleaning on demand

Operating unit and plastic front

Comply with the following specifications for cleaning:

- ▶ Use a soft, non-fuzzing cloth.
- ▶ Use a chlorine- and vinegar-free cleaning fluid or a plastics cleaning agent.
- ▶ Check the material compatibility before application!
- ▶ Never use solvents or benzene.
- ▶ Use surface disinfectants which are suitable for plastics. Observe the manufacturer's information on the respective surface disinfectant.

Pump pit and non-return valve

If the rinse water has not been removed entirely after a program, the non-return valve must be cleaned.

1. Remove the coarse and fine sieves and remove the residue and deposits from the pump pit.
2. Remove the non-return valve upwards by pulling on its grip and pull it out of the pump pit.



3. Clean the non-return valve under running water. Do not use any dishwashing detergent.
4. Replace the non-return valve and the fine and coarse sieve in the pump pit.
5. Start the "Rinse" program.

Avoiding staining

Stains on the instruments or the device can develop from poor water quality. In particular, heavy metals or chloride deposits can result in the development of stains and rust. To avoid the development of stains on the instruments or the washing chamber, we recommend a final rinse with deionised water (DI water). All water-conducting parts of the device are made of non-rusting materials. This excludes the formation of stains or rust caused by the device. Often a rust-producing instrument already suffices for third-party rust to form on the other instruments or in the device. Further information is provided in the document *How to optimise cleaning and protect instruments*.

Changing the filter in the drying fan

Exceeding the permissible level of blockage can result in a worsened drying outcome. For this reason, the device checks the degree of blockage automatically. Exceeding the tolerances results in the issue of the relevant display message.



DANGER

Careless touching of the HEPA filter may damage the lamellae.

This may impair filter services and recontamination of the disinfected instruments.


- Only touch the HEPA filter by the frame and do not damage the lamellae.




PLEASE NOTE

For reasons of hygiene, the pre-filter and HEPA filter must be changed during maintenance.

1. Pull the process agent drawer forwards.



2. Undo the screw on the cover cap of the drying fan by hand and lift up the cover cap.



3. Pull out the pre-filter upwards and replace it.

4. Take out the HEPA filter upwards and replace it.



5. Close the cover cap and turn the screw hand-tight.

Maintenance



NOTICE

Continuing operation beyond the maintenance interval can result in malfunctions in the device.

- Maintenance should only be performed by trained and authorized customer services technicians, or stockist technicians.
 - Maintain the specified servicing intervals.
-

Regular maintenance is vital to ensure reliable operation and value retention of the device. All function and safety-relevant components and electrical units are checked during maintenance and replaced where necessary. The maintenance is performed by an authorized customer services / stockist technician working in accordance with the maintenance instructions. If the device is freely accessible, the maintenance takes around 3 hours plus test run and possible works not included in the regular maintenance plan.

Maintenance is to be performed after every 1000 cycles or after 18 months at the latest. Please note any restrictions arising from the process agents used (see [Approved process agents](#) [▶ page 11]).

(Process) validation

A reproducible cleaning and disinfection outcome can only be ensured via correct operation (inc. use of suitable accessories). The practice operator is responsible for ensuring reproducibility through the use of batch, routine and / or periodic inspections (e.g. validation).

This requirement is made (in Germany) by e.g. the *Medizinprodukte-Betreiberverordnung* (§ 8 Abs. 2 MPBetreibV); DGKH, DGSV and AKI directives and the recommendations from the Robert-Koch-Institut. This requirement is also made in international regulations. This is based on DIN EN ISO 15883, which is also valid in Germany.

Please comply with all valid national regulations and specifications. In case of doubt, please consult the relevant professional association.

- Only use the loading pattern specified and approved within the scope of the validation.
- We cannot provide a guarantee for non-MELAG accessories, even if they are in possession of validation.
- The validator and technical service can download a "Recommendation for the validation of MELAtherm 10" (doc.: AS_001-17_10DT_EN) from the MELAG service portal .

11 Pause times

Pause times between individual programs are not necessary. Instruments can be re-arranged and then cleaned and disinfected directly after a program run or interruption.

Run the "Rinse" program twice before treatment following pause times longer than two days e.g. following a weekend.

Given an ophthalmic application, run an empty batch (Ophthalmo-Program) twice before treatment following pause times longer than two days in order to flush the DI cartridge

Long operating pauses (longer than two weeks)

- Decommission the device if you plan to have a pause of over two weeks (see [Decommissioning](#) [▶ page 52]).

Decommissioning

Preparation for transport

Decommissioning in preparation for transport outside the practice should only be undertaken by MELAG-authorized individuals.

Following longer operating pauses



DANGER

Danger of acid burns by irritant substances!

Improper handling of the process agents may cause caustic burns and health damage.

- Protect your eyes, hands, clothing and all surfaces from contact with the process agents.
- Comply with the information from the manufacturer of the process agents.
- Every type of fluid (e.g. in the drawer, in the device floor tank or fluid emerging from the device) issued from the device as the result of damage could potentially contain aggressive process agents.

When decommissioning the device for a long pause (e.g. due to holiday), proceed as follows:

1. Remove the suction lances from the container and place them in a canister of water. The suction lances should be immersed to a minimum of 80%.
2. Run the "air-removal" program to free the metering system from process agents.
3. Return the suction lances to the process agents and screw them on tightly.
4. The interior of the washing chamber must be dry.
5. Switch off the device at the power switch.
6. Remove the mains plug from the socket after use.
7. Turn off the water inflow.

Recommissioning



NOTICE

Run the Air removal program twice before running the first decontamination process. Then start the usual disinfection program with an unloaded basis basket.

- ▶ Comply with the specifications in chapter [First steps](#) [▶ page 19] when performing the recommissioning.

Transport within the practice



CAUTION

Danger of injury from incorrect carrying.

Lifting and carrying too heavy a load can result in spinal injury. Failure to comply with these provisions can result in crushing.

- Please ensure that you comply with the transport instructions.
-

Comply with the following provisions pertaining to transport within the practice:

- ▶ Empty the device entirely
- ▶ Remove the inserts and the basis basket.
- ▶ Seal the water inflow hose
- ▶ Close the door before moving the device.
- ▶ Avoid strong vibrations.

Frost protection

Operate, store and transport the device in a generally frost-free environment. Should any residual fluids freeze in the device, the device should be held at room temperature for a minimum of two hours so that they can thaw.

Re-startup after change of locality

When recommissioning after a move, proceed as with the first commissioning (see chapter [First steps](#) [▶ page 19]).

12 Malfunctions

Warnings are marked in the display with a **W** and malfunction messages with an **M**. Follow all the instructions relating to a warning or malfunction message issued by the device display.

General events

General events serves to inform and provides assistance in the operation of the device. Malfunction-free operation of the device is still possible.

Warning

A warning helps to ensure malfunction-free operation and recognition of undesirable situations. React to a warning quickly to prevent the resulting malfunction.

Malfunction message

Malfunction messages are issued when it is not possible to ensure safe operation or cleaning and disinfection. These can appear on the display shortly after switching on or during a program run. If a malfunction occurs during a program run, the program will be aborted and will be classed as unsuccessful.



WARNING

Danger of infection from program abort.

If a program is interrupted before the start of a drying phase, the load is taken as not having been cleaned and disinfected. This endangers the health of the patient and the practice team.

- Only abort a current program in exceptional reasons.
- Never open the door after a program start.
- Treat the instruments again after a program abort.

General events



CAUTION

Danger of injury from insufficient protective measures!

Performing work without first taking the corresponding protective measures can result in injuries.

- Comply with the working safety measures required by the respective tasks.

The following table indicates possible causes for certain events and the corresponding operating information for their remedy. Should you be unable to find the relevant event, or your efforts do not redress the problem, you can contact your nearest stockist or authorized MELAG customer service provider. Please have your device serial number and a detailed description of the malfunction contained in the malfunction message to hand.

Incident	Possible cause	Remedy
Banging or rattling noise in the washing chamber during a program run.	The rinse arm bangs against instruments or containers. The load moves in the washing chamber.	Interrupt the program and re-arrange the load. Start the program again.

Incident	Possible cause	Remedy
White layer on the instruments	The internal water softening unit has not been programmed correctly.	Check the water hardness of the tap water and reset the internal water softening unit if necessary. See Description of the device [▶ page 13].
	Water-insoluble, hardened treatment residue (e.g. dental cement or root canal disinfectants) remain on the instruments.	Remove the treatment residue manually immediately after instrument application.
	Ultrasound gel residue can be found on the instruments.	Avoid the use of Quats-based (quaternary ammonium cation) cleaning fluids and disinfectants in the manual cleaning of instruments with gels or lubricants. Gels with thickening agents, especially polyacrylic acid separate after contact with Quats. Should you prefer to change gels, products with cation-compatible thickening systems are suitable. Consult the gel manufacturer or the process agents manufacturer for further information.
Poor cleaning outcome	The basis basket / insert basket / mount are incorrectly loaded or are too full.	Ensure correct arrangement and avoid overloading.
	Load results in a spray shadow.	Ensure the correct arrangement of the instruments.
	The cleaning agent is unsuitable for this type of soiling.	Use a suitable cleaning agent for automatic cleaning
	Encrusted soiling on the instruments.	Do not allow soiling to dry on. Rinse off soiling immediately.
	Rinse arm nozzles or injector rail nozzles blocked.	Rectify the blockage in accordance with the description provided in chapter Maintenance [▶ page 47].
	Pump pit sieves are soiled.	Clean the coarse and fine sieves in accordance with with the description provided in Maintenance [▶ page 47].
Empty display	The device is not switched on.	Check that the device is connected to the power supply and is switched on.
	The fuse in domestic installation has tripped. This can be caused by operating a number of electrical devices at the same time.	Check the house fuse (see type plate for the min. fuse).
Display: Salt exhausted, refill! Then wait for signal	There regeneration salt is exhausted.	Fill the salt container with regenerating salt. The signal (a tone) informs the operator that the salt in the salt container has been recognized and that operation can be continued.
Residual moisture on and / or on the instruments.	The basis basket / insert basket / mount are incorrectly loaded or are too full.	Ensure correct arrangement and avoid overloading.
	The interior structure of the instruments is too complex / the interior volume is insufficient.	Dry the instruments with clean (medical) compressed air.

Warnings

Warning	Possible cause	Remedy
W214	The CF card was removed from the slot during a running program and re-inserted.	Once the program has been completed, select DOCU MENU in the display and output the current log. Do not remove the CF card during active logging. Logging is active when the red LED illuminates.
W215 W216 W217	The CF card does not function correctly. The system does not recognize a CF card or is unable to read it.	1. Save the log on an external data carrier. 2. Working in the display, select DOCU MENU and navigate to 15Format CF card . Format the CF card in the device.
	The CF card memory is too large (max. 4 GB).	Insert a CF card with a memory of up to 4 GB. The MELAG logo must point rightwards during insertion. MELAG recommends using original accessories only.
W218	An already-existing log has been recognized on the CF card whilst outputting the log via the DOCU MENU.	Acknowledge the message with the button 4. The existing log will not be overwritten.
W219 W220	The CF card does not function correctly. The system does not recognize a CF card or is unable to read it.	1. Save the log on an external data carrier. 2. Working in the display, select DOCU MENU and navigate to 15Format CF card . Format the CF card in the device.
	The CF card memory is too large (max. 4 GB).	Insert a CF card with a memory of up to 4 GB. The MELAG logo must point rightwards during insertion. MELAG recommends using original accessories only.
W221	The CF card memory is full. No further logs can be saved.	1. Save the log on an external data carrier. 2. Working in the display, select DOCU MENU and navigate to 15Format CF card . Format the CF card in the device.
W222 W223 W224 W225 W226 W227	The CF card does not function correctly. The system does not recognize a CF card or is unable to read it.	1. Save the log on an external data carrier. 2. Working in the display, select DOCU MENU and navigate to 15Format CF card . Format the CF card in the device.
	The CF card memory is too large (max. 4 GB).	Insert a CF card with a memory of up to 4 GB. The MELAG logo must point rightwards during insertion. MELAG recommends using original accessories only.
W228	The CF card is too slow. Either the CF card is no longer recognized following a reset or it was inserted in the slot under voltage.	1. Save the log on an external data carrier. 2. Insert a new CF card (max. 4 GB) in the card slot. The MELAG logo must point rightwards during insertion. MELAG recommends using original accessories only.
W229	The CF card was removed from the slot during a writing action.	Once the program has been completed, select DOCU MENU in the display and output the current log. Do not remove the CF card during active logging. Logging is active when the red LED illuminates.

Warning	Possible cause	Remedy
W230	The CF card does not function correctly. The system does not recognize a CF card or is unable to read it.	<ol style="list-style-type: none"> 1. Save the log on an external data carrier. 2. Working in the display, select DOCU MENU and navigate to 15Format CF card. Format the CF card in the device.
	The CF card memory is too large (max. 4 GB).	<p>Insert a CF card with a memory of up to 4 GB. The MELAG logo must point rightwards during insertion.</p> <p>MELAG recommends using original accessories only.</p>
W231	The CF card does not function correctly. There is no CF card in the slot.	<p>Insert a CF card with a memory of up to 4 GB. The MELAG logo must point rightwards during insertion.</p> <p>MELAG recommends using original accessories only.</p>
	The system does not recognize a CF card or the CF card cannot be read.	Push the CF card in the card slot until it the ejector button triggers.
W232	The CF card does not function correctly.	Acknowledge the message with the button 4.
W233	The CF card is currently being initialized or written.	
W234 W235 W236 W237	The CF card does not function correctly. The system does not recognize a CF card or is unable to read it.	<ol style="list-style-type: none"> 1. Save the log on an external data carrier. 2. Working in the display, select DOCU MENU and navigate to 15Format CF card. Format the CF card in the device.
	The CF card memory is too large (max. 4 GB).	<p>Insert a CF card with a memory of up to 4 GB. The MELAG logo must point rightwards during insertion.</p> <p>MELAG recommends using original accessories only.</p>
W238	The CF card is not functioning correctly and cannot be formatted.	<p>Insert a new CF card (max. 4 GB) in the card slot. The MELAG logo must point rightwards during insertion.</p> <p>MELAG recommends using original accessories only.</p>
W239 W240	The CF card does not function correctly. The system does not recognize a CF card or is unable to read it.	<ol style="list-style-type: none"> 1. Save the log on an external data carrier. 2. Working in the display, select DOCU MENU and navigate to 15Format CF card. Format the CF card in the device.
	The CF card memory is too large (max. 4 GB).	<p>Insert a CF card with a memory of up to 4 GB. The MELAG logo must point rightwards during insertion.</p> <p>MELAG recommends using original accessories only.</p>
W372	The internal log memory of the device is full. Not all logs have been outputted.	<ol style="list-style-type: none"> 1. Select DOCU MENU in the display and output the logs of the internal memory. 2. Start the program again. 3. If this notification is displayed repeatedly, delete the internal memory.

Warning	Possible cause	Remedy
W377	The system does not recognize an output medium. The system does not recognize a log printer, even though it is connected.	Check the settings in SETUP MENU → 02Autom. logging.
	Automatic logging is active in the SETUP MENU. However, a log printer is not connected.	<ol style="list-style-type: none"> Working in the display, select DOCU MENU and save the logs on the CF card or on a PC. Working in SETUP MENU deactivate → 02Autom. logging. The display changes from ACTIVE to INACTIVE.
W386	The internal log memory contains logs which have yet to be outputted. The memory is almost full.	Acknowledge the message with the button 4. The program starts. As soon as the program has ended, working in the display, select the DOCU MENU and output all logs from the internal memory (CF card or external data carrier).
W394	Not all logs from the internal device memory have been saved on the CF card.	Acknowledge the message with the button 4. The logs are written and saved on the CF card.
W395	Not all logs have been outputted from the internal device log memory via the EDM printer.	Acknowledge the message with the button 4. The logs are outputted and printed.
W396	Not all logs have been loaded onto the FTP server from the internal log memory of the device.	Acknowledge the message with the button 4. The logs are outputted and saved.
W397	The system is unable to locate a PC for log output. Even though the MELAG device is connected to a PC, it is unable to establish a connection. <ol style="list-style-type: none"> Even though the device is connected to a PC, it is unable to establish a connection for log output. The device is not connected to a PC, but the PC option has been activated in SETUP MENU. 	<ol style="list-style-type: none"> Check the network connection to the PC/server. Switch on the PC/server. Restart the documentation software.
	The MELAG device is not connected to a PC, but the option "Computer" is active in SETUP MENU → 02Autom. logging.	Working on the display, select SETUP MENU → 02Autom. logging and deactivate the option "Computer". The display changes from YES to NO.
W414	The rinse aid container has been exhausted.	<ol style="list-style-type: none"> Replace the rinse aid container, working in accordance with the working safety regulations. Alternatively, refill. Start the "Air removal" program. NOTICE! Use only process agents which you have used before.
W424	The neutralizer container has been exhausted.	<ol style="list-style-type: none"> Replace the neutralizer container, working in accordance with the working safety regulations. Alternatively, refill. Start the "Air removal" program. NOTICE! Use only process agents which you have used before.

Warning	Possible cause	Remedy
W425	The cleaning agent container has been exhausted.	<ol style="list-style-type: none"> 1. Replace the cleaning agent container, working in accordance with the working safety regulations. Alternatively, refill. 2. Start the "Air removal" program. <p>NOTICE! Use only process agents which you have used before.</p>
W428	There is almost no regeneration salt left	Refill regenerating salt, see Filling the regeneration salt [▶ page 21]
W447	The rinse pressure in the washing chamber is too low. Large containers with the opening pointing upwards may be sorted into the device. This diverts water from the rinsing process.	Setup the containers with their openings facing downwards.
W450	The water inflow is insufficient.	Check the water supply of the device. Open the water inflow tap completely.
W475	The HEPA filter in the drying fan is soiled.	Replace the HEPA-filter and drying fan, see Changing the filter in the drying fan [▶ page 50].
W477	The requisite pressure for the drying has not been achieved. The pre-filter in the drying fan is soiled.	Replace the drying fan pre-filter, see Changing the filter in the drying fan [▶ page 50].
	The lid of the drying fan has not been locked correctly.	Lock the lid of the drying fan correctly.
W478	The HEPA filter and the pre-filter in the drying fan are soiled.	Replace the HEPA-filter and pre-filter, see Changing the filter in the drying fan [▶ page 50].
W500	The display of date and time of the system clock are incorrect.	Working in the display, select SETUP MENU and set the correct date and time, see Setting date and time [▶ page 43].
W501	The CF card does not function correctly. There is no CF card in the slot.	Insert a CF card with a memory of up to 4 GB. The MELAG logo must point rightwards during insertion. MELAG recommends using original accessories only.
	The system does not recognize a CF card or the CF card cannot be read.	Push the CF card in the card slot until it the ejector button triggers.
W502	The system is unable to locate a PC for log output. The network connection has been interrupted.	Check the network connection to the PC/ server.
	The PC/server is not switched on.	Switch on the PC/server.
	The documentation software has not been started.	Restart the documentation software.
	A PC is not connected but the option "Computer" is active in SETUP MENU → 02Autom. logging.	Working on the display, select SETUP MENU → 02Autom. logging and deactivate the option "Computer". The display changes from YES to NO.
W533	The temperature in the washing chamber is very high. The door is blocked and cannot be unlocked immediately.	<p>CAUTION! The instruments are hot.</p> <p>Press the keys indicated in the display to acknowledge the notification. The door can be opened.</p> <p>PLEASE NOTE Take suitable safety precautions such as standing at a safe distance and wearing heat-proof gloves.</p>

Warning	Possible cause	Remedy
W534	The temperature in the washing chamber is very high. The door is blocked and cannot be unlocked immediately.	DANGER! Danger of scalding! The instruments are hot. 1. Wait until the temperature of the washing chamber has cooled. 2. Press the keys indicated in the display.
W549	The conductivity of the DI water is insufficient (greater than 15 μ S). The MELAdem 53 cartridge is exhausted.	Replace the MELAdem 53 cartridge.
	The quality of the DI water supply is not of sufficient quality.	Check the DI water supply.
W560	The maximum permissible mains voltage (270 V) has been exceeded.	Have the connection conditions checked by an electrician.
W561	The minimum permissible mains voltage (190 V) was undercut.	Have the connection conditions checked by an electrician.
W562	The maximum permissible mains frequency (63 Hz) was exceeded.	Have the connection conditions checked by an electrician.
W563	The minimum permissible mains frequency (45 Hz) was undercut.	Have the connection conditions checked by an electrician.
W575	The date and time are invalid.	Check the settings in SETUP MENU.
W622	The maximum permissible maintenance period or the maximum permissible number of cycles (1000 cycles) has been achieved since commissioning or the last maintenance.	Arrange for maintenance with an authorized customer services or a stockist technician. You can continue to start the device.
W625	The temperature during pre-cleaning is too high. The temperature during the water inflow is higher than 45 °C.	Check the water supply to the device.
W671	An insufficient conductivity value (>15 μ S und <25 μ S) was measured in the washing chamber during disinfection in the Ophthamo-Program. This could be caused by carry-over of process agents, regenerating salt or deposits. The program successfully completed despite the warning.	1. Close the lid of the salt container correctly. 2. Setup the containers in the device with their openings facing downwards. 3. Check the hollow bodies before decontamination for their free passage and correct position. 4. Clean the filter screen in the instrument connection equipment. 5. Remove and clean the coarse and fine sieve, see Regular checks and cleaning [▶ page 47]. 6. Insert the non-return valve in the pump pit correctly, see Cleaning on demand [▶ page 49]. 7. Check whether foreign bodies are in the non-return valve.

Malfunction messages

Malfunction	Possible cause	Remedy
F137	The cleaning agent metering pump is not functioning correctly. The metering system may be blocked.	<ol style="list-style-type: none"> 1. Switch off the device off and then on again. 2. Start the program again.
F139	The fan of the display is not functioning correctly.	<ol style="list-style-type: none"> 1. Switch off the device off and then on again. 2. Start the program again.
F140	The fan of the diffuser is not functioning correctly.	<ol style="list-style-type: none"> 1. Switch off the device off and then on again. 2. Start the program again.
F141	The neutralizer metering pump is not functioning correctly. The metering system may be blocked.	<ol style="list-style-type: none"> 1. Switch off the device off and then on again. 2. Start the program again.
F142	The rinse aid metering pump is not functioning correctly. The metering system may be blocked.	<ol style="list-style-type: none"> 1. Switch off the device off and then on again. 2. Start the program again.
F143	The solenoid valve for the cold water does not switch.	<ol style="list-style-type: none"> 1. Switch off the device off and then on again. 2. Start the program again.
F144	The solenoid valve for the regeneration does not switch.	<ol style="list-style-type: none"> 1. Switch off the device off and then on again. 2. Start the program again.
F145	The solenoid valve for the steam condenser does not switch.	<ol style="list-style-type: none"> 1. Switch off the device off and then on again. 2. Start the program again.
F146	The solenoid valve of the DI inflow hose does not switch.	<ol style="list-style-type: none"> 1. Switch off the device off and then on again. 2. Start the program again.
F147	The solenoid valve of the cold water inflow hose does not switch.	<ol style="list-style-type: none"> 1. Switch off the device off and then on again. 2. Start the program again.
F154	The temperature difference between the two temperature sensors (temperature control and temperature log) in the washing chamber is too high.	<ol style="list-style-type: none"> 1. Switch off the device and wait c. 30 minutes with the door open. 2. Switch on the device and re-start the program.
F155	The temperature difference between the two temperature sensors (temperature control and temperature log) in the washing chamber is too high.	<ol style="list-style-type: none"> 1. Switch off the device and wait c. 30 minutes with the door open. 2. Switch on the device and re-start the program.
F156	The temperature sensor for monitoring the drying is not functioning correctly.	<ol style="list-style-type: none"> 1. Switch off the device and wait c. 30 minutes with the door open. 2. Switch on the device and re-start the program.
F159	The collection tank pump has not been emptied correctly.	<ol style="list-style-type: none"> 1. Switch off the device off and then on again. 2. Start the program again.

Malfunction	Possible cause	Remedy
F160	The coarse or fine sieve is soiled.	<ol style="list-style-type: none"> 1. Switch off the device.. 2. Clean the coarse and fine sieve, see Regular checks and cleaning [▶ page 47]. 3. Switch on the device and re-start the program.
F161	The washing chamber pressure required for drying has not been reached.	<ol style="list-style-type: none"> 1. Switch off the device off and then on again. 2. Start the program again.
F162	The requisite rinse pressure has not been reached.	<ol style="list-style-type: none"> 1. Switch off the device off and then on again. 2. Start the program again.
F163	The cleaning agent metering pump is not functioning correctly. The metering system may be blocked.	<ol style="list-style-type: none"> 1. Switch off the device off and then on again. 2. Start the program again.
F165	The fan of the display is not functioning correctly.	<ol style="list-style-type: none"> 1. Switch off the device off and then on again. 2. Start the program again.
F166	The fan of the diffuser is not functioning correctly.	<ol style="list-style-type: none"> 1. Switch off the device off and then on again. 2. Start the program again.
F167	The neutralizer metering pump is not functioning correctly. The metering system may be blocked.	<ol style="list-style-type: none"> 1. Switch off the device off and then on again. 2. Start the program again.
F168	The rinse aid metering pump is not functioning correctly. The metering system may be blocked.	<ol style="list-style-type: none"> 1. Switch off the device off and then on again. 2. Start the program again.
F169	The solenoid valve for the cold water does not switch.	<ol style="list-style-type: none"> 1. Switch off the device off and then on again. 2. Start the program again.
F170	The solenoid valve for the regeneration does not switch.	<ol style="list-style-type: none"> 1. Switch off the device off and then on again. 2. Start the program again.
F171	The solenoid valve for the steam condenser does not switch.	<ol style="list-style-type: none"> 1. Switch off the device off and then on again. 2. Start the program again.
F172	The solenoid valve of the DI inflow hose does not switch.	<ol style="list-style-type: none"> 1. Switch off the device off and then on again. 2. Start the program again.
F173	The solenoid valve of the cold water inflow hose does not switch.	<ol style="list-style-type: none"> 1. Switch off the device off and then on again. 2. Start the program again.
F 410	The rinse aid container has been exhausted.	<ol style="list-style-type: none"> 1. Replace the rinse aid container, working in accordance with the working safety regulations. Alternatively, refill. NOTICE! Use only process agents which you have used before. 2. Start the "Air removal" program.

Malfunction	Possible cause	Remedy
F411	The neutralizer container has been exhausted.	<ol style="list-style-type: none"> 1. Replace the rinse aid container, working in accordance with the working safety regulations. Alternatively, refill. NOTICE! Use only process agents which you have used before. 2. Start the "Air removal" program.
F412	The cleaning agent container has been exhausted.	<ol style="list-style-type: none"> 1. Replace the rinse aid container, working in accordance with the working safety regulations. Alternatively, refill. NOTICE! Use only process agents which you have used before. 2. Start the "Air removal" program.
F426	No cleaning agent is pumped. The cleaning agent container has been exhausted or air is being transported.	<ol style="list-style-type: none"> 1. Replace the cleaning agent container, working in accordance with the working safety regulations. Alternatively, refill. NOTICE! Use only process agents which you have used before. 2. Start the "Air removal" program.
	The hose to the suction lance is kinked.	<ol style="list-style-type: none"> 1. Rectify the kinks or crushing on the metering hoses. 2. Start the "Air removal" program.
	Air bubbles have developed in the metering system after long immobilization times.	Start the "Air removal" program.
F427	No neutralizer is pumped. The neutralizer container has been exhausted or air is being transported.	<ol style="list-style-type: none"> 1. Replace the neutralizer container, working in accordance with the working safety regulations. Alternatively, refill. NOTICE! Use only process agents which you have used before. 2. Start the "Air removal" program.
	The hose to the suction lance is kinked.	<ol style="list-style-type: none"> 1. Rectify the kinks or crushing on the metering hoses. 2. Start the "Air removal" program.
	Air bubbles have developed in the metering system after long immobilization times.	Start the "Air removal" program.
F431	No cleaning agent is pumped. The cleaning agent container has been exhausted or is nearly empty.	<ol style="list-style-type: none"> 1. Replace the cleaning agent container, working in accordance with the working safety regulations. Alternatively, refill. NOTICE! Use only process agents which you have used before. 2. Start the "Air removal" program.
	The hose to the suction lance is kinked.	<ol style="list-style-type: none"> 1. Rectify the kinks or crushing on the metering hoses. 2. Start the "Air removal" program.
	Air bubbles have developed in the metering system after long immobilization times.	Start the "Air removal" program.

Malfunction	Possible cause	Remedy
F432	No neutralizer is pumped. The neutralizer container has been exhausted or is nearly empty.	1. Replace the neutralizer container, working in accordance with the working safety regulations. Alternatively, refill. NOTICE! Use only process agents which you have used before. 2. Start the "Air removal" program.
	The hose to the suction lance is kinked.	1. Rectify the kinks or crushing on the metering hoses. 2. Start the "Air removal" program.
	Air bubbles have developed in the metering system after long immobilization times.	Start the "Air removal" program.
F433	Water is in the pump pit after pumping out. The coarse or fine sieve is soiled.	Clean the coarse and fine sieve, see Regular checks and cleaning [▶ page 47].
	The non-return valve in the pump pit is missing or fitted incorrectly.	Insert the non-return valve in the pump pit correctly, see Regular checks and cleaning [▶ page 47]
	The non-return valve is blocked by a foreign body.	Check the non-return valve for foreign bodies and remove them if you find any.
F434	Water is in the pump pit after pumping out. The coarse or fine sieve is soiled.	Clean the coarse and fine sieve, see Regular checks and cleaning [▶ page 47].
	The non-return valve in the pump pit is missing or fitted incorrectly.	Insert the non-return valve in the pump pit correctly, see Regular checks and cleaning [▶ page 47].
	The non-return valve is blocked by a foreign body.	Check the non-return valve for foreign bodies and remove them if you find any.
	The effluent hose is kinked.	Check the installation of the waste water hose.
	The outflow or effluent hose is blocked.	Check the siphon and the effluent hose for blockage.
F440	The current program has been terminated prematurely. The load is classed as not cleaned and disinfected.	1. Acknowledge the message with the button 4. 2. Press the keys indicated in the display.

Malfunction	Possible cause	Remedy
F449	The rinse pressure in the washing chamber is too high. The water inflow is insufficient.	Check the water inflow of the device. Open the water inflow tap completely.
	The basis basket has been inserted incorrectly or not at all.	Insert the basis basket in the washing chamber correctly, see Inserting the basis basket [▶ page 21].
	Too many non-filled apertures on the injector rail.	Seal non-filled apertures on the injector rail with a sealing screw.
	The coarse or fine sieve is soiled.	Remove and clean the coarse and fine sieve, see Regular checks and cleaning [▶ page 47].
	Large containers with the opening pointing upwards may be sorted into the device. This diverts water from the rinsing process.	Sort the containers in the device with their openings facing downwards.
	Strong foam generation: The instruments have been pre-cleaned in or placed in a foam-generating solution and have then been subject to insufficient rinsing.	Rinse the instruments thoroughly before decontamination.
	Strong foam generation: Strong soiling of the filter disc in the universal adapter for transfer instruments.	Remove and replace the soiled filter disc. Clean the reusable filter screen.
	Strong foam generation: Non-approved process agents (rinse aid or cleaning agent) have been used.	NOTICE! Use only MELAG-approved process agents.
F451	The temperature difference between the two temperature sensors in the washing chamber is too great. The temperature sensors were not covered with water sufficiently. The upper rinse arm revolves too slowly.	Clean the upper rinse arm and check its ease of movement.
F462	The water inflow is insufficient. The water inflow tap has not been opened completely.	Check the water inflow of the device. Open the water inflow tap completely.
	The sieve in the cold water connection is blocked.	Remove and clean the sieve in the cold-water connection.
	The cold water inflow hose is kinked.	Check the installation of the cold water inflow hose.
F464	The water inflow is insufficient. The water inflow tap has not been opened completely.	Check the water inflow of the device. Open the water inflow tap completely.
	The sieve in the cold water connection is blocked.	Remove and clean the sieve in the cold-water connection.
	The cold water inflow hose is kinked.	Check the installation of the cold water inflow hose.
F466	Insufficient DI water inflow. The DI water supply has been interrupted.	Check the DI water supply. Check the DI water system for its correct function.
	The sieve in the DI water connection is blocked.	Remove and clean the sieve in the DI water connection.
	The DI water inflow hose is kinked.	Check the installation of the DI water inflow hose.

Malfunction	Possible cause	Remedy
F467	The water inflow is insufficient. The water inflow tap has not been opened completely.	Check the water inflow of the device. Open the water inflow tap completely.
	The sieve in the cold water connection is blocked.	Remove and clean the sieve in the cold-water connection.
	The cold water inflow hose is kinked.	Check the installation of the cold water inflow hose.
F468	Insufficient DI water inflow. The DI water supply has been interrupted.	Check the DI water supply. Check the DI water system for its correct function.
	The sieve in the DI water connection is blocked.	Remove and clean the sieve in the DI water connection.
	The DI water inflow hose is kinked.	Check the installation of the DI water inflow hose.
F471	The notification is triggered through a faults operating sequence in the DIAGNOSIS +SERVICE menu.	Switch off the device off and then on again.
F474	The HEPA filter is not recognized. A HEPA filter has not been inserted	Insert the HEPA filter correctly.
	The HEPA filter for the drying fan has not been inserted correctly.	Check whether the HEPA filter for the drying fan has been inserted correctly.
	The lid of the drying fan has not been locked correctly.	Close the lid of the drying fan correctly.
F476	The requisite pressure for the drying has not been achieved. The HEPA filter for the drying fan has not been inserted correctly.	Check whether the HEPA filter for the drying fan has been inserted correctly.
	The lid of the drying fan has not been locked correctly.	Close the lid of the drying fan correctly.
F484	The rinse pressure in the washing chamber is too high. The water inflow is insufficient.	Check the water inflow of the device. Open the water inflow tap completely.
	The basis basket has been inserted incorrectly or not at all.	Insert the basis basket in the washing chamber correctly. The injector rail should be located on the right-hand side and dock with the blind cap on the port of the rear wall, see Inserting the basis basket [▶ page 21].
	Too many non-filled apertures on the injector rail.	Seal non-filled apertures on the injector rail with a sealing screw.
	The coarse or fine sieve is soiled.	Remove and clean the coarse and fine sieve, see Regular checks and cleaning [▶ page 47].
	Large containers with the opening pointing upwards may be sorted into the device. This diverts water from the rinsing process.	Sort the containers in the device with their openings facing downwards.
	Strong foam generation: The instruments have been pre-cleaned with a foam-generating solution and have then been subject to insufficient rinsing.	Rinse the instruments thoroughly before decontamination.
	Strong foam generation: Strong soiling of the filter disc in the universal adapter for transfer instruments.	Remove and replace the soiled filter disc. Clean the reusable filter screen.
	Strong foam generation: Non-approved process agents have been used.	NOTICE! Use only MELAG-approved process agents.

Malfunction	Possible cause	Remedy
F505	The salt store has been exhausted. No new regeneration can be performed.	Refill regenerating salt, see Filling the regeneration salt [▶ page 21] A program can be started if the salt has dissolved in the water. Do not start the program until the regenerating salt has been filled and the signal tone has sounded.
F509	Liquid in the device floor trough.	<p>CAUTION! Avoid contact with liquids in the floor trough; they can contain process agents.</p> <ol style="list-style-type: none"> 1. Switch off the device.. 2. Close the water inflow tap 3. Contact the authorized customer services or stockist technician, stating the device serial number and the malfunction number displayed.
F510	Too high a water level was measured in the washing chamber during a program run.	<ol style="list-style-type: none"> 1. Press the keys indicated in the display. 2. Close the door and start the program again.
F512	The running program was interrupted by a power cut	<p>NOTICE! The load of the interrupted program has not been cleaned or disinfected.</p> <ol style="list-style-type: none"> 1. Acknowledge the message with the button 4. 2. Start the program again.
F524	The door of the device is blocked and cannot be closed correctly.	Check the door area for blockages.
F531	The emergency release on the door was actuated during a program run.	<p>NOTICE! The load of the interrupted program has not been cleaned or disinfected.</p> <ol style="list-style-type: none"> 1. Acknowledge the message with the button 4. 2. Close and lock the door correctly. 3. Start the program again.
F535	The fine sieve has been fitted incorrectly.	Insert the fine sieve correctly. The arrow on the fine sieve must point towards the left-hand corner of the washing chamber.
F536 F537	The upper/lower rinse arm is mechanically blocked.	Check the freedom of motion of the upper/lower rinse arm.
	The impulse nozzle of the upper/lower rinse arm is blocked.	Remove and clean the upper/lower rinse arm.
	The basis basket has been inserted in the incorrect position or not at all.	Insert the basis basket correctly. The injector rail must dock on to the connection fitting.
	Fine deposits on the rinse arm bearing above or on the sliding disc.	Remove and clean the upper/lower rinse arm. Clean the sliding disc with a cloth.
	The water inflow is not sufficient.	<p>Check the water inflow to the device:</p> <ol style="list-style-type: none"> 1. Remove and clean the sieve in the cold-water connection. 2. Check the installation of the inflow hose. 3. Open the water inflow tap completely.

Malfunction	Possible cause	Remedy
F538 F539	The upper/lower rinse arm is mechanically blocked.	Check the freedom of motion of the upper/lower rinse arm.
	The impulse nozzle of the upper/lower rinse arm is blocked.	Remove and clean the upper/lower rinse arm.
	The basis basket has been inserted in the incorrect position or not at all.	Insert the basis basket correctly. The injector rail must dock on to the connection fitting.
	Fine deposits on the rinse arm bearing above or on the sliding disc.	Remove and clean the upper/lower rinse arm. Clean the sliding disc with a cloth.
	The water inflow is not sufficient.	Check the water inflow to the device: <ol style="list-style-type: none"> 1. Remove and clean the sieve in the cold-water connection. 2. Check the installation of the inflow hose. 3. Open the water inflow tap completely.
F546	The MELAdem 53 cartridge has not been bled correctly. A sudden water flow causes short-term incorrect measured values.	<ol style="list-style-type: none"> 1. Bleed the MELAdem 53 cartridge of air see the separate instructions "Removing the air from the MELAdem 53/MELAdem 53 C". 2. Start the program again.
F548	The conductivity of the DI water is insufficient (greater than 60 μ S). The MELAdem 53 cartridge is exhausted.	Replace the MELAdem 53 cartridge.
	The quality of the DI water supply is not of sufficient quality.	Check the DI water supply.
F571	The program cannot be started as brine is still in the water softening unit or washing chamber. Only the program "Regeneration" may be started.	Start the "Regeneration" program.
F583	The water inflow was interrupted during the active program.	<ol style="list-style-type: none"> 1. Open the water inflow tap completely. 2. Start the program again. <p>The water inflow must be secured during the entire duration of the active program.</p>
F620	Strong foam development in the washing chamber. The instruments are pre-cleaned or placed in a foam-generating solution.	Load the instruments into the MELAtherm without pre-treatment or rinse them thoroughly after loading.
	Non-approved process agents (incorrect rinse aid or cleaning agent) have been used.	Use only MELAG-approved process agents.
	The metering concentration has been set incorrectly.	Check the settings for the metering concentration and if necessary, arrange for correction by an authorized customer service or a stockist technician.
	Strong soiling of the filters in the transfer instrument adapter.	Clean or renew the filter in regular intervals.
F624	The collection tank is not pumped out.	<ol style="list-style-type: none"> 1. Switch off the device off and then on again. 2. Start the program again.
F626	The temperature during pre-cleaning is too high.	Check the water inflow to the device.

Malfunction	Possible cause	Remedy
F632	The coarse or fine sieve is soiled.	<ol style="list-style-type: none"> 1. Remove and clean the coarse and fine sieve, see Regular checks and cleaning [▶ page 47]. 2. Switch off the device off and then on again. 3. Start the program again.
F653	The water inflow was interrupted during the active program.	<ol style="list-style-type: none"> 1. Open the water inflow tap completely. 2. Start the program again. <p>The water inflow must be secured during the entire duration of the active program.</p>
F660 F661	The power supply for the MELAtherm 10 DTA is insufficient.	<ol style="list-style-type: none"> 1. Check whether the power plug has been inserted correctly in the socket. 2. Check the fuses in the sub-distribution.
F662	The upper rinse arm is soiled.	Remove the upper/lower rinse arm and clean the nozzles, see Regular checks and cleaning [▶ page 47].
F669	The coarse or fine sieve is strongly soiled.	<ol style="list-style-type: none"> 1. Remove and clean the coarse and fine sieve, see Regular checks and cleaning [▶ page 47]. 2. Switch off the device off and then on again. 3. Start the program again.
F670	The water inflow was interrupted during the active program.	<ol style="list-style-type: none"> 1. Open the water inflow tap completely. 2. Start the program again. <p>The water inflow must be secured during the entire duration of the active program.</p>
F672	<p>An insufficient conductivity value ($\geq 25 \mu\text{S}$) was measured in the washing chamber during disinfection in the Ophthalmic-Program.</p> <p>This could be caused by carry-over of process agents, regenerating salt or deposits. The program successfully completed despite the warning.</p>	<ol style="list-style-type: none"> 1. Close the lid of the salt container correctly. 2. Setup the containers in the device with their openings facing downwards. 3. Check the hollow bodies before decontamination for their free passage and correct position. 4. Clean the filter screen in the instrument connection equipment. 5. Remove and clean the coarse and fine sieve, see Regular checks and cleaning [▶ page 47]. 6. Insert the non-return valve in the pump pit correctly, see Cleaning on demand [▶ page 49]. 7. Check whether foreign bodies are in the non-return valve.
F673	The Ophthalmic-Program does not start. A DI connection is not set in SETUP MENU.	<ol style="list-style-type: none"> 1. Connect the DI water. 2. Working on the display, select SETUP MENU → 01DI water and set the parameters to YES.

Malfunction	Possible cause	Remedy
F675	Water is in the pump pit after pumping out. The coarse or fine sieve is soiled.	Remove and clean the coarse and fine sieve, see Regular checks and cleaning [▶ page 47].
	The non-return valve in the pump pit is missing or fitted incorrectly.	Insert the non-return valve in the pump pit correctly.
	The non-return valve is blocked by foreign bodies.	Check whether foreign bodies are in the non-return valve, see Cleaning on demand [▶ page 49].

13 Technical data

Table 7: MELAtherm 10 DTA/DTB device dimensions

Device type	Free standing	Under-desk unit	Cabinet device
Dimensions (HxWxD) ²	83.6 x 59.8 x 67.83 cm	81.8 x 59.8 x 67.83 cm	124 x 59.8 x 67.83 cm
Empty weight	85 kg	79 kg	106 kg
Operating weight	113 kg	119 kg	182 kg
Max. set-up height	1500 m (it may be necessary to reduce the disinfection temperature depending on the installation height. Consult technical manual)		

Device type	MELAtherm 10 DTA	MELAtherm 10 DTB
Electrical connection	3N AC 380-415V; 50/60 Hz; 3x16 A; 9.3 kW ³	3N AC 220-240V; 50/60 Hz; 3x16 A; 3.3 kW ⁴
Washing chamber (HxWxD)	29 x 45.5 x 42.3 cm	
Volume of the washing chamber	84 l	
Max. load of the door	15 kg	
Acoustic power	Median value 68 dB(A), max.73 dB(A)	
Waste heat	0.75 kWh (2.7 MJ)	
Ambient temperature	5-40 °C (recommended max. 25 °C)	
Relative humidity	max. 80% at 31 °C, decreasing in a linear fashion up to a relative humidity of 50% at 40 °C	
Installation category	2	
Air pressure	75 kPa –106 kPa	
Connection CW / DI water	3/4" internal thread (for the connection to a standard 3/4" connection with external thread)	
Effluent water connection	DN21	
Water quality	Drinking water according to Drinking Water Ordinance (TrinkW2001) / observe local specifications	
Minimum flow pressure	150 kPa (1.5 bar)	
Recommended flow pressure	250 kPa (2.5 bar)	
Max. water pressure	1000 kPa (10 bar)	
Max effluent temperature	93 °C (<1 min, ca. 5.5 l)	
Cold water temperature	1-26 °C	
Amount of effluent water per hour	ca. 30 l (in small intervals)	
Capacity of drain pump	max. 40 l/min. (volume in effluent hose)	
Length of the inflow hoses and outlet hose	each 1.80 m (extension optionally available)	
Length of power cable	2 m	
Degree of soiling	Category 2	
Degree of protection (following IEC 60529)	IP20	
CE mark	CE 0197	

²⁾ Appropriate for a 60-cm deep working surface

³⁾ Comply with the maximum voltage range of 360-440V

⁴⁾ Comply with the maximum voltage range of 207-253V

14 Accessories and Spare Parts

Accessories

You can obtain the specified articles together with an overview of further accessories from your stockist. Information regarding the instrument decontamination accessories can be found in the current MELAG price list.

	Article	Article no.
Optionally available	Floor unit (HxWxD 40 x 59.8 x 59.8 cm)	11020
	Stainless steel cover plate (HxWxD 1.8 x 59.8 x 59.8 cm)	65310
Water treatment	MELAdem 53	01038
	MELAdem 53 C	01036
For documentation:	MELAflash CF card	01043
	MELAflash card reader	01048
	MELAprint 44 log printer	01144
	Ethernet adapter for MELAprint 42/44	40295
Process agents	5 l container for process agents	64010
	1 l container for rinse aid	60910
Others	Pre-filter	68130
	HEPA filter	51240
	Feed funnel	68200
	Container tap for a 5 l and 10 l container	70100

15 Documentation and approval

			Program/ load	Process successful?	Process approval?	Approval instruments?	Remarks	Signature
				<input type="checkbox"/> yes <input type="checkbox"/> no	<input type="checkbox"/> yes <input type="checkbox"/> no	<input type="checkbox"/> yes <input type="checkbox"/> no <input type="checkbox"/> Partially		
				<input type="checkbox"/> yes <input type="checkbox"/> no	<input type="checkbox"/> yes <input type="checkbox"/> no	<input type="checkbox"/> yes <input type="checkbox"/> no <input type="checkbox"/> Partially		
				<input type="checkbox"/> yes <input type="checkbox"/> no	<input type="checkbox"/> yes <input type="checkbox"/> no	<input type="checkbox"/> yes <input type="checkbox"/> no <input type="checkbox"/> Partially		
				<input type="checkbox"/> yes <input type="checkbox"/> no	<input type="checkbox"/> yes <input type="checkbox"/> no	<input type="checkbox"/> yes <input type="checkbox"/> no <input type="checkbox"/> Partially		
				<input type="checkbox"/> yes <input type="checkbox"/> no	<input type="checkbox"/> yes <input type="checkbox"/> no	<input type="checkbox"/> yes <input type="checkbox"/> no <input type="checkbox"/> Partially		
				<input type="checkbox"/> yes <input type="checkbox"/> no	<input type="checkbox"/> yes <input type="checkbox"/> no	<input type="checkbox"/> yes <input type="checkbox"/> no <input type="checkbox"/> Partially		
				<input type="checkbox"/> yes <input type="checkbox"/> no	<input type="checkbox"/> yes <input type="checkbox"/> no	<input type="checkbox"/> yes <input type="checkbox"/> no <input type="checkbox"/> Partially		
				<input type="checkbox"/> yes <input type="checkbox"/> no	<input type="checkbox"/> yes <input type="checkbox"/> no	<input type="checkbox"/> yes <input type="checkbox"/> no <input type="checkbox"/> Partially		
				<input type="checkbox"/> yes <input type="checkbox"/> no	<input type="checkbox"/> yes <input type="checkbox"/> no	<input type="checkbox"/> yes <input type="checkbox"/> no <input type="checkbox"/> Partially		

Glossary

A0-value

The A0 value represents a standard for the elimination of microorganisms and the deactivation of viruses in the disinfection procedure with damp heat. The A0 value depends on temperature and time.

AKI

Abbr.: working group instrument preparation ("Arbeitskreis Instrumentenaufbereitung")

Authorized personnel

An authorized personnel can be medical stockists, depot technicians or MELAG-authorized customer services trained by MELAG.

Batch

Collection of sterilization material which has been processed together in the same sterilization program.

Cleaner

The cleaning agent serves the removal of organic material (e.g. blood). The maximum pH value of the working solution of a mildly-alkaline cleaning agent is 11.

DI water

Demineralized water (also known as deionized water) is water (H₂O) which does not contain any of the minerals (salts and ions) present in normal spring or tap water.

Empty batch

Program run without a load or accessories (only with a basis basket).

Neutralizer

In addition to neutralizing the cleaning agent, the neutralizer serves to protect the instruments and brightens their stainless steel surface.

pH Value

The pH value is a measure of the strength of the acid or alkali effect of a watery solution.

Rinse aid

The rinse aid serves the subsequent rinsing of the instruments before drying. The rinse aid dries the load faster and without stain accretion.

VDE

Abbr. (German): "Verband der Elektrotechnik, Elektronik und Informationstechnik e.V." (Alliance of the Electronics, Electrotechnical and IT Industry).

MELAG Medizintechnik oHG

Geneststraße 6-10
10829 Berlin
Germany

Email: info@melag.com
Web: www.melag.com

Responsible for content: MELAG Medizintechnik oHG
We reserve the right to technical alterations

Your stockist

