

# Vcheck cProgesterone

BIONOTE Marketing team

Feb. 2020




# 01 Progesterone

- What is Progesterone?
- Application
- Guidelines

# 01 Progesterone

## What is Progesterone?

- **Progesterone** is a steroid hormone and part of the group of sex hormones. **The reproduction cycle** is the results of a rise and drop of different sex hormones.
- **Anestrus** : Resting ovary
- **Proestrus** : Ovarian follicles growing, Estrogen ↑  
Edema of the vulva, sanguineous discharge
- **Estrus** : Preovulatory luteinization of follicles. **Progesterone** ↑  
LH peak →(2 days) → **Ovulation**  
Clearer vaginal discharge, less oedematous vulva  
Allow mating of the male dog.
- **Diestrus** : Corpus luteum, Continuous rise in PRG



# 01 Progesterone

## Application

Species	Purpose
Canine	<ul style="list-style-type: none"><li>▪ <b>To determine optimal breeding dates. (Primary purpose)</b></li><li>▪ To predict parturition dates or time a Cesarean section.</li><li>▪ To detect reproductive disorders such as split heats, delayed puberty, silent estrus or hypoluteidism.</li></ul>

# 01 Progesterone

## Application

- Start **serial progesterone measurements** 3 to 4 days after a bloody vulvar discharge is first noted.
- Continue progesterone testing every other day until determination of ovulation and onset of the fertile period.


< 1.0 ng/mL (< 3.18 nmol/L)	1.0 – 1.99 ng/mL (3.18 – 6.33 nmol/L)	2.0 – 2.99 ng/mL (6.36 – 9.51 nmol/L)	3.0 – 4.99 ng/mL (9.54 – 15.87 nmol/L)	5.0 – 12.0 ng/mL (15.90 – 38.16 nmol/L)	> 12.0 ng/mL (> 38.16 nmol/L)
Anestrus or proestrus	Pre-LH surge	LH surge	Post-LH surge, Pre-ovulation	Ovulation (It may vary with breed and size.)	Post-ovulation

\* 1 ng/ml is equal to 3.18 nmol/L.

# 01 Progesterone

## Guidelines

Peak fertility typically occurs 4-7 days after the LH surge (or 2-5 days after ovulation).



# 01 Progesterone

## Guidelines

- **Natural breeding** – ideally breed every other day while the female is showing signs of standing heat. If only 2 breedings will be performed, attempt to breed 4 and 6 days after the progesterone predicted LH surge.
- **Fresh or chilled semen** – ideally inseminate 3 and 5 days after the progesterone predicted LH surge.
- It may vary with individual dogs.

# 02 Product Introduction

## Vcheck cProgesterone

- Specifications
- Key Features
- Test Procedure
- Reference Range
- Performance



# 02 Product Introduction

## Vcheck cProgesterone

- Specifications


- ✓ Species : Dog
- ✓ Sample : Serum 50 µl
- ✓ Testing Time : 15 minutes
- ✓ Measurement : Quantitative
- ✓ Measurement Range : 1.0 – 30.0 ng/mL
- ✓ Storage Condition : 2- 8 °C



# 01 Progesterone

## Performance

- Vcheck cProgesterone demonstrates excellent correlation with IMMULITE that used in reference laboratories.
- Vcheck cProgesterone can provide an accurate and reliable results in-house.



# 01 Progesterone

## Precaution

- Do not use SST (serum-separating tube)
- Hemolyzed samples or serum with precipitate such as fibrin may give erroneous results
- Sample storage Refrigerated (2 - 8°C): 7 days  
Frozen(-20°C): 6 months
- 1 ng/ml is equal to 3.18 nmol/L

# 01 Progesterone






## Precaution

- Cross reactivity with other steroid hormone

Cross-Reactivity substance		Cross(%)
<b>Estrone</b>	10ng/ml	<1%
<b>Estriol</b>	100ng/ml	<1%
<b>b-Estradiol</b>	20ng/ml	<1%
<b>17a-hydroxyprogesterone</b>	100ng/ml	5.64%
<b>Cortisol</b>	1000ng/ml	<1%
<b>Cortisone</b>	500ng/ml	<1%
<b>Corticosterone</b>	100ng/ml	<1%
<b>21-Deoxycortisol</b>	1000ng/ml	<1%
<b>Danazol</b>	1000ng/ml	<1%
<b>5-Pregnen-3<math>\beta</math>-ol-20-one</b>	200ng/ml	3.94%
<b>5<math>\alpha</math>-Pregnane-3,20-dione</b>	200ng/ml	5.14%
<b>17<math>\alpha</math>,20<math>\beta</math>-Dihydroxy- 4-pregnen-3-one</b>	200ng/ml	<1%
<b>Pregnenolone</b>	200ng/ml	6.33%

# 01 Progesterone

## Competitor


	BioNote	Eurolyser	Fuji	Boditech	Healvet
					
<b>Analyzer</b>	Vcheck V200	Solo/Cube-vet	Dri-chem immuno au10v	Vet chroma	HV-FIA3000
<b>Species</b>	Canine	Canine	Canine	Canine	Canine
<b>Specimen</b>	Serum 50 ul (Plasma will be add)	Serum/plasma 40/5 ul	Serum/plasma 100 ul	Serum/plasma 100 ul	Serum/plasma 75 ul
<b>Measurement time</b>	15 min	8 min	10 min	15 min	15 min
<b>Dynamic range</b>	1 – 30.0 ng/ml	1.2 – 8.0 ng/ml 6.0 – 20.0 ng/ml	0.2 – 40.0 ng/ml	1 – 25 ng/ml	1 – 50 ng/ml
<b>Storage</b>	2-8 °C	2-8 °C	2-8 °C	2-8 °C	4-30 °C
<b>Size</b>	10T	6T	10T	6T	10T

# 01 Progesterone

## Competitor – Eurolyser

- Eurolyser canine progesterone demonstrates **low accuracy** in comparison with reference method (Immulite progesterone) in internal study.

Correlation between Eurolyser progesterone and Immulite



$$y = \text{slope } 0.6301x + 0.2084$$
$$R^2 = 0.9278$$

Slope has to be within 0.8 - 1.2 according to Eurolyser's evaluation report.

# Thank you

BIONOTE Marketing team

Jan. 2020

