

Vcheck Canine NT-proBNP

BIONOTE Marketing team

Aug. 2020



Vcheck Canine NT-proBNP

01 NT-proBNP

02 Heart diseases in dogs


02 Product Introduction

- Vcheck Canine NT-proBNP

01 NT-proBNP

What is NT-proBNP?

- **B-type natriuretic peptide (proBNP)**
 - Produced in the muscle cells of the heart
 - Increases with excessive stretching of the cells
- proBNP is cleaved into BNP and NT-proBNP*
(*NT-proBNP: N-terminal pro-B type natriuretic peptide)



NT-proBNP is stable and has a long half-life, making it a more desirable biomarker.

⇒ Thus, it is used to assess the magnitude of cardiac muscle stretching

01 NT-proBNP

NT-proBNP in Dogs

- What NT-proBNP levels tell us -


- **Screens for asymptomatic DCM*** in large breeds
*DCM: Dilated Cardiomyopathy
 - ✓ Highly sensitive and specific for detecting DCM
- **Differentiates severe MMVD** patients from low-stage patients**
**MMVD: Myxomatous Mitral Valve Degeneration
 - ✓ Differentiates dogs with and without heart failure among the MMVD patients
- **Distinguishes cardiac from respiratory disease**
 - ✓ Useful in dogs with dyspnea requiring an emergency care



▲ A dog presenting difficulty breathing (dyspnea)

DCM (Dilated Cardiomyopathy)

- Dilated Cardiomyopathy (DCM)



The most common cause of heart disease in large breeds

(10% of all heart disease)



- Breed

- **Large breeds** (Boxer Dogs, Doberman Pinschers, Great Danes, Irish Wolfhounds, Saint Bernards);
- Rarely develops in small breeds

- Doberman Pinscher


- Very high risk of developing cardiomyopathy
- **Recommendation:** Echocardiogram + ECG (Holter) (Annually from 3-4 years of age)

- Clinical signs

- Lethargy, weakness, weight loss, collapse
- [Heart failure] dyspnea, increased respiratory rate,

DCM (Dilated Cardiomyopathy)

- Dilated Cardiomyopathy (DCM)



* Myocardial failure of the left side is most common.

Ventricular wall thinning; Heart muscle loses its contraction ability



The pressure of the blood inside the heart causes the wall to stretch

👉 Increasing NT-proBNP levels



Results in an enlarged heart



Can cause heart failure (pulmonary edema)

DCM (Dilated Cardiomyopathy)

✓ Diagnosis

Clinical examination

- **Symptoms:** Exercise intolerance, lethargy, weakness...
- **Auscultation:** Systolic murmur(-/+), gallop rhythm (-/+)

Electrocardiography (ECG)

- The best screening test (**VPC***: a common finding in early DCM)
- **Holter monitor:** 24 hour electrocardiogram (> 50 (100) VPCs/24 hours)
- **5-minute ECG:** > 1 VPC/5 min.

*VPC: Ventricular Premature Constriction



Radiograph


- Not a good screening test in the early stages
- **Severe disease:** Enlarged left atrium, pulmonary edema...



▲ DCM patient with heart failure

Echocardiogram

- Systolic dysfunction (FS < 20%)
- Enlarged left atrium / Increased LVIDs, LVIDd



▲ DCM patient (**Decreased FS**)

DCM (Dilated Cardiomyopathy)

✓ Diagnosis

Cardiac Biomarker: NT-proBNP

- A useful additional diagnostic test
- **Clinical utilities**
 - ① Screens for DCM in Even Asymptomatic Dogs
 - ② Evaluates the Survival Time
 - ③ Increases the Diagnostic Accuracy When Used With a Holter Monitor



DCM (Dilated Cardiomyopathy)


✓ Diagnosis

Cardiac Biomarker: **NT-proBNP**

- A useful additional diagnostic test
- **Clinical utilities**
 - ① [Screens for DCM in Even Asymptomatic Dogs](#)
 - ② Evaluates the Survival Time
 - ③ Increases the Diagnostic Accuracy When Used With a Holter Monitor

Evaluation of N-terminal pro-B-type natriuretic peptide as a diagnostic marker of various stages of cardiomyopathy in Doberman Pinschers

Gerhard Wess, DVM, Dr med vet, Dr habil; Verena Butz, DVM, Dr med vet; Monia Mahling, Dip Stat; Katrin Hartmann, DVM, Dr med vet, Dr habil



NT-proBNP Conc. for detecting all stages of DCM

Cut-off (400 pmol/L) Sensitivity 81%, Specificity 75%

Annual screening tests are recommended to determine if DCM is present even if your dog appears to be healthy.

DCM (Dilated Cardiomyopathy)

✓ Diagnosis


Cardiac Biomarker: **NT-proBNP**

- A useful additional diagnostic test
- **Clinical utilities**
 - ① Screens for DCM in Even Asymptomatic Dogs
 - ② Evaluates the Survival Time
 - ③ Increases the Diagnostic Accuracy When Used With a Holter Monitor

J Vet Intern Med 2012;26:1330–1336

Prospective Evaluation of NT-proBNP Assay to Detect Occult Dilated Cardiomyopathy and Predict Survival in Doberman Pinschers

G.E. Singletary, N.A. Morris, M. Lynne O'Sullivan, S.G. Gordon, and M.A. Oyama



In asymptomatic Dobermans with DCM

DCM (Dilated Cardiomyopathy)

✓ Diagnosis

Cardiac Biomarker: **NT-proBNP**

- A useful additional diagnostic test
- **Clinical utilities**
 - ① Screens for DCM in Even Asymptomatic Dogs
 - ② Evaluates the Survival Time
 - ③ Increases the Diagnostic Accuracy When Used With a Holter Monitor

J Vet Intern Med 2012;26:1330–1336

Prospective Evaluation of NT-proBNP Assay to Detect Occult Dilated Cardiomyopathy and Predict Survival in Doberman Pinschers

G.E. Singletary, N.A. Morris, M. Lynne O'Sullivan, S.G. Gordon, and M.A. Oyama


Criteria	Sens.	Spec.	Accuracy
NT-proBNP > 457	70 %	81 %	75.5 %
Holter monitoring	76.7 %	-	-
Holter or NT-proBNP > 457	94.5 %	87.8 %	91.0 %

▲ Diagnostic utility of NT-proBNP assay and Holter for the detection of occult DCM (In 155 asymptomatic Dobermans)

The addition of **NT-proBNP assay** to Holter recording as a diagnostic test improved sensitivity of detection!

DCM (Dilated Cardiomyopathy)

✓ Diagnosis (Algorithm)



MMVD (Myxomatous Mitral Valve Degeneration)

© 2016 Boehringer Ingelheim Vetmedica, Inc.

- **Myxomatous Mitral Valve Degeneration (MMVD)**



- ✓ **Most common heart disease**

(> 70% of all cardiovascular disease in dogs)

- **Age:** Older dogs
- **Breeds at elevated risk:** Small or toy breeds (< 9 kg)
 - Cavalier King Charles (Very high incidence)


- ✓ **Pathophysiology** (left heart)

- ① As dogs age, the valve becomes thickened
- ② It starts to leak backwards → **MMVD**
- ③ Abnormally high pressure
- ④ Could lead to pulmonary edema → **Heart failure**



MMVD (Myxomatous Mitral Valve Degeneration)

- ACVIM Consensus Classification



MMVD (Myxomatous Mitral Valve Degeneration)

• Diagnosis


Clinical examination

- **Symptoms:** Exercise tolerance, dyspnea, continuous panting, cyanosis ...
- **Auscultation:** Murmur, crackle sound


Radiograph

- **Indicator of heart size**
 - VHS (Vertebral Heart Size) > 10.5 vertebrae
 - VLAS (Vertebral Left Atrial Size) > 3 vertebrae
- **Pulmonary patterns**
 - Assessment of pulmonary edema


1. Stacey Fox-Alvarez, DVM, MPH, University of Florida
2. Malcolm, E. L., et al. 2018
3. Borgeat K. *In Practice* 2020



▲ Way to measure VHS in dogs



▲ Way to measure VLAS in dogs




▲ Heart failure signs caused by mitral valve disease

MMVD (Myxomatous Mitral Valve Degeneration)



• Diagnosis

Echocardiogram

- The best way to diagnose heart diseases
- **2D mode**
 - Regurgitation to LA (at systole)
 - LA:Ao ≥ 1.6
 - Epeak ≥ 150 mm/s
- **M mode**
 - LVIDd ≥ 25 mm (3kg)



▲ 2D mode



▲ 2D mode

MMVD (Myxomatous Mitral Valve Degeneration)

- **Diagnosis**

NT-proBNP (Biomarker)

- A useful clinical test for assessing the severity of heart disease
- **Clinical utilities**
 - ① Helps Discriminate Cardiac from Noncardiac Causes of Dyspnea
 - ② Discriminates Between Asymptomatic Dogs and Dogs With Heart Failure
 - ③ Monitors the Effectiveness of Medical Treatment for Heart Failure
 - ④ Predict Survival in Dogs With Heart Failure




MMVD (Myxomatous Mitral Valve Degeneration)

- **Diagnosis**

NT-proBNP (Biomarker)


- A useful clinical test for assessing the severity of heart disease
- **Clinical utilities**
 - ① [Helps Discriminate Cardiac from Noncardiac Causes of Dyspnea](#)
 - ② Discriminates Between Asymptomatic Dogs and Dogs With Heart Failure
 - ③ Monitors the Effectiveness of Medical Treatment for Heart Failure
 - ④ Predict Survival in Dogs With Heart Failure

Journal of Veterinary Internal Medicine  **Open Access**

J Vet Intern Med 2015;29:171–179

Relationship of Plasma N-terminal Pro-brain Natriuretic Peptide Concentrations to Heart Failure Classification and Cause of Respiratory Distress in Dogs Using a 2nd Generation ELISA Assay

P.R. Fox, M.A. Oyama, M.J. Hezzell, J.E. Rush, T.P. Nguyenba, T.C. DeFrancesco, L.B. Lehmkuhl, H.B. Kellihan, B. Bulmer, S.G. Gordon, S.M. Cunningham, J. MacGregor, R.L. Stepien, B. Lefbom, D. Adin, and K. Lamb



MMVD (Myxomatous Mitral Valve Degeneration)

- **Diagnosis**

NT-proBNP (Biomarker)

- A useful clinical test for assessing the severity of heart disease
- **Clinical utilities**
 - ① Helps Discriminate Cardiac from Noncardiac Causes of Dyspnea
 - ② [Discriminates Between Asymptomatic Dogs and Dogs With Heart Failure](#)
 - ③ Monitors the Effectiveness of Medical Treatment for Heart Failure
 - ④ Predict Survival in Dogs With Heart Failure


Veterinary Clinical Pathology ISSN 0275-6382

ORIGINAL RESEARCH

The diagnostic relevance of NT-proBNP and proANP 31–67 measurements in staging of myxomatous mitral valve disease in dogs

Johanna Wolf¹, Nicola Gerlach¹, Karin Weber¹, André Klima², Gerhard Wess¹

¹Clinic of Small Animal Medicine and ²Statistical Consulting Unit, LMU University, Munich, Germany



MMVD (Myxomatous Mitral Valve Degeneration)

- **Diagnosis**

NT-proBNP (Biomarker)

- A useful clinical test for assessing the severity of heart disease
- **Clinical utilities**
 - ① Helps Discriminate Cardiac from Noncardiac Causes of Dyspnea
 - ② Discriminates Between Asymptomatic Dogs and Dogs With Heart Failure
 - ③ Monitors the Effectiveness of Medical Treatment for Heart Failure
 - ④ Predict Survival in Dogs With Heart Failure


Received: 16 August 2017 | Revised: 17 April 2018 | Accepted: 3 May 2018
DOI: 10.1111/jvim.15228

Journal of Veterinary Internal Medicine **ACVIM**
American College of Veterinary Internal Medicine
Open Access

STANDARD ARTICLE

Effect of prespecified therapy escalation on plasma NT-proBNP concentrations in dogs with stable congestive heart failure due to myxomatous mitral valve disease

Melanie J. Hezzell¹ | Chloë L. Block¹ | Danielle S. Laughlin¹ | Mark A. Oyama¹




MMVD (Myxomatous Mitral Valve Degeneration)

- **Diagnosis**


NT-proBNP (Biomarker)

- A useful clinical test for assessing the severity of heart disease
- **Clinical utilities**
 - ① Helps Discriminate Cardiac from Noncardiac Causes of Dyspnea
 - ② Discriminates Between Asymptomatic Dogs and Dogs With Heart Failure
 - ③ Monitors the Effectiveness of Medical Treatment for Heart Failure
 - ④ Predict Survival in Dogs With Heart Failure

Journal of Veterinary Cardiology (2012) 14, 399–408



ELSEVIER




Journal of
Veterinary
Cardiology
ESVC

www.elsevier.com/locate/jvc

Lowered N-terminal pro-B-type natriuretic peptide levels in response to treatment predict survival in dogs with symptomatic mitral valve disease

Johanna Wolf, MRCVS ^a, Nicola Gerlach, DVM ^a, Karin Weber, DVM ^a, André Klima, Dip Stat ^b, Gerhard Wess, DVM, Dr. habil ^{a,*}



NT-proBNP 7–30 days after initial diagnosis of CHF

- **< 965 pmol/l**: had significantly longer cardiac survival times

02 Product Introduction

Vcheck Canine NT-proBNP

- Specifications
- Key Features
- Test Procedure & Reference Range

Vcheck Canine NT-proBNP

Specifications



- ✓ **Species** : Dog
- ✓ **Sample** : Serum 100 µl
- ✓ **Testing Time** : 15 minutes
- ✓ **Measurement** : Quantitative
- ✓ **Measurement Range** : 500 – 10,000 pmol/L
- ✓ **Storage Condition** : 2 - 8 °C

Key Features

- **Quantitative Analysis**
Numerical results, not qualitative
- **Proven Accuracy and Reproducibility**
Correlated against a laboratory ELISA method
- **Rapid and accurate results**
Simple procedure and quick results within 15 min.


Product No.	Product Name	Product Type	Packing Unit
VCF132DC	Vcheck Canine NT-proBNP	Device	5 Tests/Kit

Vcheck Canine NT-proBNP

Test Procedure & Reference range

The samples should be tested immediately after collection.


If serum samples are not tested immediately, freeze the serum for storage (-20 °C or colder).



Add 100 µl of the sample to the assay diluent tube



Mix well 5-6 times



Add the mixed sample (100 µl)

* Stored samples should be placed at room temperature 30 min. before use.

< 900 pmol/L	900 – 1,800 pmol/L	> 1,800 pmol/L
Normal	Suspected* Additional diagnostics are recommended	Abnormal* Additional diagnostics are recommended

* 'Abnormal' or 'Suspected' NT-proBNP test results should always be interpreted in combination and other diagnostic findings, such as an echocardiogram.

** Concentration over 735 pmol/L in Doberman Pinschers indicates an increased risk for occult DCM.

Thank you
Any Questions?

BIONOTE Marketing team

Aug. 2020

