

# Vcheck D-dimer

## Product Information

SEP 2019

BIONOTE MARKETING DIVISION






# Overview

## D-dimer Diagnosis

# D-dimer Diagnosis

## What is D-dimer?



**D-dimer** is a specific degradation fragment of cross-linked fibrin.

→ **High D-dimer**



- ① intravascular fibrin formation
- ② plasmin-mediated fibrinolysis

**Diagnosis** DIC (Disseminated intravascular coagulation)

Thromboembolism disease



# D-dimer Diagnosis

## When Is It Done?



# D-dimer Diagnosis

## DIC in Dogs



◀ petechiae

### Clinical signs

- Bleeding from the nose
- Blood in the eye
- Small pin-point bruising of the skin (called petechiae)
- Larger bruises (called ecchymosis)
- Increased heart rate
- Difficulty breathing
- Increased respiratory rate

# D-dimer Diagnosis

## DIC in Dogs

### Box 1. Clinical Conditions Potentially Associated with DIC

#### Systemic bacterial infections/sepsis

- Gram-negative bacteria (endotoxin)
- Gram-positive bacteria (bacterial coat mucopolysaccharide, enzymes)

#### Viral diseases

- Canine parvovirus
- Infectious canine hepatitis
- Feline panleukopenia
- FIP

#### Canine parasitic and rickettsial infections

- Babesiosis (*Babesia canis rossi*)
- Monocytic ehrlichiosis (*Ehrlichia canis*) and Rocky Mountain spotted fever (*Rickettsia rickettsii*)
- Leishmaniasis (*Leishmania infantum*)
- Spirocerosis (*Spirocerca lupi*)
- Heartworm disease (*Dirofilaria immitis*)
- Caval syndrome (*Angiostrongylus vasorum* and *Dirofilaria immitis*)

#### Neoplasia

- Solid tumors (e.g., mammary tumors)
- Lymphoid leukemia
- Myeloproliferative disorders
- Lymphoma
- Hemangiosarcoma (canine)
- Pulmonary adenocarcinoma

#### Immunologic disorders

- Immune-mediated hemolytic anemia
- Hemolytic transfusion reaction
- Transplant rejection

#### Vascular disorders

- Aortic aneurysm
- Vasculitis
- Hemangioma

#### Massive tissue injury

- Heatstroke and hyperthermia
- Gastric dilatation–volvulus (canine)
- Burns
- Head trauma
- Fat embolism
- Surgery (especially extensive and orthopedic)
- Severe mechanical trauma

#### Reaction to toxins

- Snakebite

#### Miscellaneous

- Pancreatitis
- Polycythemia
- Hepatic failure

# D-dimer Diagnosis

## DIC in Dogs


### At least 3 of the following criteria

- ❖ Abnormal aPTT, PT, or TCT value
- ❖ Low plasma fibrinogen concentration
- ❖ Low plasma AT III activity
- ❖ High serum FRA concentration
- ❖ Low platelet count



# D-dimer Diagnosis

## Thromboembolism in Dogs



Abnormalities in blood flow



Stasis of blood



Increased contact between platelets and coagulation factors with the endothelium



**Coagulation**


### KNOWN RISK FACTORS (Thromboembolism)

- ✓ Protein-losing disease (nephropathy or gastrointestinal)
- ✓ Cancer
- ✓ Sepsis
- ✓ Pancreatitis
- ✓ Congestive heart failure
- ✓ Immune-mediated disease
- ✓ Endogenous, exogenous corticosteroids




# D-dimer Diagnosis

## Thromboembolism in Dogs



### Venous Thromboembolism

- **Pulmonary thromboembolism** (cat, dog)




## Pulmonary Thromboembolism

- ✓ **Acute onset** of dyspnea, hypoxemia and hypocapnia & Tachycardia
- ✓ **Commonly have:** underlying cardiac disease, neoplasia, corticosteroid administration, DIC
- ✓ **Clinical signs:** nonspecific
  - lethargy
  - anorexia
  - weight loss

# D-dimer Diagnosis

## Algorithm for Pulmonary Thromboembolism (PTE) in dogs



*J Am Anim Hosp Assoc 2005;41:145-149.*

# D-dimer Diagnosis

**Table 1. Range and median D-dimer concentrations for all groups**

Group (n)	D-dimer (mg/l)	
	Range	Median
Group 1 (26)	0.1-0.5	0.2
Group 2 (9)	0.5-8.0	2.0
Group 3 (58)	0.1-6.9	1.5

Group 1: clinically healthy adult dogs

**Group 2: clinically ill dogs with TE/DIC**

Group 3: clinically ill dogs without evidence of TE/DIC

**Table 3. Distribution of D-dimer concentrations among groups investigated**

Group	D-dimer concentration (mg/l)					Median (mg/l)
	0.1 to 0.5	>0.5 to <1	1 to <2	2 to <3	>3	
Healthy	26					0.2
Neoplastic	3	3	6	4	3	1.7
Immune mediated	1	2	1	1		0.6
Inflammation	3	3	6	1	1	1.0
Postoperative			3		1	1.6
Haemorrhage	6		1	1	1	0.4
Miscellaneous	2	3		1	1	0.7
TE/DIC	1	1	2	2	3	2.0

DIC Disseminated intravascular coagulation, TE Thromboembolic disease

**Table 2. Plasma D-dimer concentrations in control, clinical illness, and TE dogs.**

D-dimer (ng/mL)					
	<250	250-500	500-1,000	1,000-2,000	>2,000
Control	30				
Neoplasia	9	1	4	1 <sup>a</sup>	1 <sup>a</sup>
Heart failure	8		1		
Liver disease	4	4	4	1 <sup>a</sup>	
Renal failure	5	2	1		
Postsurgery	11	3	6		
TE			4	9	7

TE, thromboembolic disease.

<sup>a</sup> Indicates 1 patient in category with hemoabdomen.

# D-dimer Diagnosis

**Table 1. Laboratory Screening Tests for DIC\***

Parameter/ Test	Early Hypercoagulable Phase	Clinical Manifestation Phase
Platelet count	=↓	↓
Schistocytosis	None	↑
PT	=↓	↑
aPTT	=↓	↑
Activated clotting time	=↓	↑
AT activity	=↓	↓
Fibrinogen	↓=↑ <sup>a</sup>	↓
FDP	=↓	↑
D-dimer	=↓	↑
Total protein C	=↓	↓
TAT	↑	↓=↑
PAP	↑	=↓


There is no single specific laboratory test for the diagnosis of **DIC**.

**D-dimer is still considered a good screening test for TE/DIC.**

# D-dimer Diagnosis

- A Good **Screening Test** For
  - **DIC (Disseminated intravascular coagulation)**
  - **Acute Thromboembolic Disease**
- **Assessment** of Pulmonary Thromboembolism
- **Monitoring** of **Antithrombotic therapy**
- Prediction of Survival Prognosis **after Surgery**

*AJVR, 64(12), 1562-1569, December 2003*



# D-dimer Diagnosis

## Factors which increase D-dimer

- ✓ **Neoplastic** (Ex. lymphoma)  
∴ a coagulopathy or hemorrhage from damaged vessels
- ✓ **Inflammatory**
- ✓ **Haemorrhagic**
- ✓ **Immune mediated**
- ✓ **Postoperative**
- ✓ **Miscellaneous** (idiopathic epilepsy and poisoning)

D-dimer concentration can be increased in a number of clinical conditions.

**Elevations in D-dimer concentration alone should not be used as the basis for the diagnosis of TE/DIC.**



# Product Overview

## Vcheck D-dimer

# Product Overview





# Product Overview :Vcheck D-dimer

## Specifications



- ✓ Species : Dog
- ✓ Sample : Plasma (only Sodium Citrate)
- ✓ Testing Time : 5 minutes
- ✓ Measurement Range : 0.1 – 10 µg/ml
- ✓ Storage Condition: 2 - 8 °C

# Product Overview :Vcheck D-dimer



## Test Procedure

Coding



Select [Standard Test] and insert a test device

Sample




Add sample 5µl into an assay diluent tube

Mix



Mix well by using a 100µl pipette

Apply



Load all of the mixed sample

# Product Overview :Vcheck D-dimer



## Reference Ranges


< 0.25 µg/ml	0.25 - 1.0 µg/ml	> 1.0 µg/ml
Normal	Elevated (Check other evidence of TE/DIC*)	Thromboembolic disease probable

\* TE: Thromboembolism, DIC: Disseminated intravascular coagulation

# Product Overview :Vcheck D-dimer



## Performance



## Good Clinical Utility

Researched by Haemaru Small Animal Clinical Research Institute & Referral Animal Hospital

- ✓ Stronger correlation with clinical sign
- ✓ High correlation with product 'N' ( $R^2=0.854$ )

*Thank you!*