

**SHELF LIFE:**

The **SureKlot™ Gauze** has a shelf life of 5 years.

**HOW SUPPLIED:**

7.6x50cm - SK750

7.6 x 100cm - SK7100



STERILE R



2460



EC REP

EMERGO EUROPE  
Prinsessegracht 20  
2514 AP The Hague  
The Netherlands

If you have any questions or comments about the product, please contact Axio Biosolutions Helpline at [care@axiobio.com](mailto:care@axiobio.com) or +91 8860786067

|           |                               |
|-----------|-------------------------------|
| STERILE R | Sterilized using irradiation  |
|           | Keep away from sunlight       |
|           | Keep dry                      |
|           | Do not reuse                  |
|           | Do not re-sterilize           |
|           | Do not use if pack is damaged |
|           | Date of manufacture           |
|           | Expiry date                   |
| LOT       | LOT number                    |

SKGZ-FUJ.03.22 Rev00

# SureKlot™

GAUZE

## INSTRUCTIONS FOR USE

AXIO®

**Manufactured by:**

Axio Biosolutions Private Limited  
Plot 18, Gujarat Pharma Techno Park  
Sari Matoda, Sanand Taluka  
Ahmedabad-382220, Gujarat, India.  
Mfg License no: MFG/MD/2020/000102  
Customer Care: +91 88607 86067

Email: [info@axiobio.com](mailto:info@axiobio.com) | Web: [www.axiobio.com](http://www.axiobio.com)

## **SUREKLOT™**

### **DEVICE DESCRIPTION:**

The **SureKlot™ Gauze** is a sterile, single-use, bioactive non-absorbable hemostatic dressing containing chitosan. It is intended to be used to control moderate to severe external bleeding in animals within minutes of application by providing an active mechanical barrier to the bleeding site. The **SureKlot™ Gauze** also provides a moist environment for wound management. It conforms to the wound and can serve as a primary or secondary dressing. It stops bleeding and protects wounds from infection and serves as an ideal product to be included in the first aid kits for animals or pets such as dogs, cats, and horses.

### **COMPOSITION:**

**SureKlot™ Gauze** is composed entirely of chitosan.

### **ACTION/BENEFITS:**

The **SureKlot™ Gauze** consists of positively charged chitosan molecules that enhance natural hemostasis by ionically attracting negatively charged blood particles such as platelets, red blood cells, and blood proteins and speed up the blood coagulation cascade. It also has a barrier to bacteria property and protects the wound from external microbial infection.

### **INTENDED USE:**

The **SureKlot™ Gauze** is intended to be used to control moderate to severe external bleeding and protects the wounds from external microbial infection in animals.

### **INDICATION FOR USE:**

The **SureKlot™ Gauze** is indicated for the management of moderate to severe external bleeding and protects the external veterinary wounds from infection. It can also be used to maintain a moist environment at the wound site.

- Lacerations
- Abrasions
- Post-operative and surgical stitches wounds
- Traumatic wounds
- Surgical debridement sites
- Maggot-infested wounds bleeding
- Necrotic wound and Cavity wounds
- Footrot
- Pus infection, paws infection, post skin disease management
- Bite wounds
- Burn wounds

### **DIRECTIONS FOR USE:**

1. Open the package and remove the **SureKlot™ Gauze** from it.
2. Cut the **SureKlot™ Gauze** to the desired size.
3. For acute wounds and bleeding control, wrap the wound using the **SureKlot™ Gauze** with uniform pressure until bleeding stops.
4. For chronic and infected wounds, debride when necessary and irrigate/clean the wound site by following standard wound management practices.
5. For removal, irrigate with water or saline and take it off
6. Can be combined with **SureKlot™ Spray** for debridement & decontamination, and providing moisture to dry and non-exuding wounds.

### **WARNINGS:**

- For external veterinary use only
- Non-absorbable dressing. Do not Implant
- Do not reuse
- Do not resterilize
- Do not use if the pack is open or damaged
- Contains chitosan from shellfish. Natural material tends to change color. Efficacy remains

### **STORAGE CONDITIONS:**

Store at ambient temperature. Avoid direct sunlight, heat, and moisture.