# beurer heart rate monitor herzfrequenz pulsuhr

germanengineering

## PM 18

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## TABLE OF CONTENTS

<b>1. Scope of delivery.....</b>	<b>2</b>
<b>2. Important Notes.....</b>	<b>2</b>
<b>3. Interesting Facts about this HR Monitor.....</b>	<b>4</b>
Measuring your pulse .....	4
Measuring the distance when running or walking .....	5
General information for training .....	6
<b>4. Functions of the HR monitor.....</b>	<b>7</b>
<b>5. General Operation of the HR Monitor.....</b>	<b>8</b>
General Operation of the HR Monitor.....	8
Buttons on the HR monitor .....	8
Display .....	9
Power saving mode .....	10
Main menus .....	10
Pulse measurement .....	11
Basic settings/Personal data .....	12
Overview.....	12
Setting the stride length for walking and running.....	13
Setting the sensitivity of the sensor for pace counting .....	13
Training .....	14
Set training zones.....	14
Display speed/distance.....	15
Display training data .....	16
Display calorie and fat consumption.....	16
Operate stop watch .....	17
Timer .....	17
Alarm clock.....	18
<b>6. Changing the Batteries .....</b>	<b>19</b>
<b>7. Troubleshooting .....</b>	<b>20</b>
<b>8. Technical Specifications and Batteries.....</b>	<b>20</b>
<b>9. Explanation of the Abbreviations on the Display.....</b>	<b>21</b>

## 1. SCOPE OF DELIVERY

Heart rate monitor with wrist band

In addition: These operating instructions and separate short instructions

## 2. IMPORTANT NOTES

Please read these operating instructions carefully, keep them for future reference and also make them available to others. Pay attention to all information in these operating instructions!

### Training

- Did you know that doing daily at least 10,000 paces or exercising for half an hour (less for older persons, more for children) is good for your health.
- This HR monitor is not a therapeutic device. It is a training instrument which was developed for the measurement and representation of the human heart rate.
- With high-risk sports, please note that the use of the HR monitor can represent an additional source of injury.
- If in doubt or in the case of illnesses, please consult your doctor about your ideal upper and lower training heart rate, as well as the recommended duration and frequency of the training. In such a way, you can achieve optimal results when training.
- **WARNING:** People suffering from cardiovascular diseases or those with pacemakers should only use this heart rate monitor after having consulted their doctor.

## Application

- The heart rate monitor is only intended for private use.
- The heart rate monitor may only be used in the manner outlined in the operating instructions. The manufacturer is not liable for damages that are caused by improper or incorrect use.
- The heart rate monitor supports you in many types of training (for example even when cycling), however all the functions can only be used during running training.
- This HR monitor is waterproof (50m) and thus suitable for use when swimming. However you must not press any buttons in the water.

## Cleaning and care

- From time to time, clean the HR monitor carefully with soapy water. In addition, once in a while the back of the monitor in particular should be wiped with a damp cloth. Then rinse the complete HR monitor with clean water. Carefully dry it with a soft towel.
- Do not use skin cream if you would like to use the HR monitor.
- Have the battery changed by a specialized dealer (CR2032 3V lithium battery).

## Repair, Accessories and Disposal

- Prior to use, you have to ensure that the HR monitor is free of damage. If in doubt, do not use the HR monitor and consult your sales representative or the customer service address provided.
- Clean your skin and fingers with soap to ensure a better transmission of the signal when measuring your pulse.
- Transmission problems can occur in the case of very dry or thick skin.

- If your pulse is not measured correctly or not at all: Wet the tips of your fingers or use a conductive paste.
- Avoid twisting and bending your wrist when measuring.
- Repairs may only be conducted by customer service or authorized dealers. Please do not in any case attempt to repair the HR monitor! In the case of complaints, please consult your dealer or the service address.
- If you have a complaint, please replace the battery at a specialized dealer beforehand.
- The warranty expires if the monitor is opened or due to improper use.
- Use the HR monitor only in combination with Beurer original accessories.
- Avoid contact with suntan lotions or the like, since these can cause damage to the imprint or plastic components.
- Please dispose of the HR monitor and used batteries only in accordance with the Waste Electric and Electronic Equipment Directive 2002/96/EC – WEEE. For inquiries, please contact the municipal authority responsible for disposal in your area.


### **3. INTERESTING FACTS ABOUT THIS HR MONITOR**

#### **Measuring your pulse**

With this HR monitor (without chest strap) you can monitor your heart rate and set a training zone corresponding to your training objectives with an upper and lower limit value. An exact ECG heartbeat is calculated and displayed without having to wear the annoying chest strap.

To do this, strap the monitor to your wrist and touch the metal application on the front of the monitor with your index finger or thumb. Your heart rate can be determined ECG-exact within a few seconds by measuring from the metal application via the arm, upper part of the body, other arm and on to the sensor plate on the back of the HR monitor. The monitor has to have direct contact with your skin so that the sensor plate on the back of your HR monitor can determine your pulse.

This HR monitor is thus suitable in particular for people who cannot or do not want to wear a chest strap and for sports such as hiking or walking. Together with the running speed, the pulse measurement can also calculate the calories consumed as well as the fat burned during training.

Tip: Measure your heart rate at regular intervals – in particular if your heart rate increases or decreases -, or at the end of each training section!

## Measuring the distance when running or walking

Sufficient movement is required for healthy living. At least 7,000 - 10,000 paces or 30 minutes exercise a day is optimal (Source: acc. to WHO, Dr. Catrine Tudor-Locke, Dr. David R. Bassett jr.).

With this HR monitor you have the possibility to check your daily amount of exercise. An integrated activity sensor provides an overview of your daily movements and calculates the number of paces, amount of time, the distance covered and the speed. On existing training monitors a special sensor had to be attached to your foot to calculate the distance. Now no external sensor is required because of the integrated activity sensor.

You can specify the typical stride length in the basic settings of the monitor so that the HR monitor can calculate the distance trained and

your running speed from the stride length and number of paces. You can also adjust the sensitivity of the integrated sensor so that the result does not depend on your personal movement activity.

## General information for training

This heart rate monitor serves primarily to measure the human pulse. By means of various settings, you can support your individual training program and monitor your pulse. You should monitor your pulse during exercise because you should prevent your heart from excessive strain on the one hand and on the other in order to achieve the optimal training results. The table provides you with some information on the selection of your ideal training zone. The maximum heart rate reduces with advancing age. Information on the training heart rate must therefore always be in relation to the maximum heart rate. The following rule of thumb helps to identify the maximum heart rate:

### **220 – age = maximum heart rate**

The following example applies to a 40-year-old person:  $220 - 40 = 180$

	Heart/health programme	Fat burning zone	Fitness programme	Endurance training programme	Anaerobic training
<b>Max. heart rate range</b>	50-65%	55-75%	65-80%	80-90%	90-100%
<b>Display on monitor</b>	HEALTH	FAT-B	AROBIC	–	–
<b>Effect</b>	Strengthens the heart and the circulatory system	The body burns the highest percentage of calories from fat. Strengthens the heart and the circulatory system, enhances fitness	Improves the respiratory and circulatory system. Ideal to promote the basic endurance	Improves speed maintenance and increases basic speed.	Overloads specific muscles. High injury risk for athletes at amateur level, risk to the heart in case of disease.
<b>Suitable for</b>	Ideal for beginners	Weight control and loss	Athletes at amateur level	Ambitious athletes at amateur level, professional athletes	Only professional athletes
<b>Training</b>	Regenerative training		Fitness training	Endurance training	Fitness training addressing specific zones

Three pulse ranges have already been preset and stored in the HR monitor (HEALTH/FAT-B/AROBIC). See chapter “Setting your training zone” to select and set your pulse limits.

## 4. FUNCTIONS OF THE HR MONITOR

### HR functions

- Exact ECG heart rate measurement
- The heart rate symbol flashes during measurement
- Individual training zone adjustable
- Acoustic and visual alarm for exit from the training zone
- Calorie consumption in Kcal (based on the walking or running speed)
- Fat burning in g (grams)
- Training zone proposal

### Activity functions

- Number of paces
- Total time of activity
- Daily deletion of training/activity values at 3:00 hours

### Speed functions

- Distance
- Current speed
- Average speed

### Time functions

- Time of day
- Calendar (up to 2099) and weekday
- Date
- Alarm clock
- Hour signal
- Pushbutton signal
- Stop watch
- Countdown timer

### Settings

- Gender, weight, height, year of birth
- Unit of measure selection
- Maximum heart rate
- Stride length for walking and running
- Sensitivity of sensor for pace counting

## 5. GENERAL OPERATION OF THE HR MONITOR


### General Operation of the HR Monitor

#### Buttons on the HR monitor

1	LIGHT	Pressing quickly illuminates the display for 3 seconds
2	MODE	Switches between the different menus and setting displays. If you keep the MODE button pressed for at least 5 seconds, you get into the setting mode for a menu. Exit the setting mode also by keeping the MODE button pressed for 5 seconds.
3	START/ STOP/ RESET	<ul style="list-style-type: none"><li>■ Start/stops functions, sets values. Keep pressed to accelerate input.</li><li>■ Counts values forwards.</li><li>■ Starts and stops the stop watch.</li><li>■ Resets the stop watch function and the results (Kcal etc.).</li></ul>
4	SELECT	<ul style="list-style-type: none"><li>■ Sets values.</li><li>■ Counts values backwards.</li><li>■ Display in training mode.</li></ul>
5	SENSOR	Metal application on the front of the HR monitor. Starts the pulse measurement automatically when touching the metal surface.

## Display

1 Top line

2 Bottom line

3 **The most important symbols:**

Hour signal is activated

Heart rate is being received

Training zone has been exceeded

Training zone has been underrun

Alarm clock is activated

Heart rate is within the pulse limits


To continue getting exact results, the battery has to be replaced as soon as possible after this symbol has been displayed.

An explanation of the display abbreviations can be found in the list of abbreviations at the end of the operating instructions.


## Power saving mode

Hold the **SELECT** button in real time mode for approx. 5 seconds to enter the power saving mode. The display is turned off in this mode, so that the battery can last longer. Please note that the functions (e.g. Time) are still working and the display is just turned off. Press any button to reactivate the display again.

## Main menus

The operation of the HR monitor is organized into three main menus. From each of the main menus, you can then switch into the relevant setting menus. Press the **MODE** button to switch between the main menus. You get into the setting menus by keeping the **MODE** button pressed for about 5 seconds. In the same way you can also exit the setting menus again and return to the relevant main menu.

## Pulse measurement

The key element of your HR monitor is to record the heart rate. The simple operation makes it easy for you during training to determine your current heart rate in between.

Put your index finger or thumb on the metal surface of the casing. Pay attention that the back of the monitor has direct contact with your skin.


Press your finger slightly on the metal surface of the casing for 3-7 seconds continuously. The heart symbol and the two horizontal lines in the middle flash during measurement.

As soon as your heart rate has been determined, it is displayed in beats per minute.

In addition, on the right you can see a percentage value which indicates where your pulse is in comparison to the recommended/set maximum pulse. An arrow appears next to the ZONE symbol, it indicates whether the heart rate exceeds or falls short of the set values.

As soon as you remove your finger or press a button, you return to the mode last used.


Please note: Do not measure your pulse under water!

## Basic settings/Personal data

### Overview

In this menu you can set the time, specify basic functions of the HR monitor and input personal settings.

From the time display, press the MODE button for approx. 5 seconds. You get into the setting menu with a number of settings; move on to the next setting by pressing the MODE button. The default settings flash if you can change them. Change these settings with the START/STOP/RESET and SELECT buttons.

Exit the setting menu by keeping the MODE button pressed for approx. 5 seconds.

**Tip:** If you keep the START/STOP/RESET and SELECT buttons pressed for more than 2 seconds, then the relevant setting counts forwards or backwards automatically until you let the button go.

## **Setting the stride length for walking and running**

An important presetting for the distance and running speed function is to enter your individual stride length as accurately as possible.

1. Run or walk a distance where you know the exact length (reference distance). For example the 100 meter track in the athletics stadium.
2. Count the paces you require for it.
3. Divide the reference distance by the number of paces, and you will get the stride length (e.g.  $100 \text{ m} / 125 \text{ paces} = 0.8 \text{ m} = 80 \text{ cm}$ ).
4. Record this value in the Stride length menu for walking or running, depending on the relevant type of pace.

If you change the type of pace during training, then this has an effect on the accuracy of the speed and distance recordings.

## **Setting the sensitivity of the sensor for pace counting**

To achieve an optimum result, the sensitivity of the sensor should be individually set for the user. Five sensitivity levels can be selected here (1=Low sensitivity, 5=High sensitivity). To get the optimal setting, we recommend test runs before the first training. The slower your running speed, the higher you should select the sensitivity level.

## **Using the monitor when jogging/walking**

Walk or run a distance in your training speed (at least 50 paces). Do not change your normal running arm movements. Compare the counted paces with the pace display of the HR monitor, see Chapter "Display training data". If the monitor does not display enough paces, then you have to increase the sensitivity of the sensor. If too many paces are displayed, then reduce the sensitivity. Once you have ad-

justed the sensitivity of the sensor to your pace, then you can achieve 95-100% accuracy.

Please note: A change in your running style or running speed can affect the accuracy.

### **Use as activity sensor**

To prevent the sensor from mistakenly interpreting your hand movements as paces, you should select a lower sensitivity level of the sensor. The sensitivity needs to be adjusted depending on the day's activities.

Please note: The HR monitor should be removed for activities, where only hand movements are carried out (e.g. cleaning windows) to avoid incorrect results.

## **Training**


First of all you should set your training zone. During training you can switch between various displays.

### **Set training zones**

You can set your individual upper and lower heart rate limit here by selecting your optimal training load. To do this, you can choose from three default values:

- optimal health [**HEALTH**] (50-65 % of the max. heart rate)
- fat burner [**FAT-B**] (55-75 % of the max. heart rate)
- fitness [**AROBIC**] (65-80 % of the max. heart rate)

Or you can specify your own limits for the heart rate. In addition, you can also set the acoustic alarm, if your heart rate is outside your pulse limits during training.

Press the MODE button in the main menu Training Mode for approx. 5 seconds. You get into the setting menu; move on to the next setting by pressing the MODE key. The settings flash if you can change them. Change the settings with the START/STOP/RESET and SELECT buttons. Confirm with the MODE button.

Then leave the setting menu again by keeping the MODE button pressed for approx. 5 seconds.

Tip: If you keep the START/STOP/RESET and SELECT buttons pressed for more than 2 seconds, the relevant setting automatically counts forwards or backwards.

## Display speed/distance

In the bottom line of the display, you can generally select between the two main displays (speed and distance) as well as automatically switch between these main displays. If you prefer to use one particular display, then you can keep this for your entire run.

If you press the SELECT button in training mode for 2 seconds, you can switch the display between speed (SP) and distance (DI). If you

press the SELECT button for 4 seconds, you enter the automatic switching mode between speed and distance (SCAN). In SCAN mode the display changes every 2 seconds. To stop the SCAN mode press the SELECT button for 2 seconds.

## Display training data

In the top line of the display, you can have different results displayed with the SELECT button:

- Training time
- Kilocalories (kC) in kilocalories [kcal]
- Fat burning (g) in grams [g]
- Distance (DI) in kilometres [km]
- Momentary speed (SP) in kilometres per hour [km/h] or miles per hour [mi/h]
- Average speed (SP with AVG symbol)
- Number of paces (ST=Steps)
- Activity time (M=Motion)
- Time of day

After pressing SELECT again, the display returns from the last screen to the first screen.

## Display calorie and fat consumption

In the training data display, you can also display the number of calories consumed and the fat burning during training when scrolling with the SELECT button

**Please note:** The calculation of the calorie and fat burning values is based on your running speed during training. The calorie and fat burning values cannot be measured exactly when cycling.

## Operate stop watch

In the Stop Watch mode, there is a zero “000000” in the top line.


- Start the stop watch with the START/STOP/RESET button.
- Stop the stop watch by pressing the START/STOP/RESET button again.
- To reset the stop watch and to delete all stored training and activity data, keep the START/STOP/RESET pressed for at least 2 seconds. **Reset all** and then **Reset done** is displayed.

## Timer


The timer has a stop watch that counts time backwards. The HR monitor emits a short acoustic alarm as soon as the counter has less than 10 seconds left to count. A continuous acoustic alarm sounds when it gets to zero, it can be switched off by pressing any button.

### Display the Timer mode

The current time and the remaining time of the timer are displayed in this mode.

### Settings of the Timer mode

Press the MODE button for approx. 5 seconds. The settings flash if you can change them. Change the settings with the START/STOP/RESET and SELECT buttons. Switch between the Hours, Minutes and Seconds settings with the MODE button, at the same time the new values are stored.

Leave the setting menu again by keeping the MODE button pressed for approx. 5 seconds.

Tip: If you keep the START/STOP/RESET and SELECT buttons pressed for more than 2 seconds, then the relevant setting automatically counts forwards or backwards.

You can start and stop the timer at any time with the START/STOP/RESET button. To reset the timer settings press the SELECT button for 2 seconds.

## Alarm clock

Your HR monitor has an alarm clock, which emits an acoustic signal when it reaches the wake-up time. Press any button (apart from LIGHT) to stop the alarm.

### Display the alarm clock


Press the START/STOP/RESET button when you are in the time display

In this mode, the set wake-up time is displayed (in the picture: 12 o'clock) and whether the alarm clock is activated (in the picture it has been switched off "OFF").

As soon as you release the START/STOP/RESET button, you get back into the time display.

## Setting the alarm clock

Press the START/STOP/RESET button in the time display for at least 2 seconds. You get into the setting menu; move on to the next setting by pressing the MODE key. The settings flash if you can change them. Change the settings with the START/STOP/RESET and SELECT buttons.


Then leave the setting menu again by keeping the MODE button pressed for approx. 2 seconds.

Tip: If you keep the START/STOP/RESET and SELECT buttons pressed for more than 2 seconds, then the relevant setting automatically counts forwards or backwards.

## 6. CHANGING THE BATTERIES

Have the battery changed by a specialized dealer (CR2032 3V lithium battery).

Used batteries do not belong in the household garbage. Dispose of these through your electronics dealer or your local collecting point for recoverables. You are legally required to do this. We will find this symbol on batteries containing toxic chemicals: Pb = battery contains lead, Cd = battery contains cadmium, Hg = battery contains mercury.

## 7. TROUBLESHOOTING

### Dry skin or very hairy arms

Use a conductive paste or wet your fingers and wrist.

### Your finger does not have proper contact with the sensor

Your finger (not only the finger tips!) has to lie flat and firmly on the metal casing and the monitor has to be strapped skin-tightly to your wrist.

### Irregular heartbeat

A correct heart rate measurement is very difficult in this case, irregular measurement periods can be expected. The result will not be reliable in all cases!

### Movement during the pulse measurement

Avoid moving intensively during the pulse measurement to get a reliable result.

## 8. TECHNICAL SPECIFICATIONS AND BATTERIES

HR monitor	<p>Waterproof to 50 m (suitable for swimming)</p> <p> The impermeability of the monitor to moisture cannot be guaranteed when jumping into water (high pressure load) and also when pressing the buttons in the rain.</p>
Battery	<p>3V Lithium battery, CR2032 type.</p> <p>Battery life approx. 12 months (depending on the number of pulse measurements). The battery already inserted upon purchase is for the test run and can be flat within a year.</p>

## 9. EXPLANATION OF THE ABBREVIATIONS ON THE DISPLAY

Abbreviations on the display in alphabetical order:

Abbreviation	Explanation
[cm]	Centimetre
[IN]	Inch
[km/h]	Kilometres per hour
[mi/h]	Miles per hour
AL	Alarm clock (OFF=switched off, ON=switched on)
AROBIC	Fitness: 65-80 % of the max. heart rate
BEEP	Pushbutton signal
BIRTH	Year of birth (4-digit, e.g. "1978")
CHIME	Hour signal
DAILY	Daily
DI	Distance in kilometres [km]
FAT-B	Fat Burner: 55-75 % of the max. heart rate
Fri	Friday
g	Fat burning in grams [g]
Gender	Gender (M=male, F=female)
HEALTH	Optimal health: 50-65 % of the max. heart rate
HGT	Height in inches [IN] or [CM] centimetres
IMPERI	Imperial units of measure
kC	Calorie consumption in kilocalories [kcal]
M	Activity time with the format [Hour:Minute:Second]
METRIC	Metric units of measure
Mon	Monday
RUN	Stride length for running in centimetres [CM] or inches [IN]
Sat	Saturday
SCAN	Automatic display change between speed and distance
SENSE	Sensitivity of the activity sensor
SP	Momentary speed in [km/h] or [m/h]
SP AVG	Average speed in [km/h] or [m/h]

Abbreviation	Explanation
ST	Number of paces
Sun	Sunday
Thu	Thursday
TR	Timer Mode
Tue	Tuesday
UNIT	Unit of measure in Imperial [IMPERI] or Metric [METRIC]
WALK	Stride length for walking in centimetres [CM] or inches [IN]
Wed	Wednesday
WGT	Weight in kilograms [kg] or pounds [LB]
ZONE	Training limits

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Beurer GmbH, Söflinger Str. 218, 89077 Ulm, Germany