


How to Wear?

- Please ensure you order the right Support as per the Size Chart below
- Open all the straps and lay the support in front of you with the Oppo logo facing upwards
- First slip on the sleeve onto the injured shoulder and settle it
- Adjust the hook and loop closure around the biceps
- Roll the two elastic straps around the front and back of the shoulder - seek assistance if you find this difficult due to the injury
- Fasten the two straps together under the armpit of the healthy shoulder to desired compression and restrained mobility
- Readjust the fastener on the injured shoulder once again



REMEMBER

- Please do not wear this support continuously except unless recommended by your Physician
- Don't bathe with the support - it will reduce the elasticity and compression quality of the support

- Be careful while fastening the straps, don't leave the hook section exposed to avoid any rashes over the day
- Never share your support - can lead to a transfer of infection or allergy

P.S.: Kindly consult your Physician or Orthopedic Specialist before using this product