


How to Wear?

- Make sure you order the right fit as per Size Chart below
- Insert foot into sleeve
- Pull up towards the thigh and stop when it covers the upper and middle thigh
- Check for compression and remove the support if it feels too sweaty to dry the skin and fabric
- To remove, pull the support from the bottom and gradually slide it
- **REMEMBER:** Do not roll or fold the support while wearing or removing it. This might damage the structural integrity of the product over time



Wear the Thigh Support as per above directions. To remove, simply reverse by pulling from the lower part of the support.

P.S.: Kindly consult your Physician or Specialist before using this product