

# QUICK START GUIDE

## THE DEVICE -

- Untie the velcro strap and press the power button 'ON' to light up the sensor
- Put your finger through the velcro strap close to the sensor receiving and transmitting tubes
- Adjust the velcro tape to the most suitable fitness and fasten the velcro tape
- Wait for a few moments till the measurement value is stable and can be steadily read on the screen
- The oximeter will automatically shut down after 10 seconds when no signal is generated



## THE APP-

### For Android OS:

- Search OxyCare app on Google Play
- Install and launch the app
- On the search field, provide the name of the Oxymeter you are using and tap on the search result displayed
- Once the device is connected via the app measurements will be displayed on the mobile screen
- Record the readings share it with your physician if needed

### For iOS:

- Search Oxycare app on Apple app store
- Install and launch the app and add a User name
- On the search field, provide the name of the Oxymeter you are using and select the search result displayed
- Once the device is connected via the app measurements will be displayed on the mobile screen
- Record the readings share it with your physician if needed

## P.S.

- For the best result make sure the hand is not trembling during measurement
- Index finger joining the first section is recommended as a suitable monitoring position

## REMEMBER

- Do not use oximeter in a magnetic resonance environment

- The device might not work well under low battery condition
- For more information please refer to the product package list