EN Compact Nebulizer

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DE Kompakt Vernebler

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EN English

Introduction

Thank you for purchasing a NB80 Compressor Nebulizer. With proper care and use, your nebulizer will provide you with many years of reliable treatments. This unit operates on standard AC power. Treatments are delivered quickly, safely and conveniently making this unit ideal for all ages. We encourage you to thoroughly read this guidebook to learn about the features of your nebulizer. Your compressor nebulizer should be used under the guidance of your physician and a respiratory therapist. Together with your physician and/or respiratory therapist, you can feel comfortable and confident knowing that you are obtaining the most effective inhalation treatments for your respiratory condition.

NOTE: Your nebulizer is intended for use in treatment of asthma, COPD and other respiratory ailments in which an aerosolized medication is required during therapy. Please consult with your physician and/or pharmacist to determine if your prescription medication is approved for use with this nebulizer. For type, dose and regime of medication follow the instructions of your doctor or respiratory therapist.

This device fulfills the provision of the EC directive 93/42/EEC (Medical Device Directive) and the European Standard EN 13544-1:2007+A1:2009 (Equipment - Part 1: Nebulizer systems and their components).

Please read this manual carefully before use and be sure to keep this manual.

Cautions

Please use general safety precautions when operating your nebulizer. This unit should be used only for its intended purpose as described in this guidebook and with medications only under the supervision and instruction of your physician. Do not use the device in anesthetic or ventilator breathing circuits.

Product cautions

RECOMMENDED USE BEFORE USING

- To avoid electrical shock, keep unit away from water.
- Do not immerse the unit in liquid.

Do not reach for a unit that has fallen into water immediately unplug the unit.

Do not use the unit if it has any damaged parts (including plug), if it has been submerged in water or dropped. Promptly send the unit for examination and repair.

The nebulizer should not be used where flammable gas, oxygen or aerosol spray products are being used.

Keep the air vents open.

Wipe the compressor daily using a soft cloth.

NOTE: Any other form of cleaning or cleaning agents may damage the finish of the unit.

Cleaning the compressor

Wipe the compressor daily using a soft cloth.

NOTE: Any other form of cleaning or cleaning agents may damage the finish of the unit.

Changing the air filter

It is important to change the air filter approximately when the air filter turns gray. It is recommended to change air filter every 2 months.

The unit should not be used until unfiltered.

Do not tilt or shake the unit when in operation.

Replace with a new, clean air filter.

Securely replace the air filter cover to the unit.

NOTE: Air filters cannot be cleaned or washed. Only NB80 air filters can be used. Do not substitute alternate material such as cotton. Do not operate without an air filter.

Troubleshooting

If any abnormality occurs during use, please check and correct the following:

1. Unit does not operate when power switch is pressed: Check the AC connection to the outlet.

2. No mist or low rate of misting:

• Check that there is medication in the nebulizer cup.

• Clean the main unit if there is any physical damage.

• Check the position of the nozzle inside the nebulizer.

• Make sure that air tube and all components are properly attached.

• Check the air filter and replace if necessary.

Protection against electric shock

• Do not store the unit in direct sunlight; high temperature or humidity.

• Always keep the unit away from small children.

• Cleaning caution

• Check air filter, nebulizer, mouthpiece and any other optional component before each use. Dirty or worn parts should be replaced.

• Do not immerse the unit in water. It may damage the unit.

• Disconnect the unit from the electrical outlet before cleaning.

• Clean all necessary parts after each use as instructed in this guidebook.

• Always dispose of any remaining medication in the medication cup after each use.

• Do not use the air tube/mouthpiece or medication remaining in the air tube. This could result in infection as a result of bacteria.

MEDICAL DISCLAIMER:

This manual and product are not meant to be a substitute for advice provided by your doctor or other medical professionals.

Don't use the information contained herein or this product for diagnosing or treating a health problem or prescribing any medication. If you have or suspect that you have a medical problem, promptly consult your doctor.

Durable periods are as follows, provided the product is used to nebulize 2ml of medication 2 times a day for 8 minutes each time at room temperature (25°C). Durable period may vary depending on usage environment.

• Main unit 5 years

• Nebulizer Kit 1 year

• Air Tube, Mouthpiece 1 year

• Air Filter 60 days

• Adult and child masks 1 year

Product specifications

AC Adapter DC12V, 1A

AC100~240V, 50~60Hz

Power Consumption ≤ 12 W

Sound Level ≤ 53 dBA (1 meter away from NB80)

Compressor Pressure Range ≥ 16 psi (110 kPa)

Operating Pressure Range ≥ 8 psi (55 kPa)

Operating Flow Range ≥ 3.0 lpm

Operating Temperature Range 50°F to 104°F (10°C to 40°C)

Operating Humidity Range ≤ 90% RH

Operating Atmospheric Pressure Range 700-1000 hPa

Storage Temperature Range -4°F to 140°F (-20°C to 60°C)

Storage Humidity Range ≤ 90 % RH

Dimension (L x W x H) 145mmx109mmx73mm (5.71x4.29x2.87 inches)

Weight 600g (without accessories)

Medication Capacity 5ml(cc)

Particle Size (MMAD) ≤ 2.5μm

Average Nebulization Rate Fully open Valve > 0.3ml/min (0.9% Saline Solution)

Closed Valve ≥ 0.1ml/min (0.9% Saline Solution)

Standard Accessories Nebulizer Kit, Air Tube, Mouthpiece, Filters (5pcs), Adult and child masks

*Subject to technical modification without prior notice.

*Performance may vary with drugs such as suspensions or high viscosity. See drug supplier's data sheet for further details.

A Product identification

1. Nebulizer Kit 2. Nozzle 3. Angled Mouthpiece

4. Air Filter 5. Child Mask 6. Adult Mask

7. AC Adaptor 8. Air Tube 9. Filter Slot

10. Air Output 11. Power Switch 12. AC Adaptor Jack

B Valve Adjustable Technology

The proprietary adjustable valve is able to deliver medications of different viscosity level according to every user's conditions and needs at ease. Our VA technology allows users to adjust different levels of nebulization rate ranging ≥ 0.1 ~ 0.3 ml/min at a consistent particle size.

The adjustable valve can be adjusted by the user in a very easy way without exchanging parts. Higher nebulization rate is open for higher viscosity medications and higher breathing capacity user while lower nebulization rate with closed valve will be more appropriate for kids / infants with lower breathing capacity.

C Assembling your nebulizer

Following the cleaning instructions in this guidebook under "Cleaning technique" prior to using your nebulizer for the first time or it has been stored for an extended period of time.

REMEMBER: Always unplug the compressor and make sure the power-switch is turned to the "OFF" position before cleaning, assembling and before or after each use.

1. Place the compressor on a flat, stable surface within reach.

2. Gently twist and pull straight up the lid of the nebulizer to separate into two parts (medication cup and cover).

3. Make sure the nozzle is properly installed on the upper cover. The stem inside the medication cup fits into the tube of the nozzle.

4. Add the prescribed amount of medication to the medication cup.

5. Reassemble the nebulizer by carefully twisting the medication cup and cover together. Make sure that the two parts fit securely.

WARNING: The symbol on this product means that it's an electronic product and following the European directive 2012/19/EU the electronic products have to be disposed on your local recycling centre for safe treatment.

D Operating your nebulizer

The nebulizer is operable up to a 45° angle. If the angle is greater than 45°, no aerosol will be generated.

1. Attach one end of the air tube connector to the air output.

2. Carefully attach the opposite end of the air tubing connector to the stem at the base of the nebulizer kit.

3. Attach the angled mouthpiece or mask to the top of the nebulizer.

4. The capacity of the medication cup is 2~5 ml.

NOTE: A 30-minute interval is recommended after each use.

Operating your nebulizer

After every use:

1. Unplug the unit from the power source.

2. Turn the unit to completely cool.

3. Carefully detach the air tubing from the nebulizer and pour out any remaining medication.

4. Follow the cleaning procedures provided in this guidebook.

E Using the AC Adapter

1. Check that the power switch is in the "OFF" position.

2. Connect the AC adapter to the AC adapter jack in the back of the unit.

3. Plug the AC adapter into the socket. (AC adapters with re-quired voltage and current indicated near the AC adapter jack.)

4. Use only the authorized AC Adapter with this nebulizer.

5. Press the power-switch to the "ON" position.

F Cleaning procedures

Rinsing technique (performed after each treatment or before first use):

1. Make sure that the power-switch has been turned to the "OFF" position and the unit has been disconnected from the power source.

2. Disconnect the air tube from the nebulizer device.

3. Gently twist and pull up the cover of the nebulizer kit to open and separate.

4. Rinse the nebulizer kit and components with hot tap water.

5. Dry clean components or completely dry air dry.




6. Reassemble the nebulizer kit.

NOTE: For the first time cleaning or after the unit has been stored for an extended period of time, thoroughly clean all components, except the air tube. The nebulizer kit is dishwasher safe.

G Cleaning the compressor

Wipe the compressor daily using a soft cloth.

NOTE: Any other form of cleaning or cleaning agents may damage the finish of the unit.

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