



QUICK START GUIDE

- Make sure you order the right-sized product from the Size Chart below
- Sit in a comfortable position and place your leg on a raised surface
- Place the product at the backside of the knee with wider side up
- Carefully match the axis of the hinge and the knee during application
- Ensure that the natural flexion movement (bending) of the knee joint is maintained
- Fasten the hook loop rejoiners on the fabric flap to a comfortable compression
- Fasten the hook loop bands from bottom to top to provide added grip
- Re-adjust in case of any discomfort



REMEMBER

- Do not roll or bend the support while using or storing it as it might decrease the elasticity and effectiveness

- Discontinue use and consult the doctor in case of impaired sensation, impaired blood flow, allergic reactions on skin and persistence or increase of pain

P.S.: Kindly consult your Physician or Orthopaedic Specialist before using this product