


How to Wear?

- Make sure you order the right-sized belt from the Size Chart below
- Stand in the neutral posture with your back straightened as recommended by the Physician
- Wrap the belt around the waist
- Keep the belt upright and in correct position as per directions of the “up inside” label
- Ensure that the lower edge of the belt aligns with the base of the spine
- Fasten the hook tape with the loop tape to a comfortable compression
- Re-adjust in case of any discomfort



REMEMBER

- Discontinue use and consult the doctor in case of impaired blood flow, allergic reactions on skin and persistence or increase of pain
- Consult the doctor in case of simultaneous use with any other medical device
- Improper use of the device may lead to adverse effects or decreased effectiveness. Follow the instructions on the packaging carefully before using
- Keep it out of reach of children

P.S.: Kindly consult your Physician or Orthopaedic Specialist before using this product