


How to Wear?

- Make sure you order the right-sized product from the Size Chart below
- Stretch the affected leg forward and place the patellar opening of the knee support over the patella (knee cap)
- Match the axis (central joint) of the hinge of the splint and the knee carefully
- Check that the natural flexion movement (bending) of the knee joint is maintained
- Fasten the hook and loop fasteners (the thigh first and then the calf) to a comfortable level of compression
- Re-adjust in case of any discomfort



REMEMBER

- Do not roll or bend the support while using or storing it as it might decrease the elasticity and effectiveness

- Discontinue use and consult the doctor in case of impaired sensation, impaired blood flow, allergic reactions on skin and persistence or increase of pain

P.S.: Kindly consult your Physician or Orthopaedic Specialist before using this product