


QUICK START GUIDE

- Make sure you order the right-sized product from the Size Chart below
- Sit comfortably and keep your leg straight
- Insert your feet into the kneecap while holding the top of it
- Wear it with a simple pull on application like socks, keeping broader end upwards
- Ensure whether the natural flexion movement (bending) of the knee joint is maintained
- Re-adjust in case of any discomfort



REMEMBER

- Do not roll or bend the support while using or storing it as it might decrease the elasticity and effectiveness

- Discontinue use and consult the doctor in case of impaired sensation, impaired blood flow, allergic reactions on skin and persistence or increase of pain

P.S.: Kindly consult your Physician or Orthopaedic Specialist before using this product