


# QUICK START GUIDE

- Make sure you order the right-sized product from the Size Chart below
- Sit in a comfortable position and place your leg on a raised surface
- Pull the product on the leg and position the patellar opening on the patella
- Ensure that the flexion movement of the hinge matches the natural flexion movement of the knee joint
- Align the fulcrum of the hinge carefully to the centre of the knee
- Fasten the hook loop rejoiners to a comfortable compression
- Remember to re-adjust in case of any discomfort



## REMEMBER

- Remove the support before sleeping
- Do not roll or bend the support while using or storing it as it might decrease the elasticity and effectiveness
- Discontinue use and consult the doctor in case of impaired sensation, impaired blood flow, allergic reactions on skin and persistence or increase of pain

- Consult the doctor in case of simultaneous use with any other medical device
- Misuse or improper use of the device can lead to adverse effects or reduced effectiveness. Please follow the instructions mentioned on the box carefully

**P.S.:** Kindly consult your Physician or Orthopaedic Specialist before using this product