


How to Wear?

- Make sure you order the right-sized product from the Size Chart below
- Sit in a comfortable position and place your leg on a raised surface
- Insert the leg into the support and pull up with the broader end facing upwards
- Position the patellar opening of the support on the patella
- Close the top and bottom straps to a comfortable compression
- Remember to re-adjust in case of any discomfort



REMEMBER

- Remove the support before sleeping
- Do not roll or bend the support while using or storing it as it might decrease the elasticity and effectiveness

- Discontinue use and consult the doctor in case of allergic reactions on skin and persistence or increase of pain
- Misuse or improper use of the device can lead to adverse effects or decreased effectiveness. Please follow the instructions mentioned on the box carefully

P.S.: Kindly consult your Physician or Orthopaedic Specialist before using this product