



QUICK START GUIDE

- Make sure you order the right-sized product from the Size Chart below
- Sit in a comfortable position and place your leg on a raised surface
- Place the product on the leg around the knee joint with the broader side facing upwards
- Position the patellar opening on the patella
- Ensure the flexion movement of the hinge matches the natural flexion movement of the knee joint
- Align the fulcrum of the hinge to the centre of the knee
- Fasten the hook loop rejoiners on the neoprene panel to a comfortable compression
- Fasten the hook loop straps on the bottom and upper side one by one for firm grip

- Remember to re-adjust in case of any discomfort



REMEMBER

- Remove the product before sleeping
- Do not roll or bend the support while using or storing it as it might decrease the elasticity and effectiveness
- Discontinue use and consult the doctor in case of impaired sensation, impaired blood flow, allergic reactions on skin and persistence or increase of pain

- Consult the doctor in case of simultaneous use with any other medical device
- Misuse or improper use of the product may lead to adverse effects or decreased effectiveness. Please follow the instructions mentioned on the box carefully

P.S.: Kindly consult your Physician or Orthopaedic Specialist before using this product