


How to Wear?

- Make sure you order the right-sized product from the Size Chart below
- Sit in a comfortable position and place your leg on a raised surface
- Place the support at the back of the knee with the wider side facing upwards
- Match the axis of the hinge with the knee. Ensure the flexion movement on the hinge of the splint matches the natural flexion movement of the knee
- Fasten the hook rejoinders on the fabric flap to a comfortable compression
- Secure the hook loop bands from bottom to top to provide added grip

- Tie the oblique elasticized straps for a comfortable counterforce
- Remember to re-adjust in case of any discomfort



REMEMBER

- Understand the direction of the counterforce carefully before application
- Remove the product before sleeping
- Do not roll or bend the support while using or storing it as it might decrease the elasticity and effectiveness
- Consult the doctor in case of simultaneous use with any other medical device

- Discontinue use and consult the doctor in case of impaired sensation, constricted/limited blood flow, allergic reactions on skin and persistence or increase of pain
- Misuse or improper use of the product may lead to adverse effects or decreased effectiveness. Please follow the instructions mentioned on the box carefully

P.S.: Kindly consult your Physician or Orthopaedic Specialist before using this product