


How to Wear?

- Make sure you order the right-sized product as per the Size Chart below
- Stand in the neutral posture with your back straightened as recommended by the Physician
- Place the special compression pad of the belt at the site of the rib fracture
- Wrap the belt around the waist
- Make sure that the belt is upright and in correct position as per directions of the “up inside” label
- Fasten the hook and loop to a comfortable compression
- Re-adjust in case of any discomfort



REMEMBER

- Discontinue use and consult the doctor in case of impaired blood flow, allergic reactions on skin and persistence or increase of pain
- Consult the doctor in case of simultaneous use with any other medical device
- Improper use of the device may lead to adverse effects or decreased effectiveness. Follow the instructions on the packaging carefully before using
- Keep it out of reach of children

P.S.: Kindly consult your Physician or Orthopaedic Specialist before using this product