 tynor

**ROM Elbow Brace**

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### Free Arm Movement:

- When free flexion movement of the elbow is required. Place flexion knob at 120 degree and extension knob at 0 degree.

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### Immobilization:

- When immobilization at a particular angle is required move both, flexion as well as extension locking pointer to the same desired angle on the graduated disc. (limited up to 90°).
- For Single action locking at 90°, 100° or 110° slightly fold the forearm so that the arrow mark comes exactly in front of the marked green arrow of desired angle, at the same time push the green locking pointer to the same angle to achieve quick immobilization.

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### Applying Optional Arm Sling:

- Apply the strap having length adjustment arrangement with D buckle, provided in wrist cuff strap.
- Bring the other loose end of the sling strap over the shoulder of same elbow side and under the opposite side arm, now hook with the same D buckle.
- Remove the Y clamp hook tape fastener to adjust the sling strap, as per required comfort level.
- The sling strap can be cut shortened, if required, by removing the Y clamp.

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# ROM Elbow Brace

English

**ROM Elbow Brace** controls and restricts degree of flexion and extension of the elbow joint through sophisticated dial system to allow a range of motion or immobilization of the elbow. Universal design allows it to be used on most of the patients and for multiple orthopedic problems.

## CONSTRUCTION:

ROM Elbow Brace is constructed from a strong aluminum frame with high quality polymer semi flexible cuffs padded with 3D molded cushion pads for upper arm and forearm. It has a ROM hinge which is a sophisticated multi point graduated mechanism (MPGM). The MPGM can be adjusted to fix any angle of flexion or provide a controlled motion between any two angles by a simple pull- push locking mechanism. Length of device can be independently adjusted, both on bicep and forearm side.

## APPLICATIONS:

- Mobilizing brace in Stiff elbow
- Post-operative immobilization
- Conservative treatment of elbow luxation
- Preventive injuries in Sports
- Ligamentous injuries.
- Stabilization of elbow and upper limb fractures.
- Collateral ligament reconstruction.
- Tennis elbow, lateral condyle release.
- Elbow hyper extension correction.
- Mild to severe soft tissue injuries.

## RECOMMENDATION:

- Do not use on unstable fractures.
- Not for use on patients with circulatory or sensory problems.
- Excess pull of built in Locking mechanism for length adjustment may damage it.

## INSTRUCTIONS:

- Use this product under strict guidance of a qualified doctor.
- Adjust padding and strapping for secure and comfortable fitting.
- Discontinue use and seek guidance of a qualified doctor, in case of
  - Impaired sensation
  - Impaired blood flow
  - Allergy rash or itching
  - If pain increases or persists
- Consult your doctor in case of concurrent use with another medical device.
- Misuse, or improper use of the device can lead to adverse effects and decreased effectiveness
- Follow instructions mentioned on the box carefully.
- Keep out of reach of children.

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English

### **POLICY FOR DEFECTIVE GOODS:**

All products manufactured by Tynor Orthotics (P) Limited carry one year warranty against material or manufacturing defect. Defective units must be returned to Tynor Orthotics (P) Limited through the Super Dealer concerned for replacement, along with duly filled " Product Return Form" . Liability of Tynor Orthotics (P) Limited shall be limited to the cost of product only.

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Hand wash with mild detergent and water below 30°C without wringing.

Dry in shade on a flat surface.

Do not dry-clean.

Do not iron.

Do not bleach.

Close hook loop fasteners before washing.



Wipe the plastic parts with moist lint free cloth and mild detergent

EC REP

MT Promedt Consulting GmbH  
 Altenhofstrasse 80  
 D-66386 St. Ingbert  
 Germany  
 tel.: +49 (0) 6894 581020  
 fax: +49 (0) 6894 581021  
 e-mail: info@mt-procons.com  
 www.mt-procons.com

**Symbol Explanation :**

-  Batch code
-  Manufacturing date
-  Manufacturer
-  Catalogue number