

How to Wear?

- Assistance by a caregiver is required while wearing the ROM Knee brace
- Hold the brace gently and fix the angle of flexion on the dial to the prescribed / required degree.
- To fix the angle of extension, pull the pointer knob softly on the dial, similar to fixing the angle flexion
- Fix the other dial similarly, at the same degree of flexion and extension
- Unfastening all the straps, place the brace on the patient's leg, aligning both the dials with the centre of the knee
- Ensure that the thigh side of the brace is in upward position
- Fasten the straps, starting from the cuff side and moving upwards
- In case of excess pressure, pain or discomfort at any point, re-adjust the tightness of the straps accordingly, or consult the doctor
- Ensure the brace is snugly fixed and provides the required degree of flexion or fixation

How to Remove?

Remove the velcro straps gently starting from thigh to knee

REMEMBER

- For a secure and comfortable fitting, adjust the padding and strapping accordingly
- Discontinue use and consult the doctor in case of impaired sensation, impaired blood flow, allergic reactions on skin and persistence or increase of pain
- Consult the doctor in case of simultaneous use with any other medical device
- Improper use of the device may lead to adverse effects or decreased effectiveness. Follow the instructions on the packaging carefully before using
- Keep it out of reach of children

P.S.: Kindly consult your Physician or Orthopaedic Specialist before using this product