


How to Wear?

- Make sure you order the right-sized product from the Size Chart below
- Insert foot into the sleeve of the thigh support
- Gently pull upwards until it covers the upper and middle thigh region
- Ensure comfortable compression and grip
- Re-adjust in case of any discomfort



REMEMBER

- Do not roll or bend the support while using or storing it as it might decrease the elasticity and effectiveness

- Discontinue use and consult the doctor in case of impaired sensation, impaired blood flow, allergic reactions on skin and persistence or increase of pain

P.S.: Kindly consult your Physician or Orthopaedic Specialist before using this product