

How to Wear?

- Assistance by a caregiver is required while wearing the Taylor's brace
- Make sure you order the right-sized brace from the Size Chart below
- The user should stand in an upright posture
- Place the lower horizontal splint just over the lower edge of the spine. Now mark the position of the vertical and horizontal splint on the body surface
- Shape the splints by bending as per the anatomy of the body indicated by the positions marked on the body
- Insert the customized splints into the brace
- Place the central soft panel of the back support in between two vertical splints exactly on the vertebral column. Ensure that the lower edge of the brace coincides with the nape of the spine
- Wrap the brace around the body
- Fasten the hook-loop closures securely to ensure a snug compression
- Secure the shoulder straps by throwing them over the shoulders and bringing back to the axial region
- Pull the straps through the buckles to tighten. Fasten with the loop tape provided vertically along the elastic panel

REMEMBER

- The vertical and horizontal splints must be removed, customized and re-inserted strictly by a professional
- Discontinue use and consult the doctor in case of impaired blood flow, allergic reactions on skin and persistence or increase of pain
- Consult the doctor in case of simultaneous use with any other medical device
- Improper use of the device may lead to adverse effects or decreased effectiveness. Follow the instructions on the packaging carefully before using
- Keep it out of reach of children

P.S.: Kindly consult your Physician or Orthopaedic Specialist before using this product