


How to Adjust Height?

- First identify the height you would want to adjust it to
 - Stand straight and let your arms hang by your side comfortably
 - Measure the distance from the ground to your inner wrist
 - This is your ideal walking stick height
 - You can always visit your Physician or Therapist to get a more ideal measure and recommendation in case of advanced health issues especially back pain or abnormal posture
- To choose right vs. left handed use, always opt for the **side opposite** to the affected limb. If you have pain in your left leg, align the quadripod for using it with the right hand



- Holding the stick in your lap or on a table, loosen the metal tension screw at the bottom
- Press the small metal spring-loaded pin on the side of the stick and press it inwards
- Push or pull the base to adjust the height of the stick as necessary

- If you want to shift it from **right-handed to left-handed use**, just flip outer tube around 180° at this point
- Once adjusted, simply align the pin with the groove you want it at
- Tighten the tension screw at the bottom again
- The ferrules can be changed or pulled out for cleaning easily. Please ensure they are well settled on the base before using the tripod daily

1



2



3



4



REMEMBER

- Please do not use the stick to lift, hold or bar weights such as luggage or as tool for defence
- The coating is prone to scratches. Please do not throw the quadripod around or misuse it for any other purpose

- If you notice damage to the body of the stick, please replace it