



OOPSIE HEROES

BEDWETTING ALARM NOTES



QUICK INSTRUCTIONS

- Download the Oopsie Bedwetting App from Google Play or the Apple App Store.
- Select the Oopsie **PLUS** + version.
- Follow the instructions on the App to register your details.
- Select an alarm sound.
- Remove a sticker from your sticker sets, press it firmly and accurately in place on the back of your sensor.
- Firmly press the sensor on the outside of your child's underpants. Place the sensor in the spot which will become wet first.
- Press **GO TO SLEEP** on your App.
- When the alarm sounds during the night, tap the App to switch off the alarm.
- Remove the sensor, wipe clean and reapply.
- In the morning **REMEMBER** to press the **WAKE UP** button on your App. This will preserve your sensor battery.

SETTING UP YOUR DEVICE

1. You can use any Apple or Android device. An old cell phone or tablet is suitable.
2. Download the Oopsie Bedwetting App. Select the Oopsie **PLUS** + version.
3. Follow the instructions on the App to register your details.
4. Select an alarm sound or record your voice alarm.
5. The device should be placed in the same room as the child, 2-3 metres from the bed.
6. Set the phone to silent/do not disturb or aeroplane mode. This will block notifications.
7. Turn the volume on the phone loud enough to wake your child.
8. You can lock the screen once you have activated the sleep-mode on the App. To turn off the alarm you will need to unlock the device. On your device you will have the ability to lock all apps except the Oopsie App. On an Apple device, look for the Guided Access function (under Settings/Control Centre) and on an Android device you can use App Pinning.
9. Always make sure your device is always fully-charged.
10. In order to avoid the device battery running flat, plug in the device at night.

CARING FOR YOUR SENSOR

- **DO NOT** submerge your sensor in water.
- Clean the sensor between uses with a damp tissue. If you do not do this, crystals can build up on the pins and stop the sensor from working.
- Each morning, wipe clean the sensor and place it in the box for safe-keeping.
- In the morning, **REMEMBER** to press the **WAKE UP** button on your App. This will preserve your sensor battery. The battery in the sensor has a life of 180 days of use.

DO'S & DONT'S

- **DO NOT** use a nappy. It will block the signal and you want your child to feel the wetting sensation.
- **DO NOT** submerge the sensor in water.
- Visit the FAQ link on the Oopsie product page on our website for more details.
- Before you start using your alarm, download our **FREE Wake Up Dry** e-book from our website www.moosebaby.co.nz or www.moosebaby.com.au. This will provide tips & tricks to success.

KEEP ME