

mini Air 360用前必☒
 成品尺寸：350mm*100mm(h)
 ☒折☒尺寸：70*100mm(h)
 128g双面☒板☒ ☒色：PANTONE:2587C
 -----延曲☒☒折

正面

VocalMist

Portable Nebulizer

Quick Start Guide

IMPORTANT

- VocalMist Portable Nebulizer is to be used with VocalMist Hydration Formula ONLY (sterile Isotonic Saline 0.9% NaCl).
- Please do not use distilled water, tap/Spring water, or oil substances for nebulization.
- Please empty and allow the VocalMist to air dry after each use.
- Please keep the exterior of the device, saline cup and accessories dry.
- Please do NOT touch the nebulizing mesh (inside the nozzle) with fingers, hard objects, or soft objects.

OPERATING STEPS

Following the arrows on top of the VocalMist, rotate the grey saline cup cap counterclockwise to open.

Add VocalMist Hydration Formula (Isotonic Saline) to your nebulizer. Please do not exceed the max volume indicator.

Reattach cap, and connect either the mask or mouthpiece.

Tap the green button to start nebulizing.

Use with mouthpiece for direct vocal fold hydration, or use with mask for treatment of both throat and sinuses.

AFTER EACH USE

Rotate the mask/mouthpiece to loosen, and then remove from device.

Open the saline cup and empty all remaining liquid.

Keep saline cup open to let nebulizer air dry completely after each use. Do not store liquid in nebulizer when not in use.

If purchased, store your VocalMist safely in your VocalMist Carrying Case between uses.

Your case has a dedicated slot for your grey cap, so the VocalMist can remain open and dry while in the case.

The simple air dry method described here is all you need to do after each use to keep your VocalMist clean and working well.

If the inside of your VocalMist cup becomes visibly dirty, please refrain from use and message us at support@myvocalmist.com for more involved disinfection instructions.

Do not remove the cup from the device unless directed to by our troubleshooting steps, as this may shorten the lifetime of your device.

These statements have not been evaluated by the Food and Drug Administration. This product is not intended to diagnose, treat, cure or prevent any disease.

Your Voice Will Thank You

myvocalmist.com | vocalmist

背面