PRODUCT CARE

Hand wash with mild detergent and water below 30°C without wringing.

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

English

INTENDED USE :

Contoured Cervical pillow is anatomically shaped and designed to support the cervical spine in neutral position during sleep. Contoured Cervical Pillow keeps the neck in slightly hyper extension position while sleeping. This counters the stress and strain of the day by relaxing the muscles, where we tend to keep the neck in the opposite "Flexion" position. It cures the cervical problems while we sleep. Corrects the sleeping posture and straightens the spine.

CONSTRUCTION:

Molded out of high quality PU foam to an optimized density to provide firm but comfortable support to the neck and head. The anatomical design with two contours two different degrees of hyper extension to the neck. Choosing the right option of the two improves the comfort.

APPLICATIONS:

- Cervical spondylosis
- Back pain
- To correct snoring induced apnea.
- Waking headache
- Scapular or arm pain.
- Mild Sprain or Stiff neck.
- As a prophylactic aid for cervical problems.
- Can be used on surgical tables, recovery rooms and treatment tables.

RECOMMENDATION :

- Do not sleep on stomach if you have neck pain or cervical problems
- The use of additional removable cover is recommend to keep the pillow neat & clean.

INSTRUCTIONS:

- Use this product under strict guidance of a qualified doctor.
- Keep it under neck while asleep
- Adjust and position the pillow so as to give the neck slight hyper extension position and comfort.
- Discontinue use and seek guidance of a qualified doctor, in case of
 - Impaired sensation
 - Pain increases or persists
- Consult your doctor in case of concurrent use with any other medical device.
- Misuse or improper use of the device can lead to adverse effects and decreased effectiveness.
- Follow instructions mentioned on the box carefully.