tynor

D.V.T. Stocking
D.V.T. Medias
D.V.T. Stockage
D.V.T. جورب

One Pair

REF I-66 & 69

1

HOW TO USE:

- Put on stockings first thing in the morning before you get out of bed. Your legs have the least amount of swelling early in the morning.
- Use a little baby powder or corn starch on your legs. This may help the stockings slide up.
- For open-toe stockings, use applicator to slide the stocking over your foot.
- Slide the stockings onto your arms with hand reaching the heel area.
- Holding the extra heel material, turn the stockings inside out till the heel.
- Broaden the mouth of the fabric cave so formed and pull till it reaches the heel.
- Unfold the stockings from heel to knee, in a manner that fabric rolls on the fabric and not on bare skin.
- Pull the upper part of the stockings from knee till half of mid thigh.
- Pull the upper part of the stockings from half of mid thigh till mid thigh.
- Smoothen the wrinkles from bottom upwards with both hands, gloved or otherwise till you get absolutely stretched garment on the leg.

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Medical compression stockings provide graduated compression, which is strongest at the ankle and it gradually decreases as we move up. This controlled and graduated compression effectively stops any abnormal back flow of blood. They compress the outer walls of superficial veins to make them strong and prevent varicose veins.

CONSTRUCTION:

Woven on the state of art circular looms, from high quality, Nylon and spandex yarn, to get a four way stretching seamless stockings with anatomical shape of the calf and a gradual taper from top to bottom.

The top end has dots of silicon to make it slippage proof.

COPING WITH VARICOSE VEINS:

- Avoid prolonged standing
- Take regular walks
- Keep leg elevated while sitting if possible.
- Wear Compression stockings before you get out of the bed while your leg is still elevated.
- Loose weight if overweight.
- Avoid constricting clothes such as girdles around the leg.

APPLICATIONS:

- Tired, aching legs
- Edema of ankle. Foot or leg
- Chronic peripheral venous insufficiency
- Varicose veins
- Spider Veins (mild varicosities)
- Deep Vein Thrombosis (DVT)
- Lymphedema
- Phlebitis
- Economy Class Syndrome (ECS)
- Lipodermatosclerosis
- Pregnancy
- Aching and heavy feeling in legs
- Blood clots, post surgery or injury.

CONTRAINDICATIONS:

Patients with peripheral obstructive arterial disease, peripheral neuropathy, congestive heart failure, septic phlebitis or oozing dermatitis should not wear Varicose Stockings.

CAUTION/INSTRUCTIONS:

- If you use lotion on your legs, let it dry before you put on the stockings.
- Use a little baby powder or corn starch on your legs. This may help the stockings slide up.
- Put on rubber dishwashing gloves. Use your hands to adjust the stocking and smooth it out.
- For open-toe stockings, use a special gadget to slide the stocking over your foot.

INSTRUCTIONS:

- Use this product under strict guidance of a qualified doctor.
- Discontinue use and seek guidance of a qualified doctor, in case of
 - Impaired sensation.
 - Impaired blood flow.
 - Allergy rash or itching.
 - Pain increases or persists.
- Consult your doctor in case of concurrent use with another medical device.
- Excessive physical exertion can lead to a risk of rise in blood pressure. Take extra precaution
- Misuse or improper use of the device can lead to adverse effects and decreased effectiveness.
- Follow instructions mentioned on the box carefully.
- Keep out of reach of children.

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POLICY FOR DEFECTIVE GOODS:

All products manufactured by Tynor Orthotics (P) Limited carry six months warranty against material or manufacturing defect. Defective units must be returned to Tynor Orthotics (P) Limited through the Super Dealer concerned for replacement, along with duly filled " Product Return Form" . Liability of Tynor Orthotics (P) Limited shall be limited to the cost of product only.

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PRODUCT CARE

Hand wash with mild detergent and water below 30°C without wringing.

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

Knee High

SIZE	INCHES	CENTIMETERS
S	9.8 to 11.8	25 to 30
M	11.4 to 13.4	29 to 34
L	13.0 to 15.0	33 to 38
XL	14.6 to 16.5	37 to 42


Measure circumference around the calf approximately 6 inches below the knee.

Thigh High

SIZE	INCHES	CENTIMETERS
S	14.6 to 16.9	37 to 43
M	16.5 to 18.9	42 to 48
L	18.5 to 20.9	47 to 53
XL	20.5 to 22.8	52 to 68

Measure circumference around mid thigh, approximately 6 inches above the knee.

CAUTION: Read user instructions before using the product

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Tynor Orthotics (P) Ltd.
TYNOR ORTHOTICS PVT.LTD