

PRODUCT CARE



Hand wash with mild detergent and water below 30°C without wringing.



Dry in shade on a flat surface.



Do not dry-clean.



Do not Iron.



Do not bleach.



Close hook loop fasteners before washing.

EC REP

MT Promed Consulting GmbH
Altenhofstrasse 80
D-66386 St. Ingbert
Germany
tel.: +49 (0) 6894 581020
fax: +49 (0) 6894 581021
e-mail: info@mt-procons.com
www.mt-procons.com

Symbol Explanation :



Batch code



Manufacturing date



Manufacturer



Catalogue number

Tynor Orthotics Pvt. Ltd.

169-170, Sec 82, Ind. Area, JLPL, Mohali-160 055, (Pb) India
e-mail: marketing@tynorindia.com • www.tynorindia.com
Customer Care : 0091-172-5288111

ISO
9001



W.H.O.
G.M.P.

ISO
13485



US FDA
registered

UM-D06-18-1

Knee Cap (With Rigid Hinge)

नी कैप (विद रिजिड हिंज)

Rodillera (Con Bisagra Rígida)

Rotuliere (Avec Charniere Rigide)

غطاء الركبة (ذو المفصل الميتين)



REF D 06

www.tynorindia.com


USER MANUAL

How To Use Cómo Utilizar



कैसे उपयोग करें Comment utiliser كيفية الاستخدام

1.



- Knee Cap (with Rigid Hinge) is worn by a simple pull on application like socks, with the broader end upwards.
- नी कैप (रिजिड हिंज) को मोजे की तरह छोड़े किनारे की उपर की ओर रख कर हलका खिचें ।
- La rodillera (con bisagra rígida) se coloca con un simple tirón a la manera de un calcetín, con su extremo más amplio hacia arriba.
- Rotulière (avec charnière rigide) est portée par une simple traction sur l'application comme les chaussettes, avec l'extrémité la plus large vers le haut.

● طیسب بحد للاحنم (نیٹملا لصفملا یذ) تیکرلا عاطغ اندترا متی
ی اعلایع سولاً فرطلا لعجم بروجلا اندترا لثم

2.



- Care should be taken to match the flexion movement of the hinge with the natural flexion movement of the knee joint.
- घुटना जोड़ के प्राकृतिक मोड़ संचालन के साथ हिंजो के संचालन मोड़ का मिलान करें ।
- Se debería prestar atención en que el movimiento de flexión de la bisagra coincida con el movimiento natural de flexión de la articulación de la rodilla.
- Des précautions doivent être prises pour adapter le mouvement de flexion de la charnière au mouvement de flexion naturel de l'articulation du genou.


● اءنڈنلا تکرجم عاطغلا لصفم اءنڈنلا تکرجم قاهاضملا تیانعلا بجز
تیکرلا لصفملا تیعیطلا

How To Use Cómo Utilizar



कैसे उपयोग करें Comment utiliser كيفية الاستخدام


3.



- The fulcrum of the hinge should be aligned to the center of the knee.
- हिंजो को आलंब को घुटने के केंद्र में संरेखित (अलाइंड) किया जाना चाहिए ।
- El punto de apoyo de la bisagra debería estar alineado con el centro de la rodilla.
- Le point d'appui de la charnière doit être aligné sur le centre du genou.

● يجب محاذاة نقطة ارتكاز المفصل مع مركز الركبة

4.



- The hook loop rejoiners are fastened to a comfortable compression by reverse buckle mechanism.
- पीछे की ओर बकल युक्ति के द्वारा हुक लूप टेप को आरामदायक दबाव तक बांधा जाता है ।
- Las cintas de sujeción se aprietan hasta proporcionar una compresión confortable mediante el mecanismo de hebilla invertida.
- Les répliques de boucle de crochet sont fixées à une compression confortable par un mécanisme à boucle inverse.

● يتم ربط أشرطة إعادة ضم حلقة الربط إلى حد الضغط المريح من خلال آلية الإيزيم العكسية.

English

Knee Cap (with Rigid Hinge)

INTENDED USE :

Knee Cap (with Rigid Hinge) is a simplified but effective design which offers the advantage of compression around the knee and the support of a side splint. It allows normal flexion of knee.

CONSTRUCTION:

Woven on the state of art circular looms to provide double layered, seam less, tapered made out of nylon and cotton yarn. Interspersed rubber thread, allows uniform compression even on uneven diameters of the limb.

Strong hinges provide lateral support to the injured tendons/ligaments of the knee, while keeping the flexion movement of the knee intact, they mimic the movement of the natural knee joint.

APPLICATIONS:

- Mild to moderate instability of knee joint and support to weak knee.
- Early stages of Osteoarthritis.
- Control painful knee movement in arthritis.
- Mild meniscus injuries.
- Mild ligament injury [medial or lateral].
- Mild to moderate ACL and PCL tears.
- Psychological assurance in Geriatric care.
- Post operative knee care.
- Post Partial meniscectomy or cartilage smoothing.
- Prophylaxis in sports activities.
- Mild strain and sprain.

RECOMMENDATION :

- Use this product under strict guidance of a qualified doctor.

INSTRUCTIONS:

- Discontinue use and seek guidance of a qualified doctor, in case of
 - Impaired sensation
 - Impaired blood flow
 - Allergy rash or itching
 - Pain increases or persists
- Misuse or improper use of the device can lead to adverse effects and decreased effectiveness.
- Follow instructions mentioned on the box carefully.
- Keep out of reach of children.