
- ① Wash hands and face thoroughly before using I-LID 'N LASH<sup>®</sup> PLUS. Apply a moist warm compress to closed eyelid in order to loosen encrusted debris and secretions.
- ② If you have sensitive eyes place a drop of I-DROP<sup>®</sup> PUR or I-DROP<sup>®</sup> PUR GEL on the eye to provide a protective barrier.
- ③ Open the container and remove one I-LID 'N LASH<sup>®</sup> PLUS wipe. When finished, close the container tightly.
- ④ Before applying the I-LID 'N LASH<sup>®</sup> PLUS, close the eye lightly, without squinting, and keep it closed until the application is complete.
- ⑤ Gently apply the gel from the wipe in a circular motion starting at the eyelash roots working outwards (please refer to Figure 1).
- ⑥ If necessary, clean the area around the lid margins and roots of the eyelashes more thoroughly by wiping in short strokes in the direction of the eyelashes (e.g. starting at the roots and going away from the eye for both lower and upper eyelashes).
- ⑦ Do not rinse. Any remaining I-LID 'N LASH<sup>®</sup> PLUS will continue to cleanse and hydrate the skin.
- ⑧ Repeat steps 2 through 6 for the second eye.
- ⑨ If you experience stinging, rinse the eye and surrounding area thoroughly with water.
- ⑩ Repeat this process with a new wipe if any crusty residue remains on the lash roots or lashes.
- ⑪ Perform this procedure once or twice per day for at least four weeks, or as per your doctor's instructions.

Start here for STEP 5

For STEP 6 use motion as below


**Fig. 1:** Apply the I-LID 'N LASH<sup>®</sup> PLUS wipe in a circular motion starting at the roots of the lashes and then out towards the lids and eyebrows.

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