Male Fertility Sperm Test

User Manual

-EN-

REF RH-SP10-C

INTENDED USE

Easy@Home Male Fertility Sperm Test is a rapid at-home test to detect the specific protein 10 (SP10) in semen qualitatively. This simple test will quickly let you know whether your sperm count is considered within normal range.

For *in vitro* diagnostic use only. For self-testing use.

This test is a quick screening test that will give you either a positive (normal sperm count) or negative (low sperm count) result. An explanation of how to read and interpret the test results is given in the "HOW TO READ THE RESULTS" Section. This test provides a preliminary screening only. Consult your physician with any results that are concerning or you have questions, particularly when preliminary screen tests show negative.

HOW DOES IT WORK

SP10 is a protein that is only expressed in testing and has been verified as an analyte for Sperm Concentration test. The concentration of SP10 is proportional to the concentration of sperm. When the sperm concentration is more than $15 \times 10^6 / \text{mL}$, the result is positive. This sperm concentration test can be used as an aid in diagnosis and treatment observation of male infertility.

When the sample is added into the specimen well, capillary action carries the specimen to migrate along the membrane. After administering the test, the C line should be present to indicate the test is valid. The absence of a C line at the read time would indicate an invalid result. If invalid, test again with the test cassette.

It may indicate the oligozoospermia and male factor infertility. For further detail you should consult your physician. A negative result does not mean that a pregnancy cannot occur and it should not be used as a reason for avoiding use of contraception.

7. Do not open the test cassette sealed pouch until specimen is collected and ready to be tested.

8. Do not touch the result window and sample well of the cassette to avoid contamination.

9. Do not mix different samples and dilution of different batches. The sample should be liquefied; then the semen liquefaction process should be no more than 60 minutes. If the sample is not completely liquefied within 60 minutes, it is considered abnormal.

10. The results should be interpreted by physicians along with clinical features and other laboratory test results. Any questions or suggestions, please contact the manufacturer.

QUESTIONS & ANSWERS

1. What should I do if the test result is invalid?

If invalid, you should retest with a new test kit. Positive test result cannot prove fertility but only prove that sperm concentration in the semen sample is not less than 15 million per milliliter, which is below that of most fertile men. However, the result of this test alone does not prove infertility, as there are other factors affecting fertility. However, some men with sperm count below this level are still able to father children naturally. It is also important to understand that your sperm count can vary from day to day. If there is any suspicion of the test result, it is recommended to collect sample again after an interval of more than seven days but less than three weeks to retest. If the result is still negative, it is recommended to consult your physician.


4. My semen sample did not become a thin liquid after 60 minutes. Can I still perform the test?

Semen samples will not liquefy as quickly or as fully as others. Even if your sample is not completely liquefied, the test kit may still show an accurate result. However, the premise is that you have leave the semen sample standing for at least 20 minutes and mixed it thoroughly as instructed, using any solid material when adding semen to the Sample Dilution Tube. If your sample has not liquefied at all or if your dropper gets clogged with solid or stringy material, you should discard the sample and collect another one 48 hours later. If the same problem occurs again, please consult your physician for a semen evaluation.

3. What does it mean if the test result is negative?

Negative results show that the sperm concentration $< 15 \times 10^6 / \text{mL}$.

a new test. When the level of sperm in specimen is at or above $15 \times 10^6 / \text{mL}$, there is a visible colored band in the Test Region(T) and indicates a positive result. When the level of sperm in specimens are lower than $15 \times 10^6 / \text{mL}$, absence of this colored band in the Test region (T) suggests a negative result.

CONTENTS OF THE TEST KIT

Materials Provided

One test requires:

1. One test pouch. One pouch contains a test cassette and a desiccant. The desiccant is for storage purposes only, and is not used in the test procedures.
2. A Sperm Transfer Device
3. A Sample Dilution Tube with 2 mL buffer solution
4. A Semen Collection Cup

Package insert and user instruction.

Materials Required but Not Provided

Timer (watch or clock)

When the sample is added into the specimen well, capillary action carries the specimen to migrate along the membrane. After administering the test, the C line should be present to indicate the test is valid. The absence of a C line at the read time would indicate an invalid result. If invalid, test again with the test cassette.

5. The test cassette is a disposable one-time use test. Do not reuse the test cassette.

6. Package contents may be disposed of in normal household waste after use.

7. Do not open the test cassette sealed pouch until specimen is collected and ready to be tested.

8. Do not touch the result window and sample well of the cassette to avoid contamination.

9. Do not mix different samples and dilution of different batches.

The sample should be liquefied; then the semen liquefaction process should be no more than 60 minutes. If the sample is not completely liquefied within 60 minutes, it is considered abnormal.

10. The results should be interpreted by physicians along with clinical features and other laboratory test results. Any questions or suggestions, please contact the manufacturer.

LIMITATIONS

1. This test kit cannot be used to identify paternity DNA testing OR as an indicator for contraception.

2. This test does not protect against sexually transmitted diseases.

3. This is a qualitative assay and the results cannot be used to determine the exact number of sperm.

4. It does not detect all fertility issues. As with all diagnostic tests, the test result must always be consistent with clinical findings.

5. The number of sperm varies widely so conclusion cannot be made after only one time of the test. Once abnormal results appear, the sample should be collected again for test after an interval of 2-7 days. Usually, 2 to 3 times re-tests are recommended.

6. This test kit is only used to evaluate sperm concentration through the detection of sperm marker protein SP10. The lowest concentration at which the hook effect appears is 423 million/mL.

7. Sperm concentration is only one factor that could impact a man's fertility status and time to pregnancy. For a complete assessment of male reproductive health, the user should consult a physician.

8. The following situations may cause inaccurate results.

- (1) Not collecting the entire ejaculate, especially the first drops of the ejaculate.
- (2) The semen sample was collected less than 48 hours or more than seven days after your last ejaculation. This may lead to inaccurate results.
- (3) Failure of ejaculate liquefaction.
- (4) The semen in the Sample Collection Cup was not mixed well enough before it was added to the Sperm Solution Tube.

PERFORMANCE CHARACTERISTICS

- Analytical Sensitivity (Limit of Detection): $15 \times 10^6 / \text{mL}$

- Analytical Specificity:

Although a number of drugs, disease organisms, and non-semen body fluids that might be found in the male reproductive tract have been tested and found no cross-reaction with the test, other substances or medical conditions that have not been tested may have an effect on the test results. Since illness can also temporarily reduce your sperm count, you should not take the test if you are (or recently have been) sick.

9. WHO Laboratory Manual for Human Semen (Version 5).

10. National Institute for Health and Care Excellence, NICE: Fertility: assessment and treatment for people with fertility problems. NICE Clinical Guidelines. 2013:274-.

11. WHO: Examination and Processing of Human Semen. Volume Edition, F; 2010:286.

12. Hirsch A. Male subfertility. BMJ. 2003;327(7416):669-72.

13. Counsel M, Bellinger R, Burton P. Vitality of oligozoospermic semen samples is improved by both swim-up and density gradient centrifugation before cryopreservation. J Assist Reprod Genet. 2004;21(5):137-42.

14. Rajmil O, Fernández M, Rojas-Cruz C, Sevilla C, Musquera M, Ruiz-Castañeda E. [Azospermia should not be given as the result of vasectomy]. Arch españoles Urol 2007; 60:55-8.

15. Dhar NB, Bhatt A, Jones JS: Determining the success of vasectomy. BJU Int 2006; 97:773-6.

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