BioCane™ Type CK50900

Cryobiological Storage Vessel

Operation Manual and Parts List LT509X17 (7000807) Rev. 10

Visit us online to register your warranty
www.thermoscientific.com/labwarranty

Thermo
SCIENTIFIC

Catalog Number	Model Number	Name
CK509X2	807	BioCane 20
CK509X3	808	BioCane 34
CK509X4	809	BioCane 47
CK509X6	820	BioCane 73

Manual Number LT509X17 (7000807)

10	ECNM-000831	02/25/21	Removed Warranty caution from preface page ii and 1-1. Added Safety Precautions in Preface page iii, WEEE Compliance in page 6-1 and removed EUT statements in page 2-1. Updated Specifications for Environmental Conditions in page 2-1.	TFI
9	41864	11/07/17	Removed CE mark, number	bpg
8	41238	12/20/16	Combined languages to one file	slb
7	40578	2/23/16	Clarified Do Not Roll warning	ccs
6	--	10/27/15	Updated symbols and labels, added cautions	ccs
5	30559	11/12/13	Added warnings and cautions to pg 1-2, environmental conditions to pg 2-1	ccs
4	27797	9/21/11	Added manufacturer icon and euro rep	ccs
3	27547	7/28/11	Updated specs - BioCane 73 vials from 2880 to 2016	ccs
2	27457/SI-10618	7/21/11	Added CE number to pg 2-1	ccs
1	--	2/24/11	Added part numbers for replacement canisters	ccs
0	--	5/28/10	Transfer to Marietta (was LT509X17 8/12/09)	ccs
REV	ECR/ECN	DATE	DESCRIPTION	By

Important Read this instruction manual. Failure to read, understand and follow instructions in this manual may result in damage to the unit, injury to operating personnel, and poor equipment performance. ▲

Caution All internal adjustments and maintenance must be performed by qualified service personnel. ▲

Warning Nitrogen gas can cause suffocation without warning. Store and use liquid nitrogen only in a well-ventilated place. As the liquid evaporates, the resulting gas displaces the normal air in the area. (The cloudy vapor that appears when liquid nitrogen is exposed to the air is condensed moisture, not the gas itself. The issuing gas is invisible.) In closed areas, excessive amounts of nitrogen gas reduces the concentration of oxygen and can result in asphyxiation. Because nitrogen gas is colorless, odorless, and tasteless, it cannot be detected by the human senses. Breathing an atmosphere that contains less than 19.5% oxygen can cause dizziness and quickly result in unconsciousness and death. Therefore, the use of oxygen monitoring equipment is strongly recommended. ▲

Warning Contact of liquid nitrogen with the skin or eyes may cause serious (freezing) injury. ▲

Warning Do not roll. Always keep liquid nitrogen vessel in an upright position. Exercise caution to prevent spilling and splashing. ▲

Material in this manual is for information purposes only. The contents and the product it describes are subject to change without notice. Thermo Fisher Scientific makes no representations or warranties with respect to this manual. In no event shall Thermo be held liable for any damages, direct or incidental, arising out of or related to the use of this manual.

Intended Use:

The Thermo Fisher Scientific cryogenic storage vessels are intended to store biological samples at cryogenic temperatures. These products are comprised of a vessel and a lid, with or without a level monitor.

Unintended Use:

These products are not intended for use within the patient environment.

When this manual is translated into other languages, the US English version is binding.

© 2021 Thermo Fisher Scientific. All rights reserved. Printed in U.S.A.

Important operating and/or maintenance instructions. Read the accompanying text carefully.

Extreme temperature hazards. Only qualified persons should perform procedures associated with this symbol.

Potential biological hazards. Proper protective equipment and procedures must be used.

This symbol when used alone indicates important operating instructions which reduce the risk of injury or poor performance of the unit.

CAUTION: This symbol, in the context of a CAUTION, indicates a potentially hazardous situation which if not avoided could result in minor to moderate injury or damage to the equipment. ▲

CAUTION: This indicates a situation which may result in property damage. ▲

Before installing, using or maintaining this product, please be sure to read the manual and product warning labels carefully. Failure to follow these instructions may cause the product to malfunction, which could result in injury or damage.

WARNING: Unauthorized repair of your freezer will invalidate your warranty. Contact Technical Service at 1-800-438-4851 for additional information. ▲

Do You Need Information or Assistance on Thermo Scientific Products?

If you do, please contact us 8:00 a.m. to 6:00 p.m. (Eastern Time) at:

1-740-373-4763	Direct
1-800-438-4851	Toll Free, U.S. and Canada
1-877-213-8051	FAX
http://www.thermofisher.com	Internet Worldwide Web Home Page
service.led.marietta@thermofisher.com	Tech Support Email Address
www.unitylabservices.com	Certified Service Web Page

Our **Sales Support** staff can provide information on pricing and give you quotations. We can take your order and provide delivery information on major equipment items or make arrangements to have your local sales representative contact you. Our products are listed on the Internet and we can be contacted through our Internet home page.

Our **Service Support** staff can supply technical information about proper setup, operation or troubleshooting of your equipment. We can fill your needs for spare or replacement parts or provide you with on-site service. We can also provide you with a quotation on our Extended Warranty for your Thermo Scientific products.

Whatever Thermo Scientific products you need or use, we will be happy to discuss your applications. If you are experiencing technical problems, working together, we will help you locate the problem and, chances are, correct it yourself...over the telephone without a service call.

When more extensive service is necessary, we will assist you with direct factory trained technicians or a qualified service organization for on-the-spot repair. If your service need is covered by the warranty, we will arrange for the unit to be repaired at our expense and to your satisfaction.

Regardless of your needs, our professional telephone technicians are available to assist you Monday through Friday from 8:00 a.m. to 6:00 p.m. Eastern Time. Please contact us by telephone or fax. If you wish to write, our mailing address is:

Thermo Fisher Scientific (Asheville) LLC
401 Millcreek Road, Box 649
Marietta, OH 45750

International customers, please contact your local Thermo Scientific distributor.

Table of Contents

Section 1	Safety Information	1-1
	Please Read	1-1
Section 2	General Specifications	2-1
	Unpacking	2-1
Section 3	Operation	3-1
	Measuring Liquid Nitrogen Quantity	3-2
Section 4	Replacement Parts	4-1
Section 5	Warranty Information	5-1
Section 6	WEEE Compliance	6-1

Section 1 Safety Information

This manual contains important operating and safety information. The user must carefully read and understand the contents of this manual prior to the use of this equipment.

Your Thermo Scientific cryogenic vessel has been designed with function, reliability, and safety in mind. It is the user's responsibility to install it in conformance with local electrical codes. For safe operation, pay attention to the warnings and cautions throughout the manual.


Please Read

All cryopreservation vessels, regardless of who manufactures them, use a vacuum to provide for the super insulative properties needed to keep nitrogen in a liquid form. They are very similar, in fact, to the thermos vessels you may have used to store coffee, soup or milk.

Remember how fragile they were? Since no vacuum is perfect or will last forever, we suggest that you monitor the consumption of liquid nitrogen used by your vessel on a regular basis, i.e. every 1-2 days. By monitoring, you may be able to anticipate subtle changes in consumption and possibly react to problems before they arise. You may opt to use a level monitor for this purpose.

Also, if your samples are super-critical, one-of-a-kind, or irreplaceable, consider allocating your samples into 2 separate vessels.

Remember the old adage about placing all your eggs into one basket...

Warnings

- This unit is intended to be used by trained personnel in a well-ventilated environment.
- Liquid nitrogen is extremely cold; it boils at -196°C .
- Follow all necessary precautions and warnings relevant to applicable hazards such as frostbite, contact burn, asphyxiation, etc., in cryogenic areas.

- To avoid injury due to frostbite or ruptured vials, use extreme care whenever handling liquid nitrogen, liquid nitrogen storage or transfer vessels or any objects which have come in contact with liquid nitrogen.
- Leave no areas of skin exposed.
- Wear adequate PPE where necessary: face shield, cryogenic gloves, cryogenic apron.
- Only use the power supply provided with the unit.
- Do not tightly seal liquid nitrogen containers or prevent nitrogen gas from escaping.
- Always handle vials with tongs.
- Use extreme care to prevent spilling and splashing liquid nitrogen during transfer and removal of storage contents and holders.
- Immediately remove any clothing or safety attire on which liquid nitrogen has been spilled or splashed .
- Get immediate medical attention for any frostbite injuries due to liquid nitrogen.
- Never overfill liquid nitrogen vessels.
- Use this product only in the way described in the product literature and in this manual. Before using it, verify that this product is suitable for its intended use. If the equipment is used in a manner not specified by the manufacturer, the protection provided by the equipment may be impaired.

When inserting or removing cans, be careful not to come in contact with the neck tube area of the vessel. Remove or insert cans slowly, in a vertical manner. Scratches on the neck tube can cause premature vacuum failure.

Do not tamper with, or remove vacuum port (covered by black plastic cap on side of vessel); this will void warranty.

Do not spill liquid nitrogen on vacuum port - this can cause vacuum failure.

Always keep liquid nitrogen in an upright position.

Section 2 General Specifications

	BioCane 20	BioCane 34	BioCane 47	BioCane 73
Height (Inches)	25.70"	26.60"	26.50"	27.00"
Diameter (Inches)	14.50"	18.20"	20.00"	22.00"
LN2 Capacity	20.50 L	34.80 L	47.40 L	73.0 L
Static Holding Time (Days)	227	193	121	121
Static LN2 Evaporation Rate*	0.1 L/day	0.18 L/day	0.40 L/day	0.60 L/day
Liquid Phase Capacity 1.2 ml	180 vials	720 vials	1512 vials	2016 vials

** Actual working time may vary due to current atmospheric conditions, container history, manufacturing tolerances, and any individual patterns of use.*

Environmental Conditions

Operating temperature: -29°C to +60°C

10% to 95% relative humidity, non-condensing.

Altitude limit: 2,000 meters.

Storage: -29°C to +60°C; 10% to 95% relative humidity.

Intended for indoor use only.

Unpacking

Before using your new vessel, carefully inspect it. Check for signs of damage which may have occurred in shipment. It is advisable to fill (see Filling Instructions) all new units with liquid nitrogen and watch liquid nitrogen loss rate for a few days. If there are any problems, call Technical Services as soon as possible.

The on-line warranty registration address is listed on the front of this manual. The unit should be registered within 30 days to make warranty effective. This information must include the serial number which is located on the handle of the vessel.

Note The most prevalent cause of failure of liquid nitrogen storage vessels is mechanical. The vessel neck tube supports the full weight of the inner shell and all liquid nitrogen it contains. A side or corner blow to the vessel causes the inner shell to swing in a pendulum motion causing the neck tube to be damaged. Any storage vessel which has been exposed to an accident, has been dropped or lowered to hit on one corner will tend to fail more rapidly than one which has not.

Caution Exercise caution when moving your BioCane™ vessel. BioCane cryobiological storage systems are not transportation vessels. Transport carts are designed for mobility within the lab or lab to lab only. Moving full vessels long distances, over cracks in floor, thresholds, on inclined ramps or in elevators can cause premature vacuum failure.

Note If samples must be transported under cryogenic conditions, consider the Thermo Scientific Arctic Express line of Dry Shippers.

Section 3 Operation

To avoid damage to your cryogenic storage vessel which may result in premature vacuum loss, it is important that the following procedure be used during the addition of liquid nitrogen to a warm vessel:

1. Add only a small amount of liquid nitrogen (2 - 10 liters) to new or warm vessels.
2. Allow this small amount of liquid nitrogen to sit in the covered vessel for a minimum of 2 hours. This will limit stress caused by the sudden temperature change associated with adding liquid nitrogen to a warm vessel.
3. Add an additional 2 - 10 liters LN₂ to vessel.
4. Allow vessel to sit for 48 hours and monitor liquid nitrogen consumption.
5. Fill vessel as desired. Remember to allow for displacement of liquid nitrogen when cannisters and canes are inserted.
6. Insert and remove cannisters slowly. Allow liquid nitrogen to run out of cannisters.

Caution Never overfill liquid nitrogen vessels.

When inserting or removing cans, be careful not to come in contact with the neck tube area of the vessel. Remove or insert cans slowly, in a vertical manner. Scratches on the neck tube area can cause premature vacuum failure.

Do not spill liquid nitrogen on vacuum port - this can cause vacuum failure. ▲

Measuring Liquid Nitrogen Quantity

1. Use a wooden yardstick to measure liquid nitrogen level. Never use a hollow tube or plastic dipstick to measure liquid nitrogen level.
2. Level will be indicated by frost line which develops when dipstick is removed and waved in a back and forth motion away from the user.

Section 4 Replacement Parts

Part Number	Description
AY509X13	Nalge Cryovials 1.2 ml, 500/case
AY509X12	Nalge Cryovials 2.0 ml, 500/case
AY509X20	CryoCane (aluminum) - each cane holds 5 vials for storage in BioCane cannisters (48/case)
AY509X21	CryoSleeve enclosures - Clear plastic sleeve fits over CryoCane for added security during handling or storage (100/case)
CS509X17	Nalgene System 100-ampule box accommodates most internally threaded 1.2 and 2.0 ml vials
CS509X18	Nalgene System 100-ampule box accommodates most internally threaded 1.2 and 2.0 ml vials
AY509X32	Nalgene System 100 Cryogenic Vials - 1 ml (500/case)
AY509X33	Nalgene System 100 Cryogenic Vials - 1.5 ml (500/case)
AY509X34	Nalgene System 100 Cryogenic Benchtop Vial Holder (1/package)
AY509X35	Nalgene System 100 Cryogenic Benchtop Vial Holder (1/package)
CryoCane coders - Aluminum identification tabs that slip over the end of a CryoCane to identify the cane written on for further identification (100/case)	
AY509X22	White
AY509X23	Yellow
AY509X24	Blue
AY509X25	Green
AY509X26	Red
Cryovial closure color coders - Plastic disks that fit into the top of Nalge Cryovials to identify individual (100/package)	
AY509X27	White
AY509X28	Yellow
AY509X29	Blue
AY509X30	Green
AY509X31	Red

Part Number	Description
HR509X25	BioCane 20 (CK509X2) [807]) replacement canister
HR509X26	BioCane 34 (CK509X3 [808]) replacement canister
HR509X27	BioCane 47 (CK509X4 [809]) replacement canister
HR509X28	BioCane 73 (CK509X6 [820]) replacement canister

Ordering Procedures

Refer to the Specification Plate for the complete model number, serial number, and series number when requesting service, replacement parts or in any correspondence concerning this unit.

All parts listed herein may be ordered from the Thermo Scientific dealer from whom you purchased this unit or can be obtained promptly from the factory. When service or replacement parts are needed, check first with your dealer. If the dealer cannot process your request, then contact our Technical Services Department.

Prior to returning any materials, contact our Technical Services Department for a “Return Materials Authorization” number (RMA). Material returned without an RMA number will be refused.

THERMO FISHER SCIENTIFIC STANDARD PRODUCT WARRANTY (LN₂ Vacuum)

The Warranty Period starts two weeks from the date your equipment is shipped from our facility. This allows for shipping time so the warranty will go into effect at approximately the same time your equipment is delivered. The warranty protection extends to any subsequent owner during the first year warranty period.

During the first year, component parts proven to be non-conforming in materials or workmanship will be repaired or replaced at Thermo's expense, labor included. *LN₂ Vacuum Integrity is covered for five years.* Installation and calibration are not covered by this warranty agreement. The Technical Services Department must be contacted for warranty determination and direction prior to performance of any repairs. Expendable items, glass, filters and gaskets are excluded from this warranty.

Replacement or repair of components parts or equipment under this warranty shall not extend the warranty to either the equipment or to the component part beyond the original warranty period. The Technical Services Department must give prior approval for return of any components or equipment. At Thermo's option, all non-conforming parts must be returned to Thermo postage paid and replacement parts are shipped FOB destination.

THIS WARRANTY IS EXCLUSIVE AND IN LIEU OF ALL OTHER WARRANTIES, WHETHER WRITTEN, ORAL OR IMPLIED. NO WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE SHALL APPLY. Thermo shall not be liable for any indirect or consequential damages including, without limitation, damages relating to lost profits or loss of products.

Your local Thermo Sales Office is ready to help with comprehensive site preparation information before your equipment arrives. Printed instruction manuals carefully detail equipment installation, operation and preventive maintenance.

If equipment service is required, please call your Technical Services Department at 1-800-438-4851 (USA and Canada) or 1-740-373-4763. We're ready to answer your questions on equipment warranty, operation, maintenance, service and special application. Outside the USA, contact your local distributor for warranty information.

Rev. 0 5/10

Section 5

Warranty Information

Section 6 WEEE Compliance

WEEE Compliance

WEEE Compliance. This product is required to comply with the European Union's Waste Electrical & Electronic Equipment (WEEE) Directive 2012/19/EU. It is marked with the following symbol. Thermo Fisher Scientific has contracted with one or more recycling/disposal companies in each EU Member State, and this product should be disposed of or recycled through them. Further information on our compliance with these Directives, the recyclers in your country, and information on Thermo Scientific products which may assist the detection of substances subject to the RoHS Directive are available at www.thermofisher.com/

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Spanish

thermoscientific.com

© 2021 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries. Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details.

Thermo Fisher Scientific (Asheville) LLC

401 Millcreek Road
Marietta, Ohio 45750
United States

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