## **User manual / *Manuel d'utilisation***

**INTEGRAL HSV Vitrification Unit  
for CBS™ High Security Vitrification straws /  
*Bac de vitrification INTEGRAL HSV pour  
palette de Vitrification Haute Sécurité CBS™***

**Read this manual carefully before using the INTEGRAL HSV Vitrification Unit  
*Lire attentivement ce manuel avant d'utiliser le Bac de Vitrification INTEGRAL HSV***

Contains: EN (p. 3), FR (p. 15), DA (p. 27), DE (p. 39),  
ES (p. 51), FI (p. 63), IT (p. 75), NL (p. 87), SV (p. 99).

Edition: 10/04/2018

**USER MANUAL**

Read this manual carefully before using the INTEGRAL HSV Vitrification Unit

## Contents

1	WARNING .....	4
2	INTRODUCTION .....	4
	<b>2.1 Definition</b> .....	4
	<b>2.2 Applications</b> .....	4
	<b>2.3 Protection and safety</b> .....	5
3	INTEGRAL HSV VITRIFICATION UNIT AND COMPONENTS .....	6
	<b>3.1 General view</b> .....	6
	<b>3.2 Cover</b> .....	7
	<b>3.3 Components and spare parts</b> .....	7
4	FIRST USE.....	7
5	OPERATING METHOD .....	8
	<b>5.1 CBS™ High Security Vitrification straws</b> .....	8
	<b>5.2 Opening CBS™ High Security Vitrification straws</b> .....	11
6	MAINTENANCE.....	13
	<b>6.1 Emptying the unit</b> .....	13
	<b>6.2 Cleaning</b> .....	13
7	TECHNICAL DATA .....	13
	<b>7.1 General</b> .....	13
	<b>7.2 Operating conditions</b> .....	13
	<b>7.3 Working life</b> .....	13
8	IMPORTANT INFORMATION.....	14
	<b>8.1 Ordering</b> .....	14
	<b>8.2 Warranty conditions for new equipment</b> .....	14
	<b>8.3 Waiving of liability</b> .....	14
	<b>8.4 Contact</b> .....	14

# 1 Warning


Before using the equipment, we recommend reading this user manual carefully, paying special attention to the sections marked with this symbol .

The INTEGRAL HSV Vitrification Unit must be used exclusively for vitrifying and heating CBS™ High Security Vitrification straws. It must be used exclusively with liquid nitrogen by trained, qualified staff.

This user manual describes the various components of the INTEGRAL HSV Vitrification Unit and gives the corresponding references to be used when ordering spare parts (see section 8.1: Ordering).

## 2 Introduction

### 2.1 Definitions

	This product must be disposed of in a suitable recovery and recycling centre
	Name and address of manufacturer
	Read the entire manual carefully before using the product
	Medical device marked as CE by the LNE / G-MED
	Risk of cold burn
	Product reference
	Product lot number
	Date of manufacture

### 2.2 Applications

The INTEGRAL HSV Vitrification Unit is intended for vitrifying, conditioning and opening CBS™ High Security Vitrification straws in the following applications:

- Reproductive biology laboratories
- Gamete banks.

## 2.3 Protection and safety

- The device complies with the "medical devices" directive 93/42/ CEE (class IIa). The conformity certificate and CE marking certificate can be obtained from <https://www.cryobiosystem-imv.com/> or on request.
- It is not subject to safety directives 2014/35/EU and 2014/30/EU.

### a. Installation, precautions and recommendations

1. This equipment must be used only by qualified staff trained in the use of cryogenic fluids such as liquid nitrogen.
2. Using this equipment may cause cold burns.
3. This equipment must be used only with liquid nitrogen.
4. This equipment must be used only in a room that meets safety criteria with regard to the use of liquid nitrogen.
5. This equipment can only be used in safety on a stable horizontal work bench or table, or a stable laboratory trolley.
6. This equipment must not be transported or moved when the chamber is filled with liquid nitrogen or if it contains samples.
7. This equipment must never be used if it is cracked.
8. To ensure complete safety, it is necessary to use the CBS™ High Security Vitrification straw opening device ref. 027636 or 027172 and a CBS™ perforated goblet with visotubes, for example of the daisy configuration type, ref. 015144.

### b. Transport and storage

The unit and its accessories must be stored and transported in the original packaging.


The unit must not be knocked or shaken.

The unit must be stored at ambient temperature in a dry place (max. temperature range: -20°C/ +60°C)

Maximum humidity: 70%.

### 3 INTEGRAL HSV Vitrification unit and components

#### 3.1 General view

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<b>A</b>	Movable, tilting goblet support	<b>G</b>	Outer shell
<b>B</b>	Cover	<b>H</b>	Support for straws awaiting vitrification
<b>C</b>	Liquid nitrogen level indicator	<b>I</b>	Vitrification block
<b>D</b>	Handle	<b>J</b>	Open straw holding area
<b>E</b>	Accessory storage rack	<b>K</b>	Tilting support for opening device
<b>F</b>	VHS straw platform	<b>L</b>	Chamber for liquid nitrogen

## 3.2 Cover

 No straw must be inside the INTEGRAL HSV Vitrification Unit when the cover is positioned.

Similarly, vitrification blocks must be kept in the lowered position.

## 3.3 Components and spare parts

Description	Quantity	References
Unit	1	N/A
VHS straw platform	1	N/A
Accessory storage rack	1	N/A
Goblet support	1	N/A
Goblet	1	015144
Vitrification block	3	027480
Opening device support	1	N/A
Level indication float	1	N/A
Cover	1	N/A
Support for straws awaiting vitrification	1 (neutral)	027622 (neutral), 027623 (white), 027631 (pink), 027632 (yellow), 027633 (green), 027634 (blue), 027635 (red)
Opening device	2	027636 (batch of 5)
Holder	1	027172 (opening set)
Tweezers	1	026923
Cover support	2	027628

## 4 First use

Place the INTEGRAL HSV Vitrification Unit on a stable horizontal table or laboratory trolley, in a room with a ventilation system and safety system adapted to the use of a cryogenic fluid such as liquid nitrogen.

Fill the unit with liquid nitrogen up to the maximum level indicated by the float (refer to section 5.1), observing the usual protocol in force in the laboratory.


Nitrogen level insufficient to complete the protocol

Minimum nitrogen level

Nitrogen level sufficient to complete the protocol

Maximum nitrogen level

## 5 Operating method

**NB:** The photos below were taken without any liquid nitrogen in the chamber to make it easier to understand the various stages involved.

### 5.1 CBS™ High Security Vitrification straws

1. Before filling the unit with liquid nitrogen, check that the following are in place:
  - the 3 vitrification blocks in their housings **in the lowered position**
  - the goblet support in its housing
  - a goblet in the goblet support
  - an opening device in the opening area
  - an opening device and a pair of tweezers in the storage rack
  - a support for straws awaiting vitrification
  - an ultrasonic sealer to seal the CBS™ High Security Vitrification straws in the work area

 Fill the unit with liquid nitrogen starting with the goblet and make sure it is **well filled**. If it is not to be used immediately, place the cover to reduce evaporation and thus maintain a sufficient level of nitrogen (cf. section 4).

 The equipment must not be transported or moved when the chamber is filled with liquid nitrogen or if it contains any samples.

2. Position a straw support on the edge of the unit (6 optional colors are available).
3. Label the CBS™ High Security Vitrification straws and place them on the straw support, with the open end upwards (no more than 7 straws per support).

**N.B.:** make sure to position the label about 35 mm from the splayed end of the straw.

4. Place 1 to 3 straws in the unit to be pre-cooled:
  - a) Place a vitrification block in the upper position.
  - b) Remove a straw from the straw support and insert it into the vitrification block.
  - c) Replace the vitrification block **immediately** in the lower position

Repeat the operation with the other vitrification supports.


**N.B.,** Liquid air will form in the straw if it is kept in the upper position for any length of time.

**N.B.:** The CBS™ High Security Vitrification straw must be kept in the lower position for at least 30 seconds in order to cool the inside completely and effectively.



5. After loading the biological sample in the gutter of the handling rod-capillary assembly (in accordance with the laboratory protocol):
  - a) Lift the vitrification block and block it in the upper position.
  - b) Position the end of the handling rod-capillary assembly gutter carefully opposite the open end of the straw.
  - c) Insert the handling rod-capillary assembly into the straw with a **quick, continuous movement**, pressing on it with a finger until the nearest end of the colored handling rod is brought flush with the open end of the straw.
  - d) Continue inserting by using the blue insertion device. When this blocks against the straw, the handling rod-capillary assembly is properly in place.

6. Seal the splayed open end of the straw and lower the vitrification block into the lower position.

7. Tilt the goblet support and introduce the CBS™ High Security Vitrification straw into the goblet using the tweezers, making sure it remains immersed in the liquid nitrogen.

Repeat operations 6 to 8 to condition all the other biological samples.

8. Using the movable goblet support, transfer all the vitrified straws to the dedicated storage tanks.

**⚠ Make sure the nitrogen level in the visotubes is sufficiently high to maintain the straws at cryogenic temperature while they are being transferred.**


## 5.2 Opening CBS™ High Security Vitrification straws

1. Before filling the unit with liquid nitrogen, check that:
  - the goblet support is in its housing
  - the opening device is in the raised position
  - the work area is properly prepared



 **Important: the opening device is raised to limit frost and cold on the blade. It must be raised after each straw has been cut.**

 2. Fill the unit with liquid nitrogen starting with the goblet and make sure it is **well filled**. If it is not to be used immediately, place the cover to reduce evaporation and thus maintain a sufficient level of nitrogen (cf. section 4).

3. Bring the goblet containing the CBS™ High Security Vitrification straw(s) to be opened from the storage unit to the INTEGRAL HSV Vitrification Unit using a goblet support, then insert them into the housing.

4. Incline the goblet support, grasp the straw to be opened with the tweezers, without removing it from the liquid nitrogen, and bring it to the cutting area.

5. Tilt the opening device into the horizontal position.


6. Grip the straw between two fingers and pull it until it touches the blade.
7. Turn the straw in the direction shown by the arrow.


8. Once the straw catches on the blade, give it a complete turn.


9. Free the straw from the opening device, leaving it immersed in the cutting area.

10. Lift the opening device into the vertical position.

11. Bring the straw into the holding area without removing it from the liquid nitrogen

12. Keeping the straw in the liquid nitrogen, remove the cut section.  
Take the handling rod-capillary/straw in order to bring it close to the thawing area, slide out of the straw and plunge the end of the gutter into the first drop of warming solution prepared beforehand in accordance with the laboratory protocol.

## 6 Maintenance

The INTEGRAL HSV requires no special maintenance.

The opening device is guaranteed for 1000 cutting operations. The photo opposite shows how it should be replaced when necessary.

### 6.1 *Emptying the unit*

Allow the liquid nitrogen to evaporate in a room that meets safety criteria concerning the use of liquid nitrogen and in compliance with local regulations.

 **To avoid any risk of splashes and burns, it is forbidden to empty the unit by turning it upside down.**

### 6.2 *Cleaning*

Cleaning must be done once the chamber is empty and at room temperature.

- 1 Remove all accessories from the unit.
- 2 Remove the platform and dry all components.
- 3 Remove the chamber and dry foam and chamber with a clean cloth.
- 4 Clean and decontaminate all accessories and the chamber of the INTEGRAL HSV unit with the usual products (distilled water, wipes, alcohol).

 **Never use solvents such as acetone, acidic or alkaline liquids, etc.**

## 7 Technical data

### 7.1 *General*

- Dimensions
  - a. Width: 283 mm - Length: 487 mm - Height: 265 mm with cover, 222 mm without cover
- Volume of chamber: 5L
- Weight: 5 kg

### 7.2 *Operating conditions*

- Only for used indoors at room temperature
- Storage: -20°C / +60°C
- Use: 18°C / 25°C with 40 to 70% humidity

### 7.3 *Working life*

3 years maximum.

## 8 Important information

### 8.1 Ordering

To order spare parts, refer to section 3.3: "Components and spare parts" above.

### 8.2 Warranty conditions for new equipment

Cryo Bio System guarantees its equipment against malfunctions resulting from defective material, manufacture or design for a period of one year from the date indicated on the warranty certificate.

All parts and labor are free of charge during this period.

Malfunctions will only be taken into account if they appear within 12 months of the date of effect of the warranty and provided that the equipment has been used in strict compliance with this user manual.

The warranty does into apply if the malfunction is the result of:

- work performed on the equipment without authorization,
- incorrect use that does not comply with the technical specifications and this user manual,
- force majeure.

In accordance with the above, Cryo Bio System promises to repair or replace all parts that are unusable for reasons covered by the warranty throughout the warranty period.

The warranty period will not be extended on account of any work carried out under the terms of the warranty, i.e. repairs or replacements of any or all of the equipment delivered.

Any dispute regarding the interpretation and/or application of the present warranty shall be settled exclusively in accordance with French law. The courts of the city of Alençon shall have sole jurisdiction, even in the case of warranty claims or multiple defendants.

### 8.3 Waiving of liability

Cryo Bio System is not liable for any damage due to external influences or to treatment and use of the equipment that is inappropriate or in contradiction with the instructions given in this user manual.

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