



VIP treatment for high success rate embryos in IVF

Miri®

Multi-room Incubator for IVF

Affordable Quality - Excellent Stability

Just a few of the things you'll love about the Miri® Multi-room Benchtop Incubator, packed with great features.



table of contents

- Miri® Benchtop Incubator Parts and General Features4
- Quality Control Features7
- Accessories8
- General Specifications9
- Ordering Information9



Miri®

“An advanced temperature regulation system for routine/long-term embryo incubation at your fingertips”

The Miri® has six (6) chambers which are completely independent of each other. This is ideal because any disruption (e.g. temperature drop after opening the lid) has zero impact on the rest of the system. Furthermore, calibration is so much simpler because there is no crossover of heat from adjacent chambers.

Temperature regulation is thus completely independent per chamber. The Miri® features a total of twelve (12) temperature controlled points. That is two (2) points for every chamber: one on the bottom and another on the heated lid. The heated lid is another great feature of the Miri® as it prevents condensation and enhances temperature uniformity across cultured dishes.



Hand in hand to provide only the best in embryo incubation

About Esco Medical

Esco Medical is one of the divisions of the Esco Group of companies, which targets innovative technological solutions for fertility clinics and laboratories. Esco Medical is the leading manufacturer and innovator of high quality IVF equipment that are designed in Denmark and mostly made in the EU.



FEATURES:

Heated Lid

- Prevents condensation
- Enhances temperature regulation and recovery
- Excellent uniformity between the top and the bottom
 - Accuracy and Uniformity: ± 0.2 °C

Heated bottom

- Provides direct heat transfer to the cultures through the optimization plate for stable heat regulation.
- Removable Heating Optimization Plate with wide selection of inserts.
 - Accuracy and Uniformity: ± 0.2 °C

Optional SAFE Sens[®] Integration

For continuous pH monitoring. *See page 6 for more info.*

Six (6) Chambers

- Completely individual chambers for easier calibration, faster recovery, less disruption, and prevents cross-contamination.

Control Panel Buttons and LED Display

Has large LED display that can be easily seen from a distance. The simple 4-button control panel allows for easy and intuitive operation.

Mute Button

Temporarily mutes alarm messages and sound for five (5) minutes



Oxygen range: 5-20%
 Carbon dioxide range: 1.9-10%
 Temperature: 25 Ambient to 40°C
 Gas recovery: less than 3 minutes
 Temperature recovery: less than 1 minute

Common Stressors:

- Temperature fluctuations
- Gas concentration fluctuations
- Non-optimal pH
- Volatile Organic Compounds (VOCs)

i Maximize embryo growth potential by providing **VIP treatment**


Elevated O₂ concentration isn't always a good thing

While oxygen (O₂) is necessary for normal aerobic metabolism, it is a double-edged sword as it can harm the developing embryo through oxidative damage. Recent studies highlight the benefit of having suppressed oxygen levels when incubating human embryos reflecting the natural low oxygen conditions in the womb.

Shhh... Do not disturb

The Miri® has an overall design that provides cultured embryos a minimum-stress environment. The 6-chamber format prevents cross-contamination while HEPA+VOC filtration cleans the airstream. The small chamber volumes and direct heat regulation further translate to faster temperature and gas recovery.

Fast Recovery



There are many advantages to using benchtop multi-room incubators. One important benefit is the speed of recovering temperature and gas parameters after opening a chamber.


The little details count



IVF practitioners deal with precious, fragile and sensitive embryos, and often, the little details make a big difference. The Miri® has a large LED display that can be easily seen from a distance. Also, the glass lid tops, while acting as chamber insulators, can be written on — a very useful feature for organization.

Provide total control of the gas phase environment

Airflow Diagram



The built-in gas mixer and the high-performance CO₂ and O₂ sensors allow accurate control of gas phase composition in the chambers.

The Miri[®] is built with Excellent Quality Control Features



The Miri[®] has reliable gas mixing system that allows gas phase flexibility.

The gas mixer of Miri[®] gives total control over CO₂ and O₂ concentration levels while also giving flexibility over what gas input is desired.** Moreover, the HEPA+VOC filter and UV sterilization* ensure only the highest quality of air is circulated to the cultures.



Pre-mixed gas is NOT required

Miri[®] can accept pure gas (and even pre-mixed gas if desired).

The use of 100% CO₂ and 100% N₂ allows the incubator to achieve steady-state gas condition faster than premixed gases. This is important as this will stabilize the pH during embryo development.

High quality airstream via HEPA+VOC filter + UV

The filter module can be easily replaced once used.

The gas in the Miri[®] is continuously recirculated through a HEPA/VOC filter. A UV-C light (185 nm) sterilizes the airstream before passing through the filter.

A suite of IVF-essential features

BNC connection for pH monitoring**, USB communication port, and port for external alarm monitoring.

The Miri[®] can be connected to a PC to avail of data logging via the supplied software included. Connections to external alarm monitoring systems and pH measurements** are also possible.

Stress-free validation of chamber parameters



PT1000 temperature sensors are built-in, which are completely independent from the main circuitry. Gas sampling ports are likewise available for all 6 chambers.




The Miri[®] can be connected to an external device such as the Esco Miri[®] GA for gas and temperature validation.

Full-featured and user-friendly


Control panel, display, and data logging software



Complete parameters are displayed. Histories of any alarm events are logged.



The data-logger stores continuous performance data of the machine throughout its use. These can be viewed in graphs.



* Not available in US
 **pH measurement data logging, pH measurement through pH probe is not available in the US. See page 6 for optional Safe Sens integration for continuous pH measurement.
 *** Input of pure gases is recommended




Conditions that put the Miri® into alarm state are recorded. It is possible for the software to send email alerts as well.

The Miri® can be connected to an easy-to-use, feature-packed data logging software installed on any ordinary PC and connected via USB. Multiple machines can be connected and managed from a single computer. All real-time parameters of the machine can be conveniently viewed. These include the temperature of all monitored temperature and gas concentration points, gas input pressures, gas flow rates, current gas readings, and all set points.

All performance data of the machine including alarms are continuously logged and can be viewed in graphs. The datalogger also automatically generate reports weekly which makes it more convenient for the user.

Accessories



For Miri® with integrated SAFE Sens (optional), order a different plate with hole to accommodate the SAFE Sens sensor (see ordering information on page 7).

The dishes fit into the inserts so that the heat is directly transferred to the media.

Heating optimization plates

Each chamber contains a heating optimization plate to facilitate heat transfer directly to the culture dishes.

- Has inserts to fit various dish sizes
- Removable for easy cleaning



Nunc™



Falcon™



Vitrolife™



LifeGlobal / GPS Dishes™



Nippon™



Sparmed - Oosafe™

Total Capacity


Heating plates customized for several types of dishes:

- 4 × Falcon® Ø 50/60 mm
- 8 × Falcon® Ø 35 mm
- 4 × Nunc® Ø 54/60 mm
- 8 × Nunc® Ø 35 mm
- 4 × Vitrolife® Dishes
- 4 × LifeGlobal / GPS Dishes™
- 4 × Nippon™
- 4 × Sparmed - Oosafe™

SAFE Sens* Continuous pH Monitoring

The Miri® can be installed with an integrated SAFE Sens technology for fast, effective, and non-invasive continuous pH monitoring product for *in vitro* fertilization (IVF) procedures.

The SAFE Sens technology employs an optical fluorescent measurement technology, used in combination with disposable sensors, which accurately and reliably monitors the pH of small volumes of fluids such as the media used in IVF.



Key Features



Continuous pH measurement

- Reading and recording every 30 minutes (default setting - adjustable).
- Single use sensor probe for up to seven (7) days of pH readings.



Data-Logging System**

- Data Logging and user alarms.
- Each TrakStation® can be connected to multiple incubators.



Compact and Efficient

- No more unnecessary openings of your incubator for spot pH measurement.
- Only requires 100 µL of media + 150 µL of oil.

* SAFE Sens is a trademark brand of Blood Cell Storage, Inc. (BCSI). SAFE Sens integration is currently offered as a factory-installed option.

**Minimum system requirements for datalogger PC/Tablet:

• Intel Core 2 Duo or AMD Athlon X2 at 2.4 GHz processor • 4Gb RAM • 15Gb Hard Disk space • Integrated Video Card • Monitor with resolution 1024 x 768 • Windows 7 pro/8 Pro/10 OS with 64 Bit architecture • USB 3.0 port for each connected device



Miri® Tablet PC and Holder

Let go of bulky computers and save space with a tablet. Monitor the performance of your Miri® incubator with a full-featured and user-friendly control panel, display and data logging software.

*Tablet holder can fit tablets with maximum of 26mm thickness and 201mm height.



MRA-1014 - Stacking frame for 2 units



Miri® Stacking Frame with drawer

Miri® Stacking Frame

Miri® has a stacking system to maximize space in your IVF laboratory.

General Specifications



Miri® Multi-room incubator : PRODUCT CODE: MRI-6A10-_-_-

Overall Dimensions (W x D x H)	700 x 580 x 150 mm (27.6" x 22.9" x 6")
Power Supply	115 / 230V, 50/60 Hz
Power Consumption	280 W
Temperature Control Range	25 - 40°C
*Gas Consumption (CO ₂)	< 2 L/h
**Gas Consumption (N ₂)	< 12 L/h
CO ₂ Control Range	1.9 - 10%
O ₂ Control Range	5 - 20%
Input Gas Pressure (CO ₂)	0.6 bar (8.7 psi)
Input Gas Pressure (N ₂)	0.6 bar (8.7 psi)
Net Weight	35 kg (77.2 lbs)
Shipping Weight	40 kg (88.2 lbs)
Shipping Dimension	840 x 735 x 300 mm (33.1" X 29" x 11.9")

* Under normal condition (CO₂ set point reached at 5.0%, all lids closed)

** Under normal condition (O₂ set point reached at 5.0%, all lids closed)




Ordering Information

ITEM CODE	MODEL CODE	DESCRIPTION
Unit		
2070047	MRI-6A10-8	Miri® Incubator; 230V, 50/60Hz
2070048	MRI-6A10-9	Miri® Incubator; 115V, 50/60Hz
2070086	MRI-6A10-SS-8	Miri® Incubator; with SAFE Sens for pH measurement, 230V, 50/60Hz
2070087	MRI-6A10-SS-9	Miri® Incubator; with SAFE Sens for pH measurement, 115V, 50/60Hz
1320045	MRI-GA	Miri® GA CO ₂ / O ₂ & Temperature Validation Unit, 115V / 230V
Accessories		
1320140	TBA	PC Holder for Miri®
1320141	TBA	Tablet PC for Data Logging of MRI-6A10
1320142	TBA	Set of Tablet PC and Holder For Data Logging of MRI-6A10
1320191	TBA	SAFE Sens TrakStation, a tablet with SAFE Sens Software, for pH monitoring
1320011	MRA-1007	HEPA+VOC filter (recommended to be replaced every 3 months)
1320018	MRA-1014	Stacking frame for 2 units
1320226	TBA	Stacking Frame for 2 Units, With Drawer at the Bottom
1081277	TBA	SAFE Sens SV2 Sensor; Pack of 10 pieces (shelf-life 12 months)
1081278	TBA	SAFE Sens QC2 Alignment Tool

ITEM CODE	MODEL CODE	DESCRIPTION
Accessories		
132003	TBA	Insert for Falcon® Dishes
132004	TBA	Insert for NUNC® Dishes
1320070	TBA	Insert for Vitrolife® Dishes
1320099	TBA	Insert for Nipro Dishes
1320100	TBA	Insert for LifeGlobal/ GPS Dishes™
1320101	TBA	Insert Without Footprint for Plain Dishes
1320118	TBA	Insert for Sparmed - Oosafe™
1320219	TBA	Insert for Falcon® Dishes, with hole for SAFE Sens
1320220	TBA	Insert for NUNC® Dishes, with hole for SAFE Sens
1320221	TBA	Insert for Vitrolife® Dishes, with hole for SAFE Sens
1320222	TBA	Insert for Nipro Dishes, with hole for SAFE Sens
1320223	TBA	Insert for LifeGlobal/ GPS Dishes™, with hole for SAFE Sens
1320224	TBA	Insert Without Footprint for Plain Dishes, with hole for SAFE Sens
1320225	TBA	Insert for Sparmed - Oosafe™, with hole for SAFE Sens

ESCO GLOBAL NETWORK




- Esco Medical Products:**
- Anti-Vibration Table (AVT)
 - CO₂ Incubators
 - Culture and Fertilization Dishes
 - Fertilisafe® IVF Workstation
 - Miri® Benchtop Multi-room Incubators
 - Mini Miri® Humidified Benchtop Incubator
 - Miri® TL Time Lapse Embryo Incubators
 - Miri® GA (Gas and Temperature Validation Unit)
 - Miri® GA Mini (Gas Analyzer)
 - Semi Closed Environment (IVF/ICSI)

Infertility is viewed as a problem that has social, psychological, and economic impacts to the afflicted individuals and couples. It is a global concern that knows no race nor creed. It has been estimated that 1 in 6 couples would struggle with infertility at least once in their lifetime.

The vision of Esco Medical is to support Assisted Reproductive Technologies (ART), such as IVF, by developing practical and state-of-the-art technological solutions for improving clinical success rates and patient satisfaction. All Esco Medical products are designed with the IVF clinic in mind and developed with the Silent Embryo Hypothesis as a guiding principle. The Silent Embryo Hypothesis states that the less disturbed an embryo can remain, the better its developmental potential will be.

It is on these foundations that Esco Medical remains committed to providing world class ART, worldwide. At Esco Medical, life has begun.



21 Changi South Street 1 • Singapore 486 777
 Tel +65 6542 0833 • Fax +65 6542 6920
 csis-medical@escoglobal.com • www.medical.escoglobal.com

Esco Global Offices: Bangladesh | China | Denmark | Hong Kong | India | Indonesia | Italy | Japan | Lithuania | Malaysia | Philippines | Russia | Singapore | South Africa | South Korea | Taiwan | Thailand | UAE | UK | USA | Vietnam



9010137 ART Equipment, MRI, Brochure, Ad, V, 040218
 Esco can accept no responsibility for possible errors in catalogues, brochures and other printed materials. Esco reserves the right to change specifications and prices without notice. All trademarks and copyrights in this material are the property of Esco and the respective companies.

PT Esco Britain Indonesia Cert. No. 0210185
 Esco Micro Products Cert. No. 0210176
 PT Esco Britain Indonesia Cert. No. 0210185/01001