

Tested Individual: Supine **Shoulder Position:** Flexed 90 degrees **Dynamometer Position:**


Just proximal to elbow

HOGGAN SCIENTIFIC, LLC.

MUSCLE TESTING POSITIONS

Diagrams indicate muscle tests with transducer placement, proper positioning and stabilization for test





Tested Individual: Supine **Shoulder Position:** Flexed 90 degrees

Dynamometer Position: Just proximal to elbow



SHOULDER ABDUCTION

Tested Individual: Supine **Shoulder Position:** Abducted 90 degrees **Dynamometer Position:** Just proximal to elbow



SHOULDER EXTERNAL ROTATION

Tested Individual: Supine **Shoulder Position:** Abducted 45 degrees **Elbow Position:** Flexed 90 degrees **Dynamometer Position:**

Just proximal to wrist



SHOULDER INTERNAL ROTATION

Tested Individual: Supine **Shoulder Position:** Abducted 45 degrees **Elbow Position:**

Flexed 90 degrees **Dynamometer Position:** Just proximal to wrist



ELBOW FLEXION

Tested Individual: Supine **Elbow Position:** 90 degrees

Dynamometer Position: Just proximal to wrist



ELBOW EXTENSION

Tested Individual: Supine **Elbow Position:** 90 degrees

Dynamometer Position: Just proximal to wrist



FOREARM PRONATION

Tested Individual: Supine **Elbow Position:** 90 degrees **Dynamometer Position:**

Against dowel 20cm from dowel held in hand



FOREARM SUPINATION

Tested Individual: Supine **Elbow Position:** 90 degrees

Dynamometer Position: Against dowel 20cm from dowel held in hand



WRIST EXTENSION

Tested Individual: Supine **Elbow Position:** 90 degrees

Dvnamometer Position: Just proximal to metacarpal phalangeal joints of hand



HIP FLEXION

Tested Individual: Supine **Hip Position:** Flexed 90 degrees **Dynamometer Position:** Just proximal to femoral

condyles



HIP ABDUCTION

Tested Individual: Supine **Hip Position:** Extended & abducted 0 degrees

Dvnamometer Position: Just proximal to lateral knee joint line



HIP EXTENSION

Tested Individual: Supine Lower Limb Position: Knee extended with distal limb supported on block

Dynamometer Position: Just distal to malleoli on Achille's tendon.



KNEE FLEXION

Tested Individual: Sitting Lower Limb Position: Hip & knee flexed 90 degrees

Dynamometer Position: Just distal to malleoli on a Achilles tendon.



KNEE EXTENSION

Tested Individual: Sitting Lower Limb Position: Hip & knee flexed 90 degrees

Dynamometer Position: Just proximal to malleoli



ANKLE DORSIFLEXION

Tested Individual: Supine **Lower Limb Position:** Knee extended & ankle in neutral dorsiflexion

Dynamometer Position: Just proximal to metacarpal phalangeal joints



ANKLE PLANTARFLEXION

Tested Individual: Supine **Lower Limb Position:** Knee extended & ankle in neutral dorsiflexion

Dynamometer Position: Over metacarpal phalangeal joints