# Instructions for use

EN

## Dental Implant, Instruments, URIS OMNI System\_T LOC Component

Instructions for use / Valid only in EU

IFU\_TTI0013, Revision 00, Document valid as of 30-Apr-2020

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Phone: +1 714 956 1488

Fax: +1 714 956 1491

[www.truabutment.com](http://www.truabutment.com)

[www.urisimplants.com](http://www.urisimplants.com)

---

EC REP

### AeMi World e.K.

Bugenhagenstr. 8, 10551 Berlin, Germany

Phone: +49 (030) 8620 3461

Fax: +49 (030) 8620 3789

---

**PARTS IDENTIFICATION**

T LOC Component Set  
(Titanium Cap)

**No. TLC-CST(6)**

Titanium Cap  
Black LDPE Male,  
Nylon Replace Males (blue,pink,clear)  
and Block out Spacer

**T LOC Titanium Cap**

**No. TLC-TC5423**

Ti-6Al-4V ELI

**T LOC Block Out Spacer White**

**No. TLC-OS56W**

Silicon

**T LOC Black Processing  
Replacement Male**

**No. TLC-PRM56K**

LDPE(Black)

**1.5 lb T LOC Retention  
Replacement Male**

**No. TLC-RRM47B**

Nylon(Blue)

**3.0 lb T LOC Retention  
Replacement Male**

**No. TLC-RRM47P**

Nylon(Pink)

**5.0 lb T LOC Retention  
Replacement Male**

**No. TLC-RRM47T**

Nylon(Clear)

**INDICATIONS**

T LOC Components are laboratory pieces used for fabrication of dental prosthetics.

**CONTRAINDICATIONS**

Not appropriate where a totally rigid connection is required. Use on a single implant with divergence of greater than 20 degrees is not recommended.

**CAUTION**

Federal (U.S.A.) law restricts this device to sale by or on the order of a licensed dentist.

**MRI SAFETY INFORMATION**

The T LOC Component has not been evaluated for safety and compatibility in the MR environment. It has not been tested for heating, migration, or image artifact in the MR environment. The safety of the T LOC Component in the MR environment is unknown. Scanning a patient who has this device may result in patient injury.

**SINGLE-USE DEVICES**

T LOC Components: The inadvertent re-use of T LOC nylon males could cause loss of retention for the overdenture due to wear from previous use or damage during removal.

**STERILIZATION**

All components and instruments are supplied NON-STERILE.

**PROSTHETIC PROCEDURE****RETENTION MALE PLACEMENT BY THE DENTIST**

1. Place a Block Out Spacer White (No. TLC-OS56W - also contained in No. TLC-ST(6) Package) over the head of each TruAbutment DS T LOC. The spacer is used to block out the area immediately surrounding the abutment. The space created will allow the full resilient function of the pivoting T LOC titanium cap over the T LOC Black Processing Replacement Male.

**NOTE:** If the Block Out Spacer White does not completely fill the space between the tissue and the metal denture cap, it is necessary to block out any remaining undercut to prevent the added acrylic resin from locking the denture onto the abutment. This can be accomplished by stacking more Block Out Spacers.

2. Insert a Ti Cap with Black Processing Replacement Male (No. TLC-PRM56K - also contained in No. TLC-CST(6) Package) onto each T LOC, leaving the Block Out Spacer White beneath it. The Black Processing Replacement Male will maintain the overdenture in the upper limit of its vertical resiliency during the processing procedure.
3. Prepare a recess in the denture to accommodate the protruding T LOC Male. There must be no contact between the denture and the T LOC cap. If the denture rests on the T LOC Cap, excess pressure on the implant will result.
4. Use a Chairside Lightcure Acrylic Resin Syringe Kit to light cure bond the T LOC Denture Cap Processing Male into the denture, or mix a permanent self-curing acrylic and place a small amount in the recess of the denture and around the T LOC Cap of the T LOC Denture Cap Processing Male.
5. Insert the denture into position in the oral cavity. Guide the patient into occlusion, maintaining a proper relationship with the opposing arch. Maintain the denture in a passive condition, without compression of the soft tissue, while the acrylic sets. Excessive occlusal pressure during the setting time may cause tissue recoil against the denture base and could contribute to dislodging and wear of the nylon males.
6. After the acrylic resin has cured, remove the denture and discard the Block Out Spacer White. Use a bur to remove excess acrylic, and polish the denture base before changing to the final male.
7. Remove the Black Processing Replacement Male from the T LOC cap.
8. The Replacement Male must seat securely into place, level with the rim of the cap.

**T LOC CAP MALE PLACEMENT BY THE LABORATORY**

1. Insertion of the proper Truabutment DS T LOC at tissue level must be completed before beginning the following impression procedure.
2. Pour the master cast.
3. Before waxing and processing the appliance, place a T LOC Black Processing Cap

Male onto each abutment in the master cast. Make sure the Cap Processing Male Assembly is fully seated.

4. Set the teeth and wax the appliance. Proceed with the processing technique of your choice through the boil-out step.
5. After the boil-out, remove the T LOC Cap and Black Processing Cap Male. Place a Block Out Spacer White over the head of each Truabutment DS T LOC. The spacer is used to block out the immediate area surrounding the Truabutment DS T LOC. The space created will allow the full resilient function of the pivoting metal denture cap over the T LOC Nylon Male.
6. Re-insert the T LOC Cap and Black Processing Cap Male onto each Truabutment DS T LOC, leaving the Block Out Spacer White beneath it. The Black Processing Replacement Male will maintain the overdenture in the upper limit of its vertical resiliency during the processing procedure.
7. Complete the processing and discard the Block Out Spacer White. Avoid damage to the final male by polishing the denture base before changing to the final T LOC Nylon Replacement Male.

**HOW TO CHANGE THE T LOC MALE**

1. Remove the nylon male from the T LOC cap.
2. Firmly push a T LOC Replacement Male into the empty T LOC cap. The Replacement Male must seat securely into place, level with the rim of the cap. Use of multiple T LOCs (3 or more) in the same dental arch may require use of a combination of various retention Male Replacements for easier removal of the prosthesis by the patient.

**PATIENT CARE**

Good oral hygiene is vital to attachment success. The TruAbutment DS T LOC must be thoroughly cleaned each day to prevent wear of the abutments due to buildup of abrasive plaque in the socket of the abutment. The use of a soft nylon bristle or end-tufted toothbrush, and superfloss to polish the abutments should be taught. A non-abrasive gel toothpaste, and an irrigation system is recommended to keep the socket of the T LOC clean.

Patients should maintain a three to four month recall for cleaning evaluation of abutment and male attachment. Replace male attachment if worn. The inside socket of the T LOC and the sulcus area around the implant abutment are the primary areas of concern. Use plastic instruments for scaling the abutments. Do not use metal instruments which may create scratches on the abutment surface. Examine patients for signs of inflammation around the implant abutments, and for implant mobility. Use a 30N-cm torque wrench to make sure the Truabutment DS T LOC is tight before dismissal.

**CAUTION**

**THE PATIENT MUST NOT BITE** their Overdentures into place as this force will result in improper wear of the Abutments, including the Replacement Males in the Overdenture. Remove the Overdenture by placing the thumbs under the edges of the Overdenture flanges and pulling each side upward/downward simultaneously. Use of the tongue may aid in removal. Once removed, a thorough cleaning is recommended.

**Cleaning your Implant Retained Denture**




Instruct the patient to follow the protocol below to ensure the longevity of their Overdenture. Fill a washing basin with some warm water to prevent fracture of the Overdenture. Apply non-abrasive toothpaste onto the soft bristle toothbrush and thoroughly clean every surface of the Overdenture. Each night, remove Overdenture and immerse in a cup of plain water.

**Further Information**

Traditional restorative protocols should be followed to process the attachments into the patient's Overdenture. Standard Overdenture care and maintenance should be followed in order to ensure the longevity of the restoration.

**LABELING SYMBOLS**

Symbols may be used on some international package labeling for easy identification.

	Do not reuse
	Batch code
	Date of manufacture
	Catalogue number
	Authorized representative
	Notified body number
	Caution, consult accompanying documents
	Manufacturer
	Consult instructions for use
	Do not use if package is damaged
	Non-Sterile

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**TruAbutment Inc.**