

Performance has a new style

Headgear

with elastic crown strap and integrated tube holder.
Self-adjusting for a more natural fit.

Small	1109307
Standard	1094082
Large	1094078

Reversible fabric frame

in soft suede or silky sateen.
Wisp silicone frame available for a more natural feeling

Silicone frame	1094080
Fabric frame	1094081

**‘Tip-of-the-nose’
cushion design,
featuring Auto Seal
technology**

Delivers a minimal contact experience with an exceptional seal

Petite	1094086
Small/ Medium	1094087
Large	1094088
Extra large	1112031



Silicone frame

Magnetic headgear clips* with adjustable tabs

Provides for quick and easy fitting and adjustment

Magnetic clip (x2)	1120434
Standard clip (x2)	1094083

360°rotating elbow

Features a quick release design with quiet exhalation

* intended for patients weighing 66 lb/30 kg and up

Helpful tips

For more fitting tips on how to adjust for leaks or cleaning instructions, please visit:

www.sleepapnea.com/products/masks/



- Before using the mask, hand wash the entire mask (see Mask Care in the Instructions for use). Wash your patient's face. Do not use moisturizer or lotion on your patient's face.
- Four cushion sizes are available.
For Wisp adult: Petite, Small/Medium, Large and Extra Large.
For Wisp youth: Petite, Small/Medium and Large.
- The most common mistake when fitting is over-tightening the headgear. Signs of overtightening include redness, sores, or bulging skin around the edges of the mask.

Fitting instructions

Achieving the right fit with Wisp adult and youth is easy. Simply follow the instructions below. Firstly, select the appropriate cushion size. There are four cushion sizes available:



1 Before putting on the mask, disconnect the headgear clips away from the frame by twisting the clips (Wisp youth) or disconnecting the magnetic clips (Wisp adult).



2 Gently hold the cushion in place over the nose and pull the headgear over your head. The headgear should lie flat against the back of your head.



3 Holding the mask in place, reconnect the clips by pushing them back onto the mask frame (Wisp youth) or reconnecting the magnetic clips (Wisp adult).



4 Pull back on the headgear tabs to evenly adjust the upper and lower headgear straps. The mask should rest comfortably against your face.



5 You may optionally slide the mask tubing through the headgear crown loop before attaching the CPAP tubing.



6 Assume different sleeping positions and move around until comfortable. If there are any excessive air leaks, make final adjustments while lying down.

© 2017 Koninklijke Philips N.V. All rights reserved. Specifications are subject to change without notice.

www.philips.com/respironics
SJ 9/20/17 PN499443

