



10310001/10310002
Hephaestus V094/V213
Vacuum Ovens
User Manual

Table of Contents

| | |
|---|-----------|
| I. Disclaimer & Warranty | 3 |
| LIMITED WARRANTY & DISCLAIMER | 3 |
| LIMITATION OF LIABILITY..... | 6 |
| II. Safety | 7 |
| III. Hephaestus V094/V213 Vacuum Ovens – Introduction..... | 9 |
| 3.1 Vacuum Oven Overview | 10 |
| IV. Technical Specifications..... | 12 |
| V. Standards and Conformity..... | 13 |
| VI. Inspection | 14 |
| 6.1 Packing List | 14 |
| VII. Operating Instructions | 15 |
| 7.1 Temperature Controller Panel | 15 |
| 7.2 Functions | 16 |
| 7.3 Priming the Vacuum Pump | 17 |
| 7.4 Priming the Vacuum Oven | 17 |
| 7.5 Setting of Temperature | 17 |
| 7.6 Start / Stop Operation | 17 |
| 7.7 Setting the Deviation Alarm and Timer..... | 18 |
| 7.8 Setting the Pre-set Start Function | 19 |
| 7.9 Setting the ENGINEER SETUP | 19 |
| 7.10 Unit, SYSTM Timer and PASSWORD Settings..... | 20 |
| VIII. Cautions..... | 21 |
| X. Fault Diagnosis..... | 22 |
| XI. Maintenance and Cleaning | 24 |

I. Disclaimer & Warranty

- The intellectual property rights for this manual belongs to Holland Green Science Corporation (hereinafter referred to as “the Company”).
- Products of the company are under the patent protection of China and other countries and regions (inclusive of obtained patents and patents currently under application).
- The company reserves the right to change product specifications without notice.
- The contents of the manual are subject to change without notice.
- The manual’s content must not be copied, extracted, or modified in any form without the Company’s permission.

LIMITED WARRANTY & DISCLAIMER

Products sold by the Company are intended solely for the use(s) specified by the manufacturer and are not to be used for any other purposes including, but not limited to, unauthorized commercial purposes. Buyer warrants and represents that Buyer will properly test, use, qualify, and/or validate each Product for intended use, and, to the extent authorized, manufacture and market any final articles made from the Company’s Products in accordance with the practices of a reasonable person who is an expert in the field, and in strict compliance with all applicable national, state and local laws and regulations. Buyer acknowledges and understands that the purchase of Products only conveys to Buyer the non-transferable right for only the Buyer to use the Products purchased in compliance with the applicable manufacturer’s Published User Statement, Limited User Statement, or Limited License, if any. While this manual contains certain representations of the potential dangers from the instrument or equipment, Buyer is solely responsible for conducting any research necessary to learn the hazards involved for any of its intended uses of the instrument or equipment purchased from the Company and to properly warn its customers, employees and other personnel who may be exposed to such instrument or equipment of any risks involved in using or handling the instrument or equipment.

The Company warrants that this product will operate or perform substantially in conformance with the Company’s published specifications and be free from defects in material and workmanship, when subjected to normal, proper and intended storage and usage by properly trained personnel, for the duration of the manufacturer’s warranty period set forth in this product documentation, published specifications or package inserts. If a manufacturer’s warranty period is not specified in the product documentation, published specifications or package inserts, the warranty period shall be one (1) year from the date of Bill of Sale to Buyer for equipment, and the earlier of any stated expiration dates or ninety (90) days for all other products (the "Warranty Period"). The Company agrees during the Warranty Period, to repair or replace, at the Company’s option, defective Products so as to cause the same to operate in substantial

conformance with said published specifications; provided that Buyer shall (a) promptly notify the Company in writing upon the discovery of any defect, which notice shall include the product model and serial number (if applicable) and details of the warranty claim; and (b) after the Company's review, the Company will provide Buyer with service data and/or a Return Material Authorization ("RMA"), which may include biohazard decontamination procedures and other product-specific handling instructions, then, if applicable, Buyer may return the defective Products to the Company with all costs prepaid by Buyer. Replacement parts may be new or refurbished, at the election of the Company. All replaced parts shall become the property of the Company. Shipment to Buyer of repaired or replacement Products shall be made in accordance with the Delivery provisions of the Terms & Conditions agreed to by the purchase of this product. Consumables are expressly excluded from this warranty.

Notwithstanding the foregoing, Products supplied by the Company but not manufactured by the Company are not warranted by the Company, and the Company disclaims all warranties associated with the accuracy of any original manufacturer or third-party supplier Product information, but the Company agrees to assign to Buyer any warranty rights in such Product that the Company may have from the original manufacturer or third-party supplier, to the extent such assignment is allowed by such original manufacturer or third-party supplier.

In no event shall the Company have any obligation to make repairs, replacements or corrections required, in whole or in part, as the result of (i) normal wear and tear, (ii) accident, disaster or event of force majeure, (iii) misuse, fault or negligence of or by Buyer, (iv) use of the Product in a manner for which they were not designed, (v) causes external to the Product such as, but not limited to, power failure or electrical power surges, (vi) improper storage and handling of the Product or (vii) use of the Product in combination with equipment or software not supplied by the Company. If the Company determines that Product for which Buyer has requested warranty services are not covered by the warranty hereunder, Buyer shall pay or reimburse the Company for all costs of investigating and responding to such request at the Company's then prevailing time and materials rates. If the Company provides repair services or replacement parts that are not covered by this warranty, Buyer shall pay the Company therefor at the Company's then prevailing time and materials rates.

ANY INSTALLATION, MAINTENANCE, REPAIR, SERVICE, RELOCATION OR ALTERATION TO OR OF, OR OTHER TAMPERING WITH, THE PRODUCT PERFORMED BY ANY PERSON OR ENTITY OTHER THAN THE COMPANY WITHOUT THE COMPANY'S PRIOR WRITTEN APPROVAL, OR ANY USE OF REPLACEMENT PARTS NOT SUPPLIED BY THE COMPANY, SHALL IMMEDIATELY VOID AND CANCEL ALL WARRANTIES WITH RESPECT TO THE PRODUCT.

THE OBLIGATIONS CREATED BY THIS WARRANTY STATEMENT TO REPAIR OR REPLACE A DEFECTIVE PRODUCT SHALL BE THE SOLE REMEDY OF BUYER IN THE EVENT OF A DEFECTIVE PRODUCT. EXCEPT AS EXPRESSLY PROVIDED IN THIS WARRANTY STATEMENT, THE COMPANY DISCLAIMS ALL OTHER WARRANTIES, WHETHER EXPRESS OR IMPLIED, ORAL OR WRITTEN, WITH RESPECT TO THE PRODUCTS, INCLUDING WITHOUT LIMITATION ALL IMPLIED WARRANTIES OF MERCHANTABILITY OR FITNESS FOR ANY PARTICULAR PURPOSE OR USE. THE COMPANY DOES NOT WARRANT THAT THE PRODUCTS ARE ERROR-FREE OR WILL ACCOMPLISH ANY PARTICULAR RESULT.

NOTHING IN THIS MANUAL OR IN THE TERMS AND AGREEMENT AGREED TO AT THE TIME OF PURCHASE SHALL BE CONSTRUED TO PROVIDE FOR THE TRANSFER, ASSIGNMENT OR ALIENATION OF THE LIMITED WARRANTIES PROVIDED BY THE COMPANY HEREIN TO ANY OTHER PARTY OTHER THAN BUYER, INCLUDING BUT NOT LIMITED TO, SUCCESSORS, HEIRS, SUBSIDIARIES OR ANY OTHER THIRD PERSON.

THE WARRANTIES IDENTIFIED ABOVE ARE THE SOLE AND EXCLUSIVE WARRANTIES WITH RESPECT TO THIS PRODUCT AND ARE IN LIEU OF ALL OTHER WARRANTIES, EXPRESS OR IMPLIED. ALL OTHER WARRANTIES ARE EXPRESSLY DISCLAIMED, INCLUDING, WITHOUT LIMITATION, THOSE OF NON-INFRINGEMENT, MERCHANTABILITY, OR FITNESS FOR A SPECIFIC PURPOSE OR USE (WHETHER ARISING FROM STATUTE, OR OTHERWISE IN LAW, OR FROM A COURSE OF DEALING, OR USAGE OF TRADE). THIS WARRANTY IS **NOT** TRANSFERABLE FROM THE ORIGINAL PURCHASER TO A SUBSEQUENT OWNER. FURTHER, THE COMPANY IS NOT LIABLE IN CASES OF DELIBERATE, NEGLIGENT OR ACCIDENTAL MISUSE OF THE PRODUCTS, USE WITH INAPPROPRIATE REAGENTS OR CONSUMABLES, DAMAGE CAUSED BY DISASTER, REPAIR OR MODIFICATIONS DONE BY ANYONE OTHER THAN THE COMPANY. THE COMPANY DOES NOT WARRANT THAT THE PRODUCT WILL NOT INFRINGE THE INTELLECTUAL PROPERTY RIGHTS OF A THIRD PARTY EITHER ALONE OR IN COMBINATION WITH OTHER PRODUCTS OR IN THE OPERATION OF ANY PROCESS. THE COMPANY'S TOTAL LIABILITY FOR BREACH OF THESE WARRANTIES SHALL BE LIMITED TO THE GROSS PURCHASE PRICE FOR THE PRODUCT AS STATED IN THE RELEVANT INVOICE OR ORDER UNDER WHICH THE ACTIVITY GIVING RISE TO LIABILITY ARISES. UNDER NO CIRCUMSTANCES SHALL THE COMPANY BE LIABLE FOR INDIRECT, INCIDENTAL, CONSEQUENTIAL, SPECIAL OR PUNITIVE DAMAGES (INCLUDING LOSS OF BUSINESS OR PROFITS) UNDER ANY THEORY OF LIABILITY HEREUNDER.

For claims under the warranty please contact your local supplier. You may also send the instrument directly to manufacturer, enclosing the invoice copy and by giving reasons for the claim.

LIMITATION OF LIABILITY

NOTWITHSTANDING ANYTHING TO THE CONTRARY CONTAINED HEREIN, THE LIABILITY OF THE COMPANY UNDER THESE TERMS AND CONDITIONS (WHETHER BY REASON OF BREACH OF CONTRACT, TORT, INDEMNIFICATION, OR OTHERWISE, BUT EXCLUDING LIABILITY OF THE COMPANY FOR BREACH OF WARRANTY (THE SOLE REMEDY FOR WHICH SHALL BE AS PROVIDED UNDER THIS SECTION) SHALL NOT EXCEED AN AMOUNT EQUAL TO THE LESSER OF (A) THE TOTAL PURCHASE PRICE THERETOFORE PAID BY BUYER TO THE COMPANY WITH RESPECT TO THE PRODUCT(S) GIVING RISE TO SUCH LIABILITY OR (B) ONE HUNDRED THOUSAND DOLLARS (\$100,000), WHICHEVER IS LESS.

With respect to any software products incorporated in or forming a part of the Product hereunder Buyer intends and agrees that such software products are being licensed and not sold, and that the words "purchase", "sell" or similar or derivative words are understood and agreed to mean "license", and that the word "Buyer" or similar or derivative words are understood and agreed to mean "licensee". Notwithstanding anything to the contrary contained herein, The COMPANY or its licensor, as the case may be, retains all rights and interest in software products provided hereunder. See the Terms & Conditions of purchase for details regarding such license.

II. Safety



IMPORTANT Safety Items to be Observed.

Safety Warnings

! DANGER (may cause serious damage to property and or casualties)

1. Please carefully read this user manual prior to operating the instrument and observe the instructions on safe operation.
2. The power source must be grounded reliably and away from any sources of electromagnetic interference.
3. Confirm that the voltage and frequency of the power supply matches the specified voltage and frequency specified prior to use.
4. The instrument should be operated on an independent power outlet.
5. Do not allow the instrument to be plugged in or unplugged whilst the power switch is in the on position.
6. Do not modify the power cord.
7. The working environment must be free of any flammable, explosive, volatile or corrosive substances.
8. The instrument is deenergised only when the power cable is disconnected from the power source.
9. This instrument is for indoor use only.
10. The instrument should be operated in a low humidity, well ventilated dust and water free environment, without direct sunshine or an intense magnetic, light or heat source.
11. Ensure that the instrument and its accessories are free of any potential defects.
12. Ensure that the instrument is secured tightly, and only standard accessories supplied by the manufacturers are used with this instrument.

! WARNING (may cause property damage or personal injury)

13. Prevent water from splashing on the electrical elements of this instrument.
14. Before assembly, disassembly, cleaning or maintenance, the power source must be disconnected.
15. It is prohibited to place flammable, explosive, volatile or corrosive substances into the oven.
16. Do not touch the chamber door, the chamber body or the surrounding surface when the set temperature is over 80°C (176°F) to avoid burns.

! ATTENTION (may affect operational performance or service life)

17. When pulling the power plug out from a power source, do not pull out using the power cord directly.
18. Only accredited and qualified professional repair technicians can open the instrument or conduct required repairs. Persons performing repairs on the instrument other than those selected or approved by the Company shall operate to void any warranty contained hereinabove for the product.
19. The 304 stainless steel liner is not acid resistant, please adhere to advised anti-corrosion measures and NEVER use acidic media within the oven to prevent damage.

III. Hephaestus V094/V213 Vacuum Ovens – Introduction

The Hephaestus V vacuum ovens are designed to efficiently process raw materials in a temperature-controlled vacuum environment. This process has the following advantages:

- The partial vacuum decreases the boiling point of water and thereby decreases the process temperature and increases the rate of moisture evaporation.
- The reduced processing temperature minimizes the oxidation rate of the product ensuring less oxidative damage.

The Hephaestus V Vacuum Oven is constructed from 304 stainless steel and the exterior surface of the oven is polymer coated. The inner chamber uses either a galvanized or a stainless-steel armour plate. Fiberglass insulation is used to thermally insulate the chamber. The cabinet door has a viewing window that allows visual inspection of the process. The window is constructed using a double-layered fortified glass design. The door uses a heavy-duty column door pin.

Temperature is monitored and controlled using a Pt100 sensor and a temperature controller. The temperature controller allows the output power with to be controlled using a PID. It also monitors the process temperature and corrects temperature fluctuations and monitors processing time(s).

3.1 Vacuum Oven Overview



- 1) vac-gauge
- 2) controller
- 3) LED switch
- 4) vac-valve
- 5) handle
- 6) LED



- 7) inlet
- 8) outlet/kf25
- 9) power switch
- 10) power cord
- 11) gas/air

The Hephaestus V094/V213 vacuum oven uses 5-15P power plug (see figure below).



3.3 After Sales Support

If problems are encountered or technical assistance support is required, when installing or using the instrument, please immediately contact the service department at ServiceUSA@hollandgreenscience.com. The Company may provide technical assistance and information regarding the instrument or equipment or service without charge at its sole discretion. Buyer assumes sole responsibility for any reliance on or use of such assistance and information, and the Company makes no warranty thereon.

Upon contact the following information is required to be provided:

- Product serial number (located on rear panel instrument nameplate)
- Warranty card
- Description of issue or problem (i.e., software or hardware)
- Method and or operating steps you have undertaken towards resolution.
- Your contact details inclusive of telephone number and email address.

For claims under the warranty please contact your local supplier. You may also send the instrument directly to manufacturer, enclosing the invoice copy and by giving reasons for the claim.

3.4 Proper Use

The instrument is designed for non-residential use e.g. in schools, laboratories or factories and to be used in conjunction only with accessories recommended within this manual and by the manufacturer.

IV. Technical Specifications

| Model Index Name | 10310001 | 10310002 |
|-------------------------------------|---|-------------|
| Voltage | AC 120V 60Hz | |
| Input Power (W) | 1350W | 1550W |
| Control Scope (°C) | RT+10~200 RT refers to Environment Temperature | |
| Fluctuating Degree (°C) | 1.0 | |
| Attainable Vacuum Degree | 133Pa | |
| Number of Shelves per Model (Layer) | 4 | 3 |
| Inner Chamber Material | 304 Stainless Steel | |
| Inner Chamber Size (mm) | 450x450x450 | 600x600x600 |
| Oven Outer Dimensions (mm) | 740x610x591 | 890x781x795 |

Note: 1. This product is produced according to Q/TIWY7-2004.
 2. All Technical specifications indexed are conducted under the situation with environment temperature 25°C, relative humidity \cong 85% and vacuum degree is \cong 0.1Mpa, and all data has been tested with a mercurial thermometer to the precision ± 0.1 °C. The mercurial head of the thermometer touching the surface of the shelf inside the cabinet.

V. Standards and Conformity

Construction in Accordance with the following Safety Standards:

UL 61010-1:2012

CSA C22.2#61010-1-12:2012

UL 61010-2-010:2015

CSA C22.2#61010-2-010:2015

Construction in Accordance with the following EMC Standards:

EN 61326-1

VI. Inspection

6.1 Packing List

Unpack the equipment carefully and check for any damage which may have arisen during transport. In the event of identified damage, please contact serviceusa@hollandgreenscience.com

| Type | Item Description | Quantity |
|------------------|------------------------------|----------|
| Main Unit | Vacuum Oven | 1 |
| Fittings | Vacuum Connection Pipe | 1 |
| Fittings | Spanner | 2 |
| Quick Setup Card | Scan to download user manual | 1 |




Note:






If you find any visible damage on the product, please do not connect the product to power supply.

VII. Operating Instructions

7.1 Temperature Controller Panel



Description of Controller Panel Keys

- 1)  Press and hold the key for 3s or more to control the run/stop of the program.
- 2)  In the setup mode, press this key to add one number and hold the key to increase numbers successively.
- 3)  In the setup mode, press this key to decrease one number and hold the key to decrease numbers successively.
- 4)  Return.
- 5)  Set the target temperature and enter the internal parameter layers.

7.2 Functions



- **F1:** first layer measurement temperature
- **F2:** second layer measurement temperature
- **F3:** third layer measurement temperature
- **F4:** fourth layer measurement temperature
- **RUN/STOP:** display of running status
- **14: 09:** time display of controller system
- **SET:** display of set value. (°C: display of temperature unit symbol)




7.3 Priming the Vacuum Pump

- Priming of the vacuum pump requires that the connecting tube provided is securely (airtight sealed) placed between the vacuum pump and the vacuum oven.
- The oven door should be securely shut, and the door handle secured in a closed position (this will prevent air from entering the chamber by causing the hole in the rubber plug and the hole in the pipe core to be 180 ° out of phase).
- Power to the vacuum pump should be turned on and then the vacuum valve should be closed by turning the handle clockwise 90° during processing.
- To stop the process, reduce the vacuum in the chamber -0.1Mpa before turning off the vacuum pump. This will ensure that no engine oil from the vacuum pump will be drawn into the chamber.



7.4 Priming the Vacuum Oven

- Turn on power to the vacuum oven (a power-on indicator light will confirm this) and after a 5 second self-test is completed; the temperature controller will automatically enter the “run” mode. This will allow the PV screen to display the process temperature and the SV screen display the set temperature.
- Once the process begins the AT indicator light will turn on and if process temperature is less than the set temperature heating lamps will turn on.


7.5 Setting of Temperature

- When the instrument is switched on, press the  button one time, and the “set temperature value” will start to blink. Set the desired temperature value (in Celsius degrees) by pressing keys.
- Use the  button to move between the digits.
- Confirm the set value by pressing the  button.


Note:

Pressing the  or  key once changes the value by 0.1.
Press and hold the keys changes the value rapidly.


7.6 Start / Stop Operation

- After setting the operating parameters, press  button for 4-5 seconds to start the heating process. When the word “STOP” appears on the upper right corner of the screen, the system is in the stopped state. When the message “RUN” appears in the

upper right corner, the system is in the running state and the screen will display values like timer, temperature measured inside the chamber.


- At any time, you can always manually stop running by pressing the  button for 4-5 seconds.

7.7 Setting the Deviation Alarm and Timer




- On the controller panel, press and hold the  key for 4-5 seconds.



- Press  key to confirm, press  or  key to move, press  key to return.




- TIMER WAIT: ON/OFF of waiting temperature. When 'ON' the waiting temperature is available. When 'OFF' no waiting temperature is available.
- ALM Temp: Setting of temperature deviation alarm: when the measured temperature > set temperature +ALM TEMP, an audible alarm will sound.
- RH Temp: Setting of the upper limit for Temperature Control.
- Timer: When the time is set as 0, timing function is not available. When other than 0 is applied the controller is in 'Timer' function.

Press the  key so, that Timer numerical value is back lit displayed, indicating the timer can now be set as required, by pressing either the  or  key.

NOTE: If no operation is performed for a long time, the system automatically returns to the home screen.









7.8 Setting the Pre-set Start Function

- Press and hold the  key and press  key, then shift the cursor to 'START SETUP', and press  key.




- 1) START TIMER: Press 'ON' to start the Pre-Set timer function settings
- 2) YEAR, MONTH, DAY, HOUR, MINUTE: Enter timer pre-set required details

7.9 Setting the ENGINEER SETUP


- Press and hold the  key and press  key, then shift the cursor to 'ENGINEER SETUP', and press  key, enter password (The default is "0") and press  to confirm. Press  key to move, press  key to confirm, and then press  or  key to adjust.

NOTE: The parameters here are factory parameters, do not adjust (except for the system date and password)



7.10 Unit, SYSTM Timer and PASSWORD Settings

- Unit: "0" °C, "1" °F
- YY: year; MM: month; DD: day; HH: hour; MM: minute; SS: second; (this is the apparatus system time setting)
- PASSWORD: for direct input and changing of new password.



VIII. Cautions

- Combustible materials and corrosive liquids or gases should not be placed in the vicinity of the drying cabinet.
- If the product is excessively damp before drying, add a “precipitation” filter between the vacuum pump and the cabinet to avoid steam entering the pump.
- If a granular product is to be dried, a 5 μm filter should be added in series before the vacuum intake.
- If the air release rubber stopper becomes difficult to turn, Vaseline may be used to lubricate the contact area.
- The vacuum chamber should always be clean. The glass on the cabinet door should only be cleaned using a soft cotton cloth.
- Whenever the vacuum cabinet is not used for a prolonged period of time, it is recommended that all exposed plated parts be cleaned by nonabrasive lubricant. It is recommended that a plastic film cover be placed over the cabinets to prevent dust from accumulating.

X. Fault Diagnosis

| Fault | Possible Causes | Solution |
|---|---|--|
| No Power | External power supply outlet has no power. | Check power outlet |
| | The power plug is not correctly inserted into the socket | Re-insert the power plug into the socket check switched 'ON' |
| | Failed fuse | Check the fuse. Replace the fuse. |
| PV Display | Temperature sensor (Pt100) is damaged | Check Pt100. (0°C, 100Ω, 0.3Ω/°C) Replace if required |
| | Temperature sensor line issues encountered. | Re-connect the lines again. |
| Temperature Not Increasing | The SET value is too low | Set temperature SV ≥ ambient +10°C |
| | Temperature controller has no output | Replace the temperature controller |
| | The heater is damaged (short circuit, or open circuit) | Replace the heater |
| | Timing function or the setting is not correct. | ST=0 or ST= (heating time +constant temperature time) |
| Temperature Issues * Out of Control * Offset or Overshoot * Error Display Tested Temperature vs. Real Temperature | The output of the temperature controller is unreliable | Replace 3041 or BTA |
| | Pt sensor not connecting well. | Reconnect PT sensor |
| Large Variance Between Test Temperature and Real Temperature. | Relevant parameters are not correctly set. | Re-set all relevant parameters, such as Ar, P etc. |
| | Insufficient vacuum | Check vacuum pump |
| | The parameters within the vacuum oven are not accurate or correct | Re-arrange Pb, and Pk parameters. |

| | | |
|--|--|---|
| Unable to maintain vacuum within chamber | The vacuum pump is not the correct size for the volume within the vacuum oven | The vacuum pumping speed should not be less than 2 Litres per second |
| | Various connecting flanges and pipes are not secured or coming loose | Check all are securely fitted if not replace applicable fittings. |
| | The vacuum gauge is damaged. | Replace the vacuum gauge. |
| | The oven door is not closing well. | Adjust the door pin distance. |
| | The airproof rubber door seal is aged and lacks elasticity. | Replace the airproof rubber door seal |
| | Air release valve and vacuum valve are not in the correct place. | Manually adjust positioning |
| | The electromagnetic valve is damaged. (Models 6090 & 6210) | Replace electromagnetic valve |
| Air Leakage Vacuum Degree Decreased to 0.092 Mpa from 0.1Mpa within 24 hours. | There is air leakage in various connect pipes and flanges | Check all and replace if required |
| | Heat incurred from the oven can distort the shape of the "O" ring causing air leakage. | Screw tight the heater seat (Located in the back of the inner bladder) or replace the "O" ring. |
| | The air release valve is not in the right place. | Manually adjust into the correct position |
| | Air leakage within the vacuum valve. | Replace the vacuum valve |
| | The electromagnetic valve cannot be closed and there is air leakage. | Replace the electromagnetic valve. |

XI. Maintenance and Cleaning

Proper maintenance and operation of the vacuum oven is required to ensure good working condition and extend the oven's service life.

Keep the vacuum oven dry and clean in routine operations. Clean the exterior of the vacuum oven with a non-abrasive cleanser and only connect the power supply when the oven surface is dry.

If liquid or moist solids are spilled the vacuum oven, please immediately disconnect the power supply and contact the service department.

- Power must be disconnected before any maintenance or cleaning.
- Surface stains on the vacuum oven should be cleaned only by a clean, soft rag and detergent.
- Keep the product clean and ensure no cleaning solution is spilled inside the vacuum oven.
- Avoid cleaning the vacuum oven with any corrosive cleaning solutions.
- If the vacuum oven is left unused for a long period, switch 'OFF' the power and store within a dry and clean area on a level surface at room temperature.