BIOS
Diagnostics

1 SECOND

EAR

THERMOMETER

Instruction Manual

235DI

1 Second Ear Thermometer

Instruction Manual

Table of Contents

- | | |
|----|---|
| 1 | Important Safety Instructions |
| 2 | Product Diagram |
| 3 | How this Thermometer Measures Temperature |
| 4 | Directions for Use |
| 5 | Setting Date and Time |
| 6 | Setting the Beeper |
| 7 | Fever Alarm |
| 8 | Changing Between Celsius and Fahrenheit |
| 9 | Memory Storage |
| 10 | Measuring Surface Temperatures |
| 11 | Care and Cleaning |
| 12 | Battery Replacement |
| 13 | Error Messages |
| 14 | Technical Specifications |
| 15 | 5 Year Warranty |

1 Second Infrared Ear Thermometer

Thank you for purchasing the BIOS Diagnostics 1 Second Ear Thermometer.

The BIOS Diagnostics 1 Second Ear Thermometer is a high quality product incorporating the latest technology and tested in accordance with international standards.

With its unique technology, this thermometer can provide safe, accurate and reliable readings with each measurement on people over 2 years of age.

Features

Multiple Uses (Wide Range Measurement)

This thermometer offers a wide measurement range feature from 0°C to 100.0°C (32.0°F to 212.0°F), meaning the unit can be used as an ear thermometer to measure body temperature or it can be used to measure surface temperature of the following:

- Milk surface temperature in a baby's bottle
- Surface temperature of a baby's bath
- Ambient temperature

Measurement in 1 Second

The innovative infrared technology allows measurement of ear temperature in only 1 second.

Multiple Reading Recall

Users will be able to recall the last 30 readings when entering the recall mode, enabling efficient tracking of temperature variations.

Fever Alarm

10 short beeps and a red LCD backlight alerts the patient that he/she may have a temperature above 37.4°C / 99.3°F.

Illuminated Probe

This thermometer includes a probe LED light which enables the user to find the correct ear position in the dark. The probe will illuminate for approximately 10 seconds when the thermometer is turned ON.

"Clean Me" Technology

After each temperature measurement "clean me" will flash on the LCD to remind users to clean the probe. See Care and Cleaning section.

Accurate and Reliable

The unique probe assembly construction incorporates an advanced infrared sensor, ensuring that each measurement is accurate and reliable.

Please read these instructions carefully before using this thermometer and keep both in a safe place.

1. IMPORTANT SAFETY INSTRUCTIONS

This instrument may only be used for the purposes described in this booklet. The manufacturer cannot be held liable for damage caused by incorrect application.


Never immerse this instrument in water or other liquids. For cleaning please follow the instructions in the «Care and Cleaning» section.

- Do not use the instrument if there are any signs of damage. Do not attempt to repair the instrument.
- Earwax in the ear canal may cause a lower temperature reading, therefore it is important to ensure the patient's ear canal is clean.
- If the measurement result is not consistent with the patient's finding or suspiciously low, repeat the measurement every 15 minutes or double check the result by another core body temperature measurement.
- This thermometer consists of high quality precision parts. Do not drop this instrument. Protect it from severe impact and shock.
- Keep the thermometer away from direct exposure to the sun and keep it in a dust-free dry area.
- If the instrument is not going to be used for a prolonged period the battery should be removed.

WARNING:

- Use of this thermometer is not intended as a substitute for consultation with your physician.
- Ensure that children do not use the instrument unsupervised; some parts are small enough to be swallowed.
- It is not recommended to use infrared thermometers on children under 2 years of age.
- Infrared thermometers should not be used to manage important health concerns, and therefore parents should not rely solely on temperature readings to evaluate the status of your child's health. If you have any concerns, consult with your doctor.

2. PRODUCT DIAGRAM

- | | |
|------------------------|------------------------------|
| 1. Measuring sensor | 4. ON/OFF button |
| 2. Display | 5. Battery compartment cover |
| 3. START button | 6. M -button (memory) |

3. HOW THIS THERMOMETER MEASURES TEMPERATURE

This thermometer measures infrared energy radiated from the eardrum and the surrounding tissue. This energy is collected through the lens and converted to a temperature value. The measured reading obtained directly from the eardrum (Tympanic Membrane) ensures the most accurate ear temperature.

Measurements taken from the surrounding tissue of the ear canal generate lower readings and may result in misdiagnosis of fever.

4. DIRECTIONS FOR USE

IMPORTANT: If using the thermometer for the first time, remove the plastic strip located in the battery compartment.

Figure 1

1. Switch ON the thermometer by pressing the **ON/OFF** button, all segments will be shown for 1 second. **(Figure 1)**

2. After one beep is heard and the °C / °F icon is flashing, the thermometer is ready for measurement. Straighten the ear canal by gently pulling the middle of the ear back and up.
3. Place the probe firmly into the ear canal, press the **START** button and keep the probe in the ear until the thermometer beeps to signal the completion of the measurement. The probe LED light will turn OFF once the **START** button is pressed to take the measurement.

NOTE:

- Always insert the thermometer in the same ear with the same direction and depth.
- Always clean the probe accordingly before and after each use. (See Care and Cleaning Section).
- It is recommended that you measure the temperature 3 times in the same ear. If the 3 measurements are different, select the highest one. (Allow 30 seconds between each measurement).
- Before measurements, remain in a stable environment and avoid exercise and bathing for 30 minutes.
- Readings from different measuring sites should not be compared as the normal body temperature varies by measuring site and time of day, being highest in the evening and lowest about one hour before waking up.

Normal body temperature ranges:

- Axillary: 34.7 – 37.3°C / 94.5 – 99.1°F
- Oral: 35.5 – 37.5°C / 95.9 – 99.5°F
- Rectal: 36.6 – 38.0°C / 97.9 – 100.4°F
- Ear: 35.4 – 37.4°C / 95.7 – 99.3°F

5. SETTING DATE AND TIME

1. Once a new battery is installed or by pressing and holding the **M**-button for 3 seconds the year will flash on the display. You can change the year by pressing the **START** button. To confirm the year press the **M**-button. The month will flash.
2. Press the **START** button to set the month digits. Press the **M**-button to confirm.

3. Follow the steps above to set the day, 12/24 hour mode, hours and minutes.

4. Once you have set the minutes and pressed the **M**-button the date and time are set and the time is displayed.

- If no button is pressed for 20 seconds, the thermometer automatically switches to the measuring screen.
- To cancel time setup, press the **ON/OFF** button during time setup. The LCD will show Date / Time icons with "--:--". Press the **ON/OFF** button to start the measurement. If no further action is taken within 60 seconds, the thermometer will automatically turn OFF.

6. SETTING THE BEEPER

1. When the thermometer is OFF, press and hold the **ON/OFF** button for 5 seconds to set the beeper. **(Figure 2)**
2. Press the **ON/OFF** button again to either turn the beeper ON or OFF. The beeper is activated when the beeper icon is not crossed out.

Figure 2

NOTE: If no button is pressed for 5 seconds, the thermometer automatically switches to measuring mode.

7. FEVER ALARM

If the thermometer detects a body temperature greater than 37.4°C/99.3°F, it will make 10 short beeping sounds and the LCD will light up RED, indicating a potential fever.

8. CHANGING BETWEEN CELSIUS AND FAHRENHEIT

This thermometer can display temperature measurements in either Celsius or Fahrenheit. To switch the display between °C and °F, simply turn OFF the unit, press and hold the **START** button for 5 seconds; after 5 seconds, the current measurement scale (°C or °F icon) will flash on the display **(Figure 3)**. Change the measurement scale between °C and °F by pressing the **START** button. When the measurement scale has been chosen, wait for 5 seconds and the unit will automatically enter the ready for measuring mode.



Figure 3

9. MEMORY STORAGE

This thermometer recalls the last 30 readings.

- **Recall mode:** Press the **M**-button to enter Recall mode when the power is OFF. The memory icon << M >> will flash. **(Figure 4)**
- **Reading 1 - the last reading:** Press and release the **M**-button to recall the last reading. Number 1 and a flashing M are displayed. **(Figure 5)**
- **Reading 30 - readings in succession:** Press and release the **M**-button consecutively to recall the last 30 readings in succession. Pressing and releasing the **M**-button after the last 30 readings have been recalled will resume the above sequence from reading 1.


Figure 4

Figure 5

10. MEASURING SURFACE TEMPERATURES

1. Press the **ON/OFF** button, after one beep sounds the thermometer is ready for use.
2. Aim the thermometer probe at the desired surface and press the **START** button.

11. CARE AND CLEANING

Use an alcohol swab or cotton tissue moistened with alcohol (70% Isopropyl) to clean and disinfect the thermometer casing and the measuring probe. Ensure that no liquid enters the interior of the thermometer. Never use abrasive cleaning agents, thinners or benzene for cleaning and never immerse the instrument in water or other cleaning liquids. Take care not to scratch the surface of the probe lens and the display.

12. BATTERY REPLACEMENT


This thermometer is supplied with one lithium battery type, CR2032.

The battery needs replacing when only the ▼ icon is shown on the display. **(Figure 6)**





Remove the battery cover by sliding it in the direction shown. Remove the battery and replace with a new one.

Figure 6

13. ERROR MESSAGES

| Display | Error |
|---|---|
|  | <p>Measured temperature too high: Displays “H” when measured temperature is higher than 100°C / 212°F</p> |
|  | <p>Measured temperature too low: Displays “L” when measured temperature is lower than 0°C / 32°F</p> |
|  | <p>Ambient temperature too high: Displays “H” in conjunction with the “▼” when ambient temperature is higher than 40°C / 104°F</p> |
|  | <p>Ambient temperature too low: Displays “L” in conjunction with the “▼” when ambient temperature is lower than 10°C / 50°F</p> |
|  | <p>Error function display: The system has a malfunction</p> |
|  | <p>Blank display: Check if the battery has been loaded correctly. Also check polarity of the battery.</p> |
|  | <p>Flat battery indication: If only the “▼” icon is shown on the display, the battery should be replaced immediately.</p> |

14. TECHNICAL SPECIFICATIONS

| | |
|-----------------------------|---|
| Type | 1 Second Ear Thermometer 235DI |
| Measurement Range | 0°C to 100.0°C/32°F to 212.0°F |
| Resolution | 0.1°C/°F |
| Measurement Accuracy | Laboratory: ± 0.2°C, 32°C to 42.2°C/ ± 0.4°F, 89.6°F to 108.0°F |
| Display | Liquid Crystal Display , 4 digits plus special icons |

| | |
|--------------------------------------|--|
| <p>Acoustic</p> | <ul style="list-style-type: none"> • The unit is turned ON and ready for the measurement: 1 short beep • Complete the measurement: 1 long beep (1 second) if the reading is less than 37.5°C/99.5°F, 10 short beeps sound, if the reading is equal to or greater than 37.5°C/99.5°F • System error or malfunction: 3 short beep sounds |
| <p>Memory</p> | <ul style="list-style-type: none"> • Auto-Display the last measured temperature • 30 readings recall in the Memory Mode |
| <p>Backlight</p> | <ul style="list-style-type: none"> • The display light will be GREEN for 1 second, when the unit is turned ON. • The display light will be GREEN for 5 seconds, when a measurement is completed with a reading less than 37.5°C/99.5°F • The display light will be RED for 5 seconds, when a measurement is completed with a reading equal to or higher than 37.5°C/99.5°F |
| <p>Operating Temperature</p> | <p>10°C to 40°C/50°F to 104°F 15% - 95% relative maximum humidity</p> |
| <p>Storage Temperature</p> | <p>10°C to 40°C/50°F to 104°F 15% - 95% relative maximum humidity</p> |
| <p>Automatic Switch-Off</p> | <p>Approximately 1 minute after last measurement has been taken.</p> |
| <p>Battery</p> | <p>CR2032 Battery 3V - at least 1,000 measurements</p> |
| <p>Dimensions</p> | <p>139 (L) x 39 (W) x 42 (H) mm</p> |
| <p>Weight</p> | <p>54.5 g (with battery), 51 g (without battery)</p> |
| <p>Reference to Standards</p> | <p>EN 12470-5; ASTM E1965; IEC 60601-1; IEC 60601-1-2 (EMC) This device complies with the requirements of the Medical Device Directive 93/42/EEC.</p> |

15. 5 YEAR WARRANTY

This product is warranted to be free from defects in material or workmanship for 5 years. It will be repaired or replaced without charge upon receipt of the unit prepaid with \$5.00 to cover handling, packaging and return postage to the address below. Please include your full name, address and daytime phone number or email address.

16975 Leslie Street
Newmarket, ON
L3Y 9A1
www.biosmedical.com

 Read the instruction manual carefully before using this device, especially the safety instructions, and keep the instruction manual for future use.

 Type BF applied part

 Batteries and electronic devices must be disposed of in accordance with the locally applicable regulations, not with domestic waste.

Using the unit in the immediate vicinity of mobile phones, microwave appliances or other devices with strong electromagnetic fields may result in impaired functioning.

IP22 Protected against solid foreign objects of 12.5 mm diameter and greater. Protected against vertically falling water drops when the device is tilted up to 15°.

Expected service life: 5 years

Do not use this device close to strong electromagnetic fields, such as mobile telephones or radio installations. Keep a distance from such devices when using this unit.

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