Project Title
Instructions

UPC (For Reference Only)

8 106841 003609 15

Date: March 3 2020
Customer: Thermor Ltd
Job Number: TH6128
Item Number: MD300C29
Designer: Alex Vranjesjevic

Reason For Project:
New Item:
New Branding:
Design Update:
N/A X

Artwork is not to be amended or scaled. If any changes are required please send through Thermor office with instructions, and we will amend here and resend artwork through. **PLEASE RETURN WITH SIGNATURE**

Colour Breakdown

Contact: Graphic's Department
graphic@thermor-us.com

PO: 24086-00
Revision: A

Thermor

Contact: Graphic's Department
graphic@thermor-us.com

Colour Breakdown

LIMITED WARRANTY

The Fingertip Pulse Oximeter has a 2 year warranty to be free of manufacturing defects in materials and workmanship under normal applications for 2 years of the original owner. If this product becomes inoperable due to defect and requires repair, return the product with all component pieces and proof of purchase to the address listed below. This warranty does not cover any shipping/transport costs. This warranty does not apply if the product is subject to misuse, neglect, rough handling or damage.

Ship the unit prepaid and insured (at owner's option) to:
Thermor Ltd.
Attn: Repair Department
16975 Leslie Street
Newmarket, ON L3Y 9A1

www.biosmedical.com
Email: support@biosmedical.com

Limitations

This warranty covers all defects encountered in normal use of the equipment and does not apply in the following cases:

- 1) If the equipment has been serviced by other than a certified center.
- 2) Damage to the equipment due to mishandling, abuse, accident or not following operating instructions.
- 3) INSTRUMENTS: Warranty does not extend to display face, batteries, carrying case or lanyard.

For troubleshooting assistance please contact Consumer Support. Our representatives are trained to provide you assistance over the phone and may be able to resolve your problem without returning. For the best possible assistance please have your unit available when calling us.

Service Phone: 1 (866) 517-2970

www.biosmedical.com

OPERATOR'S MANUAL

MD300C29 Fingertip Pulse Oximeter

General Description

Oxygen binds to hemoglobin in red blood cells when moving through the lungs. It is transported throughout the body as arterial blood. A pulse oximeter uses two frequencies of light (red and infrared) to determine the percentage (%) of hemoglobin in the blood that is saturated with oxygen. The percentage is called blood oxygen saturation, or SpO₂. A pulse oximeter also measures and displays the pulse rate at the same time it measures the SpO₂ level.

Diagram of Operation Principle

1. Red and infrared-ray Emission Tube
2. Reirs or Infrared-ray Receipt Tube

Precautions For Use

1. Before use, carefully read the manual.
2. Operation of the fingertip pulse oximeter may be affected by the use of an electrical unit (ESU).
3. The fingertip pulse oximeter must be able to measure the pulse properly to obtain an accurate SpO₂ measurement. Verify that nothing is hindering the pulse measurement before relying on the SpO₂ measurement.
4. Do not use the fingertip pulse oximeter in an MRI or CT environment.
5. Do not use the fingertip pulse oximeter in situations where alarms are required: The device has no alarms. It is not for continuous monitoring.
6. Do not use the fingertip pulse oximeter in an explosive atmosphere.
7. The fingertip pulse oximeter is intended only as an adjunct in patient assessment. It must be used in conjunction with other methods of assessing clinical signs and symptoms.
8. In order to ensure correct sensor alignment and skin integrity, the maximum application time at a single site for our device should be less than half an hour.
9. Do not sterilize the device using autoclaving, ethylene oxide sterilizing, or immersing the device in liquid. The device is not intended for sterilization.
10. Follow local ordinances and recycling instructions regarding disposal or recycling of the device and device components, including batteries.
11. This equipment complies with IEC 60601-1-2:2014 for electromagnetic compatibility for medical electrical equipment and/or systems. However, because of the proliferation of radio-frequency transmitting equipment and other sources of electrical noise in healthcare and other environments, it is possible that high levels of such interference due to close proximity or strength of a source might disrupt the performance of this device.
12. Portable and mobile RF communications equipment can affect medical electrical equipment.
13. This equipment is not intended for use during patient transport outside the healthcare facility.
14. This equipment should not be used adjacent to or stacked with other equipment.
15. It may be unsafe to:
 - disassemble, repair or modify the equipment
 - use accessories, detachable parts and materials not described in the instructions for use
 - interconnect this equipment with other equipment not described in the instructions for use
 - disassemble, repair or modify the equipment
16. These materials that contact with the patient's skin contain medical silicone and ABS plastic enclosure are all pass the ISO10993-5 Tests for invitro cytotoxicity and ISO10993-10 Tests for irritation and delayed-type hypersensitivity.

Inaccurate measurements may be caused by

1. Significant levels of dysfunctional hemoglobin (such as carboxyl - hemoglobin or methemoglobin).
2. Intravascular dyes such as indocyanine green or methylene blue.
3. High ambient light. Shield the sensor area if necessary.
4. Excessive patient movement.
5. High-frequency electro-surgical interference and defibrillators.
6. Venous pulsations.
7. Placement of a sensor on an extremity with a blood pressure cuff, arterial catheter, or intravascular line.
8. The patient has hypotension, severe vasoconstriction, severe anemia, or hypothermia.
9. The patient is in cardiac arrest or is in shock.
10. Fingernail polish or false fingernails.
11. Weak pulse quality (low perfusion).
12. Low hemoglobin.

Product Features

1. High brightness OLED display SpO₂, PR, Pulse bar and Plethysmograph.
2. Low power consumption; battery-low indicator
3. Automatically power off; 2 AAA alkaline batteries
4. When no signal or low signal is detected, the pulse oximeter will power off automatically in 8 seconds.

Contraindication

It is not for continuous monitoring

Intended Use

The Fingertip Pulse Oximeter is a handheld non-invasive device intended for spot-checking of oxygen saturation of arterial hemoglobin (SpO₂) and Pulse Rate of adult, adolescent and child patients in hospitals, hospital-type facilities and homecare.

Operation Instructions

1. Install two AAA batteries into battery compartment correctly.
2. Place clamp over finger nail as the right diagram
3. Insert one finger into rubber hole of the oximeter fully.
4. Press the switch once on the front panel.
5. Finger and body should not tremble during measuring.
6. Read correspond data from display screen.
7. Two display modes

After turning on the oximeter, each time you press the power switch, the oximeter will switch to another display mode. There are 2 display modes shown as follows:

When you press the power switch for more than one second, the brightness of the oximeter will be changed by degrees. There are 2 levels on brightness; the default is level 1.

NOTE: Please use medical alcohol to clean the rubber inside the oximeter, and clean the test finger using alcohol before and after each test. (The rubber inside of the oximeter belongs medical rubber, which has no toxin and no harmful to the skin of human being).

Product Accessories

1. One lanyard
2. Two AAA batteries
3. One user's manual
4. One carrying case

Battery Installation

1. Install two AAA batteries into the battery compartment. Match the plus (+) and minus (-) signs in the compartment. If the polarities are not matched, damage may be caused to the oximeter.
2. Slide the battery door cover horizontally along the arrow.

Notes:
Please remove the batteries if the pulse oximeter will not be used for long periods of time.

Maintenance and Storage

1. Replace the battery regularly in a timely manner when low voltage lamp is lighted.
2. Clean the surface of the fingertip oximeter before it is used in diagnosis for patients.
3. Remove the batteries if the oximeter is not operated for a long time.
4. It is best to store the product at 25°C~+70°C and 55%RH.
5. Keep in a dry place. Extreme moisture may affect oximeter lifetime and may cause damage.
6. Dispose of battery properly; follow any applicable local battery disposal laws.

Cleaning the fingertip pulse oximeter

Use 70% medical alcohol to clean the silicone touching the finger inside of oximeter with a soft cloth dampened with 70% isopropyl alcohol. Also clean the being tested finger using alcohol before and after each test. Do not pour or spray liquids onto the oximeter, and do not allow any liquid to enter any openings in the device. Allow the oximeter to dry thoroughly before reuse.

The fingertip pulse oximeter requires no routine calibration or maintenance other than replacement of batteries.

The use life of the device is five years when it is used for 15 measurements every day and 10 minutes per one measurement. Stop using and contact local service center if one of the following cases occurs:

- An error in the Possible Problems and solutions is displayed on screen.
- The oximeter cannot be powered on in any case and not the reasons of battery.
- There is a crack on the oximeter or damage on the display resulting readings cannot be identified; the spring is invalid; or the key is unresponsive or unavailable.

Disinfecting

The applied parts touching the patients' body are required to be disinfected once after each use. The recommended disinfectants include: ethanol 70%, isopropyl alcohol 70%, glutaraldehyde 0.5%, and disinfectant. Disinfection may cause damage to the equipment and is therefore not recommended for this pulse oximeter unless otherwise indicated in your hospital's servicing schedule. Clean the pulse oximeter before disinfecting it. CAUTION! Never use EO or formaldehyde for disinfection.

Specifications

1. Display Type:	OLED display
2. SpO₂	Measurement range: 70%-100% Accuracy: 70%~100%; ±2%; 0%-69% no definition Resolution: 1%
3. Pulse Rate	Measurement range: 30bpm~250bpm Accuracy: 30bpm~99bpm, ±2bpm; 100bpm~250bpm, ±2% Resolution: 1bpm
4. Probe LED Specifications	
	Wavelength Radiant Power
RED	660±3nm 3.2mW
IR	905±10nm 2.4mW

NOTE: The information about wavelength range can be especially useful to clinicians.

5. Power Requirements

Two AAA alkaline Batteries
Power consumption: Less than 40mA
Battery Life: Two AAA 1.5V, 1200mAh alkaline batteries could be continuously operated as long as 18 hours.
6. Environment Requirements
Operation Temperature: 5°C ~ 40°C
Storage Temperature: -25°C ~ +70°C
Ambient Humidity: 15%~ 93% no condensation in operation; 93% no condensation in storage/transport

Atmosphere pressure: 70kPa~106kPa

7. Equipment Response Time

As shown in the following figure (Figure 1) Response time of slower average is 8s.

8. Classification

According to the type of protection against electric shock: INTERNALLY POWERED EQUIPMENT;

According to the degree of protection against electric shock: TYPE BF APPLIED PART; (applied part: the rubber hole of the device);

According to the degree of protection against ingress of water: IP22


According to the mode of operation: CONTINUOUS OPERATION

Clinical Study Summary

The following details are provided to disclose actual performance observed in the clinical validation study of healthy adult volunteers. The ARMS value analysis statement and Bland-Altman plot of data is shown as follows:

Item	90~100	80~90	70~80
#pts	78	66	63
Bias	1.02	0.40	-0.48
ARMS	1.66	1.46	1.93

Possible Problems and Solutions

Problems	Possible reason	Solution
SpO ₂ or PR can not be shown normally	1. Finger is not inserted correctly 2. Patient's SpO ₂ value is too low to be measured	1. Retry by inserting the finger 2. There is excessive illumination 3. Try some more times. If you can make sure no problem exist in the product, please go to a hospital timely for exact diagnosis.
SpO ₂ or PR is shown unstably	1. Finger might not be inserted deep enough. 2. Excessive patient movement	1. Retry by inserting the finger 2. Be calmness
The oximeter can not be powered on	1. No battery or low power of battery 2. Batteries might be installed incorrectly 3. The oximeter might be damaged	1. Please replace batteries 2. Please reinstall the batteries 3. Please contact with local customer service centre
Indication lamps are suddenly off	1. The product is automatically powered off when no signal is detected longer than 8 seconds 2. The battery power is too low to work	1. Normal 2. Replace the batteries
"Err7" is displayed on screen	Err7 means all the emission LED or reception diode is damaged	Please contact with local customer service centre

Symbol Definitions

Symbol	Definition	Symbol	Definition
	Type BF applied part.		No SpO ₂ Alarm
	Follow instruction for use		Power switch
	Protected against dripping water.		Lot No.
%SpO2 symbol"/>	Oxygen saturation		Storage temperature and relative humidity
	Pulse Rate (bpm)		Manufacturer's information
	Low power indication		Date of Manufacture
	Attention		Conformity to WEEE Directive

- Notes:**
1. The illustrations used in this manual may differ slightly from the appearance of the actual product.
 2. The specifications are subject to change without prior notice.

Declaration

Requirement - Test	Result/Comments	Verdict
Class 7 - Emissions	---	---
Classification	---	---
Class A or B	Class B	---
Group 1 or 2	Group 1	---
CISPR 11, 14-1, 32 or ISO 7137	CISPR 11	---
Conducted RF Emissions	N/A	N/A
Radiated RF Emissions	---	P
Distance Power (if applicable)	N/A	N/A
Harmonic Distortion per IEC61000-3-2 (Class A, B, C, D)	N/A	N/A
Voltage Fluctuations and Flicker per IEC1000-3-3	N/A	N/A
Class 8 - Immunity	---	---
Electrostatic Discharges	IEC 61000-4-2	P
Radiated RF EM Fields and Proximity Wireless Fields	IEC 61000-4-3	P
Electrical Fast Transients and bursts	IEC 61000-4-4	N/A
Burges	IEC 61000-4-5	N/A
Conducted Disturbances, induced by RF fields	IEC 61000-4-6	N/A
Voltage Dips and Interruptions	IEC 61000-4-11	N/A
Rated Power-Frequency Magnetic Field	IEC 61000-4-8	P

Distributed by:
Thermor Ltd.
16975 Leslie Street
Newmarket, ON L3Y 9A1
WWW.BIOSMEDICAL.COM

Manufactured by:
Beijing Choice Electronic Technology Co., Ltd.
Room 4104, No. A12 Yuquan Road, Haidian District, Beijing P.R. China.

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


Quick Operation Guide


ChoiceMMed_MD300C29_Ver2.0

- Before using the oximeter, please remove protective sticker covering the display.


- 1 Displaying Screen
- 2 Power Key
- 3 Shaft Cover
- 4 Battery Box Cover

Brief Description of Front Panel

The pulse bar graph displays corresponding with the patient's pulse beat. The height of the bar graph shows the patient's pulse strength.

Battery Installation


Please put in the batteries according to the polarity marked inside the battery compartment.

Using the Lanyard

1. Thread thinner end of the lanyard through the loop.
2. Thread thicker end of the lanyard through the threaded end before pulling it tightly.


Warnings!

- Keep the oximeter away from young children. Small items such as the battery door, battery, and lanyard are choking hazards.
- Do not hang the lanyard from the device's electrical wire.

Operation Instruction

Place one of your fingers into the pulse oximeter to the end and press the switch button one time on front panel to turn it on.

There are 2 display modes. After turning on the pulse oximeter, each time you press the power switch, the pulse oximeter will switch to another display mode.

When you press the power switch for more than one second, the brightness of the pulse oximeter will be changed by degrees. There are 2 levels on brightness. The default level is level 1.

Keep your hands still for the reading. The pulse oximeter will power off automatically in 8 seconds if there's no finger inside.

Warnings and Notes

Warnings:

1. Keep the pulse oximeter away from young children. Small parts such as the battery door and the batteries etc. May be hazardous if swallowed.
2. The lanyard may cause strangulation in conditions that may cause it to twist around the neck.

Notes:

1. Read the manual carefully before use.
2. The illustration used in this manual may differ slightly from the appearance of the actual product.
3. Follow local ordinances and recycling instructions regarding disposal or recycling of the device and device components, including batteries.

Frequently Asked Questions

What is a Pulse Oximeter?

A pulse oximeter is a non-invasive device that indirectly monitors the blood oxygen saturation (SpO₂) and pulse rate (heart rate). It displays both blood oxygen saturation (SpO₂) & pulse rate (heart rate). Pulse oximeters provide an easy way of assessing your blood oxygen level and pulse rate.

What is SpO₂?

SpO₂ is also known as oxygen saturation. Oxygen saturation is a measure of how much oxygen the blood is carrying as a percentage of the maximum it could carry.

What is the normal range of SpO₂?

The normal range for SpO₂ is typically considered from 95%~99%. The SpO₂ measurement may be lower for people who live at high altitudes. Ask your health professional this question as it pertains to you.

What is the normal range for pulse rate?

The normal resting range for Pulse rate is typically considered from 60 ~100 beats per minute. Ask your health professional this question as it pertains to you.

What kind of conditions may cause an inaccurate reading?

Cold hands, poor circulation, very weak pulse, fingernail polish and acrylic nails may cause inaccurate results.

The SpO₂ is not changing - it's stuck?

SpO₂ does not change like pulse rate. It is slow to change.

The Pulse Rate is changing rapidly.

Your heart rate changes with emotions, excitement and exercise.

I do not see the battery light indicator.

The battery light indicator only appears when the battery is low.

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© 2019 ChoiceMed. ALL RIGHTS RESERVED