

## DrQuickLook™ Instructions

1. Remove unit from box.
2. **Turn ON DrQuickLook™ by sliding the small power switch on bottom of unit.**
3. Place DrQuickLook™ in docking/charger station on a flat surface/counter and plug in AC adapter, one end in to docking station and the plug end in to a standard AC outlet
4. Place DrQuickLook™ front side out in the docking/charger station, the charging indicator light on the front, directly below the 3.5" LCD monitor, will illuminate once contact is made indicating the unit is now in charging mode.
5. **Allow to charge for 8 to 10 hours when first removed from packaging.**

***THE AUTO SHUT OFF FEATURE*** will automatically shut off unit when placed flat down on a counter/table after one minute. The unit will also automatically shut off after one minute if held in a horizontal position. ***The unit will automatically turn back ON when returned back to a vertical position.*** AUTO SHUT OFF is instant when placed in docking/charger station. **To turn unit back on** remove from docking/charger station or pick up from counter/table and hold unit in an upright vertical position.

**NOTE:** When not in use, **ALWAYS** place DrQuickLook™ in the docking/charging station for safety **AND** to insure the unit is always charged. Batteries will last up to 4 hours of continuous use - so be sure to place unit in charger when not in use.

**Always use appropriate infection control plastic sleeves on the camera wand.**  
**Failure to do so can result in spread of contamination in the dental setting.**

***USING THE CAMERA WAND.*** Camera probe **BUTTON** features:

- A. **FREEZE MODE** Button on top of camera handle when pressed will freeze image temporarily, when pressed again image un freeze and will go back to normal video viewing mode
- B. **ZOOM MODE** Button on the side of the camera handle is used to zoom the frozen image when in *freeze mode* only. Press repeated for up to 3x *zoom mode*. Press again to return image back to normal image mode.

**Note:** When freeze image button is pressed and image is in *freeze mode*, LED lights will shut off. When freeze mode button is pressed again, image is in video mode and LED lights will be illuminated. When finished using, place camera wand back on top of the DrQuickLook unit, line up groove and allow the built-in magnets to hold in place.

Once a permanent location has been established for the DrQuickLook system, use the enclosed mounting screws to install DrQuickLook™ on the wall or cabinet surface near an AC plug.