Troubleshooting

General

■ Cannot Turn On the Device

- Be sure the power switch has been slid completely to the ON position.
- The battery is discharged. Charge the device sufficiently first.
- The battery is dead. This device does not have a serviceable battery. You may have to send the device back to the manufacturer for a replacement.
- Be sure the power supply is the one supplied with the unit. The plug has a sticker on the back indicating it is AC ADAPTER Model K25S075300U.
- Be sure the wall outlet is functional and the power supply is fully in place.
- Be sure the power supply is plugged completely into the charger base.

■ Cannot Access the SD Card

- The SD format may be incorrect. Insert the card and the unit should automatically reformat the SD Card. If this does not occur, discard the SD Card and try a fresh one.
- The SD Card may have been ejected but not removed from the device. Pop out the card and reinsert it.

■ Cannot Review Images

- Be sure the SD Card is correctly inserted into the device.
- If you have selected to eject the SD Card, you will need to remove the card and reinsert it.

■ Cannot Draw on Image

- The image is currently in Magnify mode. Press the Down button to switch off Magnify to start drawing.

■ The Device Seems Unresponsive

- The device will go to sleep after a period of no user interaction while lying flat. Place the device upright to wake it up.
- The battery may be discharged. Charge the device.

■ There is No Audio from the Speaker in Video Playback

- Volume may be set too low. Press the Up button to increase the volume.
- A headphones jack may be inserted in the device. Remove the headphones jack from the device.

Cameras

■ The Image Does Not Look Right

- The camera settings may be incorrectly set. Go to the Camera Settings to make necessary adjustments. Below are the default settings for brightness, hue and contrast.
- The viewer shows nothing but a blank screen. The intraoral camera wand may be ON but resting on the top of the viewer.

■ Cannot Capture Images

- Images cannot be transferred to the SD Card if the camera wand is in live survey mode. You must first **freeze the image** then **save to the SD Card**.
- Images cannot be saved if there is no SD Card in the device. Insert the SD Card into the device first. If you have selected to eject the SD Card, you will need to remove the card and reinsert it.

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