

SNAP N SAVE

The Snap N Save model consists of a seat with attached sliding rings, left and right leg frame, a backrest, two sliding rails, and four legs. This instruction sheet used for Models 61211, 61281, 61291, 61311, 61381, 61391, Weight capacity is 400 lbs for regular length benches only, 350 lbs. for medium length and 300 lbs. for extra long. The picture may differ from model to model but the structure is the same.



- 1 Press and hold down the push button brass locking pins on the sliding rails. Insert both sliding rails into the side leg frame (Fig. 1)
- 2 Stand it up vertically.
- 3 Determine which side of the transfer bench will face the shower head
- 4 Insert both sliding rails into sliding rings under the seat on each side. (Fig. 2)
- 5 Take the side leg frame with the handle bar and fit it onto the sliding rails as you did with the other side in Step 1. (Fig. 3)
- 6 Put the bench on the ground horizontally. Press the two brass locking pins of back brace tubes, then insert into tubes beneath the seat and push it until they click into place. (Fig. 4)
- 7 Insert each leg into the leg frame and set the same hole location on all four legs. Adjust the rubber tips on each leg to fully contact the floor, by turning them.

To Use:

- 1 Slide the bench seat all the way to the end outside of the tub (the end opposite to the handle bar). Press down from the middle of locking mechanism (located on the side of the seat) until it snaps onto the side brace. Then transfer the patient onto the seat. For added security, use the safety belt around the patient's lap. Release the locking mechanism and slide the patient with the seat all the way into the tub. Lock the seat again by pressing down the locking mechanism on the other side of the seat.
- 2 After completing the shower, unlock the locking mechanism and slide the patient with the seat out of the tub all the way to the opposite end. Push down the locking mechanism to lock the seat. Release the safety belt, then transfer the patient off the bench.

Maintenance:

To keep your sliding bench in optimal condition, we recommend to use WD-40, spray a little bit onto a damp cloth and wipe on the sliding rails. This will help keep the sliding smooth. Clean off any soap scum at least once a week and check to make sure all buttons are pop out and secure.

WARNING:

If any button does not pop out completely, or you feel the bench is not secure, do not use. Immediately contact the place of purchase.